
BrainFlow Documentation

Andrey Parfenov

Nov 13, 2020

CONTENTS

1	Supported Boards	3
1.1	Playback File Board	5
1.2	Streaming Board	5
1.3	Synthetic Board	6
1.4	OpenBCI	7
1.5	NeuroMD	13
1.6	G.TEC	15
1.7	Neurocity	16
2	Installation Instructions	19
2.1	Python	19
2.2	C#	19
2.3	R	19
2.4	Java	20
2.5	Matlab	20
2.6	Julia	20
2.7	Compilation of Core Module and C++ Binding	21
2.8	Android	22
3	User API	25
3.1	Python API Reference	25
3.2	C++ API Reference	39
3.3	Java API Reference	50
3.4	C# API Reference	65
3.5	R API Reference	80
3.6	Matlab API Reference	81
3.7	Julia API Reference	82
4	Data Format Description	83
4.1	Units of Measure	83
4.2	Generic Format Description	83
4.3	OpenBCI Specific Data	84
5	Code Samples	87
5.1	Python	87
5.2	Java	98
5.3	C#	108
5.4	C++	117
5.5	R	142
5.6	Matlab	145

5.7	Julia	148
5.8	Notebooks	153
6	Integration with Game Engines	161
6.1	Unity	161
6.2	Unreal Engine	165
7	BrainFlow Dev	167
7.1	Code style	167
7.2	CI and tests	167
7.3	Pull Requests	167
7.4	Instructions to add new boards to BrainFlow	168
7.5	Instructions to build docs locally	168
7.6	Debug BrainFlow's errors	168
7.7	BrainFlow Emulator	169
8	Ask Help	173
8.1	Contact Info, Feature Request, Report an Issue	173
8.2	Issue format	173
8.3	Contributors	173
9	Partners and Sponsors	175
9.1	OpenBCI	175
10	MIT License	177
11	Partners and Sponsors	179
12	Search	181
	Index	183

BrainFlow is a library intended to obtain, parse and analyze EEG, EMG, ECG and other kinds of data from biosensors. It provides a **uniform data acquisition API for all supported boards**, it means that you can switch boards without any changes in code and applications on top of BrainFlow are board agnostic. Also there is **powerful API to perform signal processing** which you can use even without BCI headset. Both of these two APIs are the same across bindings.

SUPPORTED BOARDS

To create an instance of BoardShim class for your board check required inputs in the table below:

Table 1: Required inputs

Board	Board Id	BrainFlow Board	BrainFlow Board	BrainFlow Board	BrainFlow Board	BrainFlow Board	BrainFlow Board	BrainFlow Board	BrainFlow Board	BrainFlow Board
Playback Board	BoardIds (-3)	PLAYBACK_FILE_BOARD	.	.	Board Id of master board	.	.	path to file for play-back		
Streaming Board	BoardIds (-2)	STREAMING_BOARD	Multicast port IP address	.	Board Id of master board	.	.	.		
Synthetic Board	BoardIds (-1)	SYNTHETIC_BOARD
Cyton	BoardIds (0)	CYTON_BOARD	Cyton serial port(COM3, /dev/ttyUSB0, /dev/cu.usbserial-xxxxxx...)
Ganglion	BoardIds (1)	GANGLION_BOARD	Ganglion serial port(COM3, /dev/ttyUSB0, /dev/cu.usbserial-xxxxxx...)	Timeout for device discovery(default 15sec)	.	.
Cyton Daisy	BoardIds (2)	CYTON_DAISY_BOARD	Cyton serial port(COM3, /dev/ttyUSB0, /dev/cu.usbserial-xxxxxx...)
Ganglion WIFI	BoardIds (4)	GANGLION_WIFI_BOARD	any local IP(default 192.168.4.1) which is free	Timeout for HTTP re-sponse(default 10sec)	.	.
Cyton WIFI	BoardIds (5)	CYTON_WIFI_BOARD	any local IP(default 192.168.4.1) which is free	Timeout for HTTP re-sponse(default 10sec)	.	.
Cyton Daisy WIFI	BoardIds (6)	CYTON_DAISY_WIFI_BOARD	any local IP(default 192.168.4.1) which is free	Timeout for HTTP re-sponse(default 10sec)	.	.
BrainBit	BoardIds (7)	BRAINBIT_BOARD	Timeout for device discovery(default 15sec)	Optional: Serial Number of BrainBit device	.
Unicorn	BoardIds (8)	UNICORN_BOARD	Optional: Serial Num	.

1.1 Playback File Board

This board playbacks file recorded using another BrainFlow board.

It can be extremely useful during development.

To choose this board in BoardShim constructor please specify:

- board_id: -3
- other_info field of BrainFlowInputParams structure, write there board_id for a board which acts like data provider(master board)
- file field of BrainFlowInputParams structure

Supported platforms:

- Windows >= 8.1
- Linux
- MacOS

In methods like:

```
get_eeg_channels (board_id)
get_emg_channels (board_id)
get_ecg_channels (board_id)
# .....
```

You need to use master board id instead Playback Board Id, because exact data format for streaming board is controlled by master board as well as sampling rate

Board Specs:

- num eeg(emg,...) channels: like in master board
- num acceleration channels: like in master board
- sampling rate: like in master board
- communication: UDP multicast socket to read data from master board

1.2 Streaming Board

BrainFlow's boards can stream data to different destinations like file, socket and so on. This board acts like a consumer for data streamed from the main process.

To use it in the first process you should call:

```
# choose any valid multicast address(from "224.0.0.0" to "239.255.255.255") and port
start_stream (450000, 'streaming_board://225.1.1.1:6677')
```

In the second process please specify:

- board_id: -2
- ip_address field of BrainFlowInputParams structure, for example above it's 225.1.1.1
- ip_port field of BrainFlowInputParams structure, for example above it's 6677
- other_info field of BrainFlowInputParams structure, write there board_id for a board which acts like data provider(master board)

Supported platforms:

- Windows >= 8.1
- Linux
- MacOS

In methods like:

```
get_eeg_channels (board_id)
get_emg_channels (board_id)
get_ecg_channels (board_id)
# .....
```

You need to use master board id instead Streaming Board Id, because exact data format for streaming board is controlled by master board as well as sampling rate.

Board Specs:

- num eeg(emg,...) channels: like in master board
- num acceleration channels: like in master board
- sampling rate: like in master board
- communication: UDP multicast socket to read data from master board

1.3 Synthetic Board

This board generates synthetic data and you dont need real hardware to use it.

It can be extremely useful during development.

To choose this board in BoardShim constructor please specify:

- board_id: -1
- you dont need to set any fields in BrainFlowInputParams structure

Supported platforms:


- Windows >= 8.1
- Linux
- MacOS
- Android

Board Specs:

- num eeg(emg,...) channels: 8
- num acceleration channels: 3
- sampling rate: 256
- communication: None

1.4 OpenBCI

1.4.1 Cyton



[Cyton Getting Started Guide from OpenBCI](#)

To choose this board in BoardShim constructor please specify:

- board_id: 0
- serial_port field of BrainFlowInputParams structure

Supported platforms:


- Windows >= 8.1
- Linux
- MacOS

On MacOS there are two serial ports for each device: /dev/tty.... and /dev/cu.... You HAVE to specify /dev/cu....

Board Spec:

- num eeg(emg,...) channels: 8
- num acceleration channels: 3
- sampling rate: 250
- communication: serial port
- signal gain: 24

1.4.2 Ganglion



[Ganglion Getting Started Guide from OpenBCI](#)

To use Ganglion board you need a [dongle](#)

To choose this board in BoardShim constructor please specify:

- board_id: 1
- serial_port field of BrainFlowInputParams structure
- mac_address field of BrainFlowInputParams structure, if its empty BrainFlow will try to autodiscover Ganglion
- optional: timeout field of BrainFlowInputParams structure, default is 15sec

To get Ganglion's MAC address you can use:

- Windows: [Bluetooth LE Explorer App](#)
- Linux: hcitool command

Supported platforms:


- Windows >= 8.1
- Linux
- MacOS

On MacOS there are two serial ports for each device: /dev/tty.... and /dev/cu.... You HAVE to specify /dev/cu....

Board Spec:

- num eeg(emg,...) channels: 4
- num acceleration channels: 3
- sampling rate: 200
- communication: Bluetooth Low Energy behind serial port from the dongle

1.4.3 Cyton Daisy



CytonDaisy Getting Started Guide from OpenBCI

To choose this board in BoardShim constructor please specify:

- board_id: 2
- serial_port field of BrainFlowInputParams structure

Supported platforms:


- Windows >= 8.1
- Linux
- MacOS

On MacOS there are two serial ports for each device: /dev/tty.... and /dev/cu.... You HAVE to specify /dev/cu....

Board Spec:

- num eeg(emg,...) channels: 16
- num acceleration channels: 3
- sampling rate: 125
- communication: serial port
- signal gain: 24

1.4.4 Ganglion with WIFI Shield



[WiFi Shield Getting Started Guide from OpenBCI](#)

[WiFi Shield Programming Guide from OpenBCI](#)

To choose this board in BoardShim constructor please specify:

- board_id: 4
- ip_address field of BrainFlowInputParams structure should contain WiFi Shield Ip address(in direct mode its 192.168.4.1), if it's empty BrainFlow will try to autodiscover WiFi Shield and in case of failure will try to use 192.168.4.1
- ip_port field of BrainFlowInputParams structure should be any local port which is free right now
- optional: timeout field of BrainFlowInputParams structure, default is 10sec


Supported platforms:

- Windows >= 8.1
- Linux
- MacOS
- Android

Board Spec:

- num eeg(emg,...) channels: 4
- num acceleration channels: 3
- sampling rate: 1600
- communication: TCP socket to read data and HTTP to send commands

1.4.5 Cyton with WIFI Shield



WIFI shield Getting Started Guide from OpenBCI

WIFI shield Programming Guide from OpenBCI

To choose this board in BoardShim constructor please specify:

- board_id: 5
- ip_address field of BrainFlowInputParams structure should contain WiFi Shield Ip address(in direct mode its 192.168.4.1), if it's empty BrainFlow will try to autodiscover WIFI Shield and in case of failure will try to use 192.168.4.1
- ip_port field of BrainFlowInputParams structure should be any local port which is free right now
- optional: timeout field of BrainFlowInputParams structure, default is 10sec


Supported platforms:

- Windows >= 8.1
- Linux
- MacOS
- Android

Board Spec:

- num eeg(emg,...) channels: 8
- num acceleration channels: 3
- sampling rate: 1000
- communication: TCP socket to read data and HTTP to send commands
- signal gain: 24

1.4.6 CytonDaisy with WIFI Shield



[WiFi Shield Getting Started Guide from OpenBCI](#)

[WiFi Shield Programming Guide from OpenBCI](#)

To choose this board in BoardShim constructor please specify:

- board_id: 6
- ip_address field of BrainFlowInputParams structure should contain WiFi Shield Ip address(in direct mode its 192.168.4.1), if it's empty BrainFlow will try to autodiscover WiFi Shield and in case of failure will try to use 192.168.4.1
- ip_port field of BrainFlowInputParams structure should be any local port which is free right now
- optional: timeout field of BrainFlowInputParams structure, default is 10sec

Supported platforms:

- Windows >= 8.1
- Linux
- MacOS
- Android

Board Spec:

- num eeg(emg,...) channels: 16
- num acceleration channels: 3
- sampling rate: 1000

- communication: TCP socket to read data and HTTP to send commands
- signal gain: 24

1.5 NeuroMD

1.5.1 BrainBit



BrainBit website

To choose this board in BoardShim constructor please specify:

- board_id: 7
- optional: serial_number field of BrainFlowInputParams structure should contain Serial Number of BrainBit device, use it if you have multiple devices
- optional: timeout field of BrainFlowInputParams structure, default is 15sec

Supported platforms:

- Windows \geq 10
- MacOS

Board Spec:

- num eeg channels: 4
- num acceleration channels: None
- sampling rate: 250

- communication: Bluetooth Low Energy

1.5.2 Callibri(Yellow)



[Callibri website](#)

Callibri can be used to record EMG, ECG and EEG, but based on signal type you need to apply different settings for device.

BrainFlow does it for you, so there are:

- CALLIBRI_EEG_BOARD (board_id 9)
- CALLIBRI_EMG_BOARD (board_id 10)
- CALLIBRI_ECG_BOARD (board_id 11)

To choose this board in BoardShim constructor please specify:

- board_id: 9, 10 or 11 based on data type
- optional: to use electrodes connected vis USB write “ExternalSwitchInputMioUSB” to other_info field of BrainFlowInputParams structure
- optional: timeout field of BrainFlowInputParams structure, default is 15sec

Supported platforms:

- Windows ≥ 10
- MacOS

Board Spec:

- num exg channels: 1
- num acceleration channels: None
- communication: Bluetooth Low Energy

1.6 G.TEC

1.6.1 Unicorn



[Unicorn website](#)

To choose this board in BoardShim constructor please specify:

- board_id: 8
- optional: serial_number field of BrainFlowInputParams structure should contain Serial Number of BrainBit device, use it if you have multiple devices

Supported platforms:

- Ubuntu 18.04, may work on other Linux OSes, it depends on dynamic library provided by Unicorn
- May also work on Raspberry PI, if you replace libunicorn.so by library provided by Unicorn for Raspberry PI

Board Spec:

- num eeg channels: 8

- num acceleration channels: 3
- sampling rate: 250
- communication: Bluetooth Low Energy

1.7 Neurosity

1.7.1 Notion 1



[Notion website](#)

[Link to Neurosity Tutorial](#)

To choose this board in BoardShim constructor please specify:

- board_id: 13
- optional: Serial Number field of BrainFlowInputParams structure, important if you have multiple devices in the same place

Supported platforms:

- Windows
- Linux
- MacOS

Board Spec:

- num eeg channels: 8
- sampling rate: 250
- communication: UDP BroadCast

1.7.2 Notion 2



[Notion website](#)

[Link to Neurosity Tutorial](#)

To choose this board in BoardShim constructor please specify:

- board_id: 14
- optional: Serial Number field of BrainFlowInputParams structure, important if you have multiple devices in the same place

Supported platforms:

- Windows
- Linux
- MacOS

Board Spec:

- num eeg channels: 8
- sampling rate: 250
- communication: UDP BroadCast

INSTALLATION INSTRUCTIONS

2.1 Python

Please, make sure to use Python 3+. Next, install the latest release from PYPI with the following command in terminal

```
python -m pip install brainflow
```

If you want to install it from source files or build unreleased version from Github, you should compile core module first and run

```
cd python-package  
python -m pip install -e .
```

2.2 C#

For C#, only Windows is currently supported.

You are able to install the latest release from [Nuget](#) or build it yourself:

- Compile BrainFlow's core module
- open Visual Studio Solution
- install required nuget packages
- build it using Visual Studio
- **make sure that unmanaged(C++) libraries exist in search path** - set PATH env variable or copy them to correct folder

2.3 R

R binding is based on [reticulate](#) package and calls Python code, so you need to install Python binding first, make sure that reticulate uses correct virtual environment, after that you will be able to build R package from command line or using R Studio, install it and run samples.

2.4 Java

You are able to download jar files directly from [release page](#)

If you want to install it from source files or build unreleased version from github you should compile core module first and run

```
cd java-package
cd brainflow
mvn package
```

2.5 Matlab

Steps to setup Matlab binding for BrainFlow:

- Compile Core Module, using instructions below
- Open Matlab IDE and open brainflow/matlab-package/brainflow folder there
- Add folders lib and inc to Matlab path
- If you want to run Matlab scripts from folders different than brainflow/matlab-package/brainflow you need to add it to your Matlab path too

2.6 Julia

Steps to setup Julia binding for BrainFlow:

- Compile Core Module, using instructions below
- Set PATH(on Windows) or LD_LIBRARY_PATH(on Unix) env variables to ensure that compiled libraries are in search path
- Install BrainFlow package locally

Example

```
# compile core module first
# set env variable
export LD_LIBRARY_PATH=/home/andreyparfenov/brainflow/installed_linux/lib/:$LD_
↪LIBRARY_PATH
cd julia-package/brainflow
julia
# type ']' to switch to pkg terminal
activate . # activate BrainFlow's env
```


2.7 Compilation of Core Module and C++ Binding

2.7.1 Windows

- Install Cmake>=3.13 you can install it from PYPI via pip
- Install Visual Studio 2017, you can use another version but you will need to change cmake generator in batch files or run cmake commands manually. Also in CI we test only VS2017
- Build it as a cmake project manually or use cmd files from tools directory

Compilation using cmd files

```
python -m pip install cmake==3.13.3
# need to run these files from project dir
.\tools\build_win32.cmd
.\tools\build_win64.cmd
```

2.7.2 Linux

- Install Cmake>=3.13 you can install it from PYPI via pip
- If you wanna distribute compiled Linux libraries you HAVE to build it inside manylinux Docker container
- Build it as a cmake project manually or use bash file from tools directory
- You can use any compiler but for Linux we test only GCC, also we test only 64bit libraries for Linux

Compilation using bash file

```
python -m pip install cmake==3.13.3
# you may need to change line endings using dos2unix or text editor for file below
# need to run this file from project dir
bash ./tools/build_linux.sh
```

2.7.3 MacOS

- Install Cmake>=3.13 you can install it from PYPI via pip
- Build it as a cmake project manually or use bash file from tools directory
- You can use any compiler but for MacOS we test only Clang

Compilation using bash file

```
python -m pip install cmake==3.13.3
# you may need to change line endings using dos2unix or text editor for file below
# need to run this file from project dir
bash ./tools/build_mac.sh
```

2.8 Android

To check supported boards for Android visit [Supported Boards](#)

2.8.1 Installation instructions

- Create Java project in Android Studio, Kotlin is not supported
- Download *jniLibs.zip* from [Release page](#)
- Unpack *jniLibs.zip* and copy it's content to *project/app/src/main/jniLibs*
- Download *brainflow-jar-with-dependencies.jar* from [Release page](#) or from [Github package](#)
- Copy *brainflow-jar-with-dependencies.jar* to *project/app/libs* folder

Now you can use BrainFlow SDK in your Android application!

Note: Android Studio inline compiler may show red errors but it should be compiled fine with Gradle. To fix inline compiler you can use *File > Sync Project with Gradle Files* or click at *File > Invalidate Cache/Restart > Invalidate and Restart*

For some API calls you need to provide additional permissions via manifest file of your application

```
<uses-permission android:name="android.permission.INTERNET"></uses-permission>
<uses-permission android:name="android.permission.ACCESS_NETWORK_STATE"></uses-
↳permission>
<uses-permission android:name="android.permission.READ_EXTERNAL_STORAGE"></uses-
↳permission>
<uses-permission android:name="android.permission.WRITE_EXTERNAL_STORAGE"></uses-
↳permission>
```

2.8.2 Compilation using Android NDK

For BrainFlow developers

To test your changes in BrainFlow on Android you need to build it using Android NDK manually.

Compilation instructions:

- [Download Android NDK](#)
- [Download Ninja](#), for Windows there is exe file in tools folder, make sure that *ninja.exe* in search path
- You can also try *MinGW Makefiles* instead Ninja, but it's not tested and may not work
- Build C++ code using *cmake* and Ninja for **all ABIs**
- Compiled libraries will be in *tools/jniLibs* folder

Command line examples

```
# to prepare project
# for arm64-v8a
cmake -G Ninja -DCMAKE_TOOLCHAIN_FILE=E:\android-ndk-r21d-windows-x86_64\android-ndk-
↳r21d\build\cmake\android.toolchain.cmake -DANDROID_NATIVE_API_LEVEL=android-19 -
↳DANDROID_ABI=arm64-v8a ..
# for armeabi-v7a
cmake -G Ninja -DCMAKE_TOOLCHAIN_FILE=E:\android-ndk-r21d-windows-x86_64\android-ndk-
↳r21d\build\cmake\android.toolchain.cmake -DANDROID_NATIVE_API_LEVEL=android-19 -
↳DANDROID_ABI=armeabi-v7a ..
```

(continues on next page)

(continued from previous page)

```
# for x86_64
cmake -G Ninja -DCMAKE_TOOLCHAIN_FILE=E:\android-ndk-r21d-windows-x86_64\android-ndk-
↳r21d\build\cmake\android.toolchain.cmake -DANDROID_NATIVE_API_LEVEL=android-19 -
↳DANDROID_ABI=x86_64 ..
# for x86
cmake -G Ninja -DCMAKE_TOOLCHAIN_FILE=E:\android-ndk-r21d-windows-x86_64\android-ndk-
↳r21d\build\cmake\android.toolchain.cmake -DANDROID_NATIVE_API_LEVEL=android-19 -
↳DANDROID_ABI=x86 ..

# to build
cmake --build . --target install --config Release -j 2 --parallel 2
```


USER API

BrainFlow User API has three main modules:

- BoardShim to read data from a board, it calls methods from underlying BoardController library
- DataFilter to perform signal processing, it calls methods from underlying DataHandler library
- MLModel to calculate derivative metrics, it calls methods from underlying MLModule library

These classes are independent, so if you want, you can use BrainFlow API only for data streaming and perform signal processing by yourself and vice versa.

BrainFlow data acquisition API is board agnostic, so **to select a specific board you need to pass BrainFlow's board id to BoardShim's constructor and an instance of BrainFlowInputParams structure** which should hold information for your specific board, check [Supported Boards](#). for details. This abstraction allows you to switch boards without any changes in code.

In BoardShim, all board data is returned as a 2d array. Rows in this array may contain timestamps, EEG and EMG data and so on. To see instructions how to query specific kind of data check [Data Format Description](#) and [Code Samples](#).

3.1 Python API Reference

3.1.1 brainflow.board_shim

```
class brainflow.board_shim.BoardIds (value)
```

```
    Bases: enum.Enum
```

```
    Enum to store all supported Board Ids
```

```
    PLAYBACK_FILE_BOARD = -3
```

```
    STREAMING_BOARD = -2
```

```
    SYNTHETIC_BOARD = -1
```

```
    CYTON_BOARD = 0
```

```
    GANGLION_BOARD = 1
```

```
    CYTON_DAISY_BOARD = 2
```

```
    GALEA_BOARD = 3
```

```
    GANGLION_WIFI_BOARD = 4
```

```
    CYTON_WIFI_BOARD = 5
```

```
    CYTON_DAISY_WIFI_BOARD = 6
```

```
BRAINBIT_BOARD = 7
UNICORN_BOARD = 8
CALLIBRI_EEG_BOARD = 9
CALLIBRI_EMG_BOARD = 10
CALLIBRI_ECG_BOARD = 11
FASCIA_BOARD = 12
NOTION_OSC_BOARD = 13
NOTION_1_BOARD = 13
NOTION_2_BOARD = 14
IRONBCI_BOARD = 15
FREEEEG32_BOARD = 17
```

```
class brainflow.board_shim.LogLevels(value)
```

Bases: enum.Enum

Enum to store all log levels supported by BrainFlow

```
LEVEL_TRACE = 0
LEVEL_DEBUG = 1
LEVEL_INFO = 2
LEVEL_WARN = 3
LEVEL_ERROR = 4
LEVEL_CRITICAL = 5
LEVEL_OFF = 6
```

```
class brainflow.board_shim.IpProtocolType(value)
```

Bases: enum.Enum

Enum to store Ip Protocol types

```
NONE = 0
UDP = 1
TCP = 2
```

```
class brainflow.board_shim.BrainFlowInputParams
```

Bases: object

inputs parameters for prepare_session method

Parameters

- **serial_port** (*str*) – serial port name is used for boards which reads data from serial port
- **mac_address** (*str*) – mac address for example its used for bluetooth based boards
- **ip_address** (*str*) – ip address is used for boards which reads data from socket connection
- **ip_port** (*int*) – ip port for socket connection, for some boards where we know it in front you dont need this parameter

- **ip_protocol** (*int*) – ip protocol type from IpProtocolType enum
- **other_info** (*str*) – other info
- **serial_number** (*str*) – serial number
- **file** (*str*) – file

exception `brainflow.board_shim.BrainFlowError` (*message: str, exit_code: int*)

Bases: `Exception`

This exception is raised if non-zero exit code is returned from C code

Parameters

- **message** (*str*) – exception message
- **exit_code** (*int*) – exit code flow low level API

class `brainflow.board_shim.BoardShim` (*board_id: int, input_params: brainflow.board_shim.BrainFlowInputParams*)

Bases: `object`

BoardShim class is a primary interface to all boards

Parameters

- **board_id** (*int*) – Id of your board
- **input_params** (*BrainFlowInputParams*) – board specific structure to pass required arguments

classmethod `set_log_level` (*log_level: int*) → `None`

set BrainFlow log level, use it only if you want to write your own messages to BrainFlow logger, otherwise use `enable_board_logger`, `enable_dev_board_logger` or `disable_board_logger`

Parameters **log_level** (*int*) – log level, to specify it you should use values from `LogLevels` enum

classmethod `enable_board_logger` () → `None`

enable BrainFlow Logger with level INFO, uses stderr for log messages by default

classmethod `disable_board_logger` () → `None`

disable BrainFlow Logger

classmethod `enable_dev_board_logger` () → `None`

enable BrainFlow Logger with level TRACE, uses stderr for log messages by default

classmethod `log_message` (*log_level: int, message: str*) → `None`

write your own log message to BrainFlow logger, use it if you wanna have single logger for your own code and BrainFlow's code

Parameters

- **log_level** – log level
- **message** (*str*) – message

classmethod `set_log_file` (*log_file: str*) → `None`

redirect logger from stderr to file, can be called any time

Parameters **log_file** (*str*) – log file name

classmethod `get_sampling_rate` (*board_id: int*) → `int`

get sampling rate for a board

Parameters **board_id** (*int*) – Board Id

Returns sampling rate for this board id

Return type int

Raises **BrainFlowError** – If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

classmethod **get_package_num_channel** (*board_id: int*) → int
get package num channel for a board

Parameters **board_id** (*int*) – Board Id

Returns number of package num channel

Return type int

Raises **BrainFlowError** – If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

classmethod **get_battery_channel** (*board_id: int*) → int
get battery channel for a board

Parameters **board_id** (*int*) – Board Id

Returns number of batter channel

Return type int

Raises **BrainFlowError** – If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

classmethod **get_num_rows** (*board_id: int*) → int
get number of rows in resulting data table for a board

Parameters **board_id** (*int*) – Board Id

Returns number of rows in returned numpy array

Return type int

Raises **BrainFlowError** – If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

classmethod **get_timestamp_channel** (*board_id: int*) → int
get timestamp channel in resulting data table for a board

Parameters **board_id** (*int*) – Board Id

Returns number of timestamp channel in returned numpy array

Return type int

Raises **BrainFlowError** – If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

classmethod **get_eeg_names** (*board_id: int*) → List[str]
get names of EEG channels in 10-20 system if their location is fixed

Parameters **board_id** (*int*) – Board Id

Returns EEG channels names

Return type List[str]

Raises **BrainFlowError** – If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

classmethod `get_eeg_channels` (*board_id: int*) → List[int]
 get list of eeg channels in resulting data table for a board

Parameters `board_id` (*int*) – Board Id

Returns list of eeg channels in returned numpy array

Return type List[int]

Raises **BrainFlowError** – If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

classmethod `get_exg_channels` (*board_id: int*) → List[int]
 get list of exg channels in resulting data table for a board

Parameters `board_id` (*int*) – Board Id

Returns list of eeg channels in returned numpy array

Return type List[int]

Raises **BrainFlowError** – If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

classmethod `get_emg_channels` (*board_id: int*) → List[int]
 get list of emg channels in resulting data table for a board

Parameters `board_id` (*int*) – Board Id

Returns list of eeg channels in returned numpy array

Return type List[int]

Raises **BrainFlowError** – If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

classmethod `get_ecg_channels` (*board_id: int*) → List[int]
 get list of ecg channels in resulting data table for a board

Parameters `board_id` (*int*) – Board Id

Returns list of ecg channels in returned numpy array

Return type List[int]

Raises **BrainFlowError** – If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

classmethod `get_eog_channels` (*board_id: int*) → List[int]
 get list of eog channels in resulting data table for a board

Parameters `board_id` (*int*) – Board Id

Returns list of eog channels in returned numpy array

Return type List[int]

Raises **BrainFlowError** – If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

classmethod `get_eda_channels` (*board_id: int*) → List[int]
 get list of eda channels in resulting data table for a board

Parameters `board_id` (*int*) – Board Id

Returns list of eda channels in returned numpy array

Return type List[int]

Raises `BrainFlowError` – If this board has no such data exit code is `UNSUPPORTED_BOARD_ERROR`

classmethod `get_ppg_channels` (*board_id: int*) → List[int]

get list of ppg channels in resulting data table for a board

Parameters `board_id` (*int*) – Board Id

Returns list of ppg channels in returned numpy array

Return type List[int]

Raises `BrainFlowError` – If this board has no such data exit code is `UNSUPPORTED_BOARD_ERROR`

classmethod `get_accel_channels` (*board_id: int*) → List[int]

get list of accel channels in resulting data table for a board

Parameters `board_id` (*int*) – Board Id

Returns list of accel channels in returned numpy array

Return type List[int]

Raises `BrainFlowError` – If this board has no such data exit code is `UNSUPPORTED_BOARD_ERROR`

classmethod `get_analog_channels` (*board_id: int*) → List[int]

get list of analog channels in resulting data table for a board

Parameters `board_id` (*int*) – Board Id

Returns list of analog channels in returned numpy array

Return type List[int]

Raises `BrainFlowError` – If this board has no such data exit code is `UNSUPPORTED_BOARD_ERROR`

classmethod `get_gyro_channels` (*board_id: int*) → List[int]

get list of gyro channels in resulting data table for a board

Parameters `board_id` (*int*) – Board Id

Returns list of gyro channels in returned numpy array

Return type List[int]

Raises `BrainFlowError` – If this board has no such data exit code is `UNSUPPORTED_BOARD_ERROR`

classmethod `get_other_channels` (*board_id: int*) → List[int]

get list of other channels in resulting data table for a board

Parameters `board_id` (*int*) – Board Id

Returns list of other channels in returned numpy array

Return type List[int]

Raises `BrainFlowError` – If this board has no such data exit code is `UNSUPPORTED_BOARD_ERROR`

classmethod `get_temperature_channels` (*board_id: int*) → List[int]

get list of temperature channels in resulting data table for a board

Parameters `board_id` (*int*) – Board Id

Returns list of temperature channels in returned numpy array

Return type List[int]

Raises BrainFlowError – If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

classmethod get_resistance_channels (*board_id: int*) → List[int]

get list of resistance channels in resulting data table for a board

Parameters board_id (*int*) – Board Id

Returns list of resistance channels in returned numpy array

Return type List[int]

Raises BrainFlowError – If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

prepare_session () → None

prepare streaming session, init resources, you need to call it before any other BoardShim object methods

start_stream (*num_samples: int = 450000, streamer_params: str = None*) → None

Start streaming data, this methods stores data in ringbuffer

Parameters

- **num_samples** (*int*) – size of ring buffer to keep data
- **parameter to stream data from brainflow, supported vals** (*streamer_params*) – “file://%file_name%:w”, “file://%file_name%:a”, “streaming_board://%multicast_group_ip%:%port%”. Range for multicast addresses is from “224.0.0.0” to “239.255.255.255”

stop_stream () → None

Stop streaming data

release_session () → None

release all resources

get_current_board_data ()

Get specified amount of data or less if there is not enough data, doesnt remove data from ringbuffer

Parameters num_samples (*int*) – max number of samples

Returns latest data from a board

Return type NDArray[Float64]

get_board_data_count () → int

Get num of elements in ringbuffer

Returns number of elements in ring buffer

Return type int

get_board_id () → int

Get's the actual board id, can be different than provided

Returns board id

Return type int

is_prepared () → bool

Check if session is ready or not

Returns session status

Return type bool

get_board_data ()

Get all board data and remove them from ringbuffer

Returns all data from a board

Return type Numpy[Float64]

config_board (*config*) → None

Use this method carefully and only if you understand what you are doing, do NOT use it to start or stop streaming

Parameters **config** (*str*) – string to send to a board

Returns response string if any

Return type str

3.1.2 brainflow.exit_codes

class brainflow.exit_codes.BrainflowExitCodes (*value*)

Bases: enum.Enum

Enum to store all possible exit codes

STATUS_OK = 0

PORT_ALREADY_OPEN_ERROR = 1

UNABLE_TO_OPEN_PORT_ERROR = 2

SER_PORT_ERROR = 3

BOARD_WRITE_ERROR = 4

INCOMING_MSG_ERROR = 5

INITIAL_MSG_ERROR = 6

BOARD_NOT_READY_ERROR = 7

STREAM_ALREADY_RUN_ERROR = 8

INVALID_BUFFER_SIZE_ERROR = 9

STREAM_THREAD_ERROR = 10

STREAM_THREAD_IS_NOT_RUNNING = 11

EMPTY_BUFFER_ERROR = 12

INVALID_ARGUMENTS_ERROR = 13

UNSUPPORTED_BOARD_ERROR = 14

BOARD_NOT_CREATED_ERROR = 15

ANOTHER_BOARD_IS_CREATED_ERROR = 16

GENERAL_ERROR = 17

SYNC_TIMEOUT_ERROR = 18

JSON_NOT_FOUND_ERROR = 19

NO_SUCH_DATA_IN_JSON_ERROR = 20

```

CLASSIFIER_IS_NOT_PREPARED_ERROR = 21
ANOTHER_CLASSIFIER_IS_PREPARED_ERROR = 22
UNSUPPORTED_CLASSIFIER_AND_METRIC_COMBINATION_ERROR = 23

```

3.1.3 brainflow.data_filter

```

class brainflow.data_filter.FilterTypes(value)
    Bases: enum.Enum

    Enum to store all supported Filter Types

    BUTTERWORTH = 0
    CHEBYSHEV_TYPE_1 = 1
    BESSEL = 2

class brainflow.data_filter.AggOperations(value)
    Bases: enum.Enum

    Enum to store all supported aggregation operations

    MEAN = 0
    MEDIAN = 1
    EACH = 2

class brainflow.data_filter.WindowFunctions(value)
    Bases: enum.Enum

    Enum to store all supported window functions

    NO_WINDOW = 0
    HANNING = 1
    HAMMING = 2
    BLACKMAN_HARRIS = 3

class brainflow.data_filter.DetrendOperations(value)
    Bases: enum.Enum

    Enum to store all supported detrend options

    NONE = 0
    CONSTANT = 1
    LINEAR = 2

class brainflow.data_filter.DataFilter
    Bases: object

    DataFilter class contains methods for signal processig

    classmethod enable_data_logger() → None
        enable Data Logger with level INFO, uses stderr for log messages by default

    classmethod disable_data_logger() → None
        disable Data Logger

    classmethod enable_dev_data_logger() → None
        enable Data Logger with level TRACE, uses stderr for log messages by default

```

classmethod **set_log_file** (*log_file: str*) → None
redirect logger from stderr to file, can be called any time

Parameters **log_file** (*str*) – log file name

classmethod **perform_lowpass** ()
apply low pass filter to provided data

Parameters

- **data** (*NDArray [Float64]*) – data to filter, filter works in-place
- **sampling_rate** (*int*) – board's sampling rate
- **cutoff** (*float*) – cutoff frequency
- **order** (*int*) – filter order
- **filter_type** (*int*) – filter type from special enum
- **ripple** (*float*) – ripple value for Chebyshev filter

classmethod **perform_highpass** ()
apply high pass filter to provided data

Parameters

- **data** (*NDArray [Float64]*) – data to filter, filter works in-place
- **sampling_rate** (*int*) – board's sampling rate
- **cutoff** (*float*) – cutoff frequency
- **order** (*int*) – filter order
- **filter_type** (*int*) – filter type from special enum
- **ripple** (*float*) – ripple value for Chebyshev filter

classmethod **perform_bandpass** ()
apply band pass filter to provided data

Parameters

- **data** (*NDArray [Float64]*) – data to filter, filter works in-place
- **sampling_rate** (*int*) – board's sampling rate
- **center_freq** (*float*) – center frequency
- **band_width** (*float*) – band width
- **order** (*int*) – filter order
- **filter_type** (*int*) – filter type from special enum
- **ripple** (*float*) – ripple value for Chebyshev filter

classmethod **perform_bandstop** ()
apply band stop filter to provided data

Parameters

- **data** (*NDArray [Float64]*) – data to filter, filter works in-place
- **sampling_rate** (*int*) – board's sampling rate
- **center_freq** (*float*) – center frequency
- **band_width** (*float*) – band width

- **order** (*int*) – filter order
- **filter_type** (*int*) – filter type from special enum
- **ripple** (*float*) – ripple value for Chebyshev filter

classmethod perform_rolling_filter()

smooth data using moving average or median

Parameters

- **data** (*NDArray[Float64]*) – data to smooth, it works in-place
- **period** (*int*) – window size
- **operation** (*int*) – int value from AggOperation enum

classmethod perform_downsampling()

perform data downsampling, it doesnt apply lowpass filter for you, it just aggregates several data points

Parameters

- **data** (*NDArray[Float64]*) – initial data
- **period** (*int*) – downsampling period
- **operation** (*int*) – int value from AggOperation enum

Returns downsampled data

Return type *NDArray[Float64]*

classmethod perform_wavelet_transform()

perform wavelet transform

Parameters

- **data** (*NDArray[Float64]*) – initial data
- **wavelet** (*str*) – supported vals: db1..db15,haar,sym2..sym10,coif1..coif5,bior1.1,bior1.3,bior1.5,bior2.2,bior2.4,bior3.7,bior3.9,bior4.4,bior5.5,bior6.8
- **decomposition_level** (*int*) – level of decomposition

Returns tuple of wavelet coeffs in format [A(J) D(J) D(J-1) ... D(1)] where J is decomposition level, A - app coeffs, D - detailed coeffs, and array with lengths for each block

Return type tuple

classmethod perform_inverse_wavelet_transform()

perform wavelet transform

Parameters

- **wavelet_output** – tuple of wavelet_coeffs and array with lengths
- **original_data_len** (*int*) – len of signal before wavelet transform
- **wavelet** (*str*) – supported vals: db1..db15,haar,sym2..sym10,coif1..coif5,bior1.1,bior1.3,bior1.5,bior2.2,bior2.4,bior3.7,bior3.9,bior4.4,bior5.5,bior6.8
- **decomposition_level** (*int*) – level of decomposition

Returns restored data

Return type *NDArray[Float64]*

classmethod perform_wavelet_denoising()

perform wavelet denoising

Parameters

- **data** (*NDArray [Float64]*) – data to denoise
- **wavelet** (*str*) – supported vals: db1..db15,haar,sym2..sym10,coif1..coif5,bior1.1,bior1.3,bior1.5,bior2.2,bior2.4,bior3.7,bior3.9,bior4.4,bior5.5,bior6.8
- **decomposition_level** (*int*) – decomposition level

classmethod perform_fft()

perform direct fft

Parameters

- **data** (*NDArray [Float64]*) – data for fft, len of data must be a power of 2
- **window** (*int*) – window function

Returns numpy array of complex values, len of this array is $N / 2 + 1$ **Return type** *NDArray[Complex128]***classmethod get_psd()**

calculate PSD

Parameters

- **data** (*NDArray [Float64]*) – data to calc psd, len of data must be a power of 2
- **sampling_rate** (*int*) – sampling rate
- **window** (*int*) – window function

Returns amplitude and frequency arrays of len $N / 2 + 1$ **Return type** tuple**classmethod get_psd_welch()**

calculate PSD using Welch method

Parameters

- **data** (*NDArray [Float64]*) – data to calc psd
- **nfft** (*int*) – FFT Window size, must be power of 2
- **overlap** (*int*) – overlap of FFT Windows, must be between 0 and nfft
- **sampling_rate** (*int*) – sampling rate
- **window** (*int*) – window function

Returns amplitude and frequency arrays of len $N / 2 + 1$ **Return type** tuple**classmethod detrend()**

detrend data

Parameters

- **data** (*NDArray [Float64]*) – data to calc psd
- **detrend_operation** (*int*) – Type of detrend operation

classmethod get_band_power (*psd: Tuple, freq_start: float, freq_end: float*) → float

calculate band power

Parameters

- **psd** (*tuple*) – psd from get_psd
- **freq_start** (*int*) – start freq
- **freq_end** (*int*) – end freq

Returns band power

Return type float

classmethod get_avg_band_powers (*data: nptyping.types._ndarray.NDArray, channels: List, sampling_rate: int, apply_filter: bool*) → Tuple
calculate avg and stddev of BandPowers across all channels

Parameters

- **data** (*NDArray*) – 2d array for calculation
- **channels** (*List*) – channels - rows of data array which should be used for calculation
- **sampling_rate** (*int*) – sampling rate
- **apply_filter** (*bool*) – apply bandpass and bandstop filters or not

Returns avg and stddev arrays for bandpowers

Return type tuple

classmethod perform_ifft ()
perform inverse fft

Parameters **data** (*NDArray[Complex128]*) – data from fft

Returns restored data

Return type NDArray[Float64]

classmethod get_nearest_power_of_two (*value: int*) → int
calc nearest power of two

Parameters **value** (*int*) – input value

Returns nearest power of two

Return type int

classmethod write_file (*data, file_name: str, file_mode: str*) → None
write data to file, in file data will be transposed

Parameters

- **data** (*2d numpy array*) – data to store in a file
- **file_name** (*str*) – file name to store data
- **file_mode** (*str*) – 'w' to rewrite file or 'a' to append data to file

classmethod read_file (*file_name: str*)
read data from file

Parameters **file_name** (*str*) – file name to read

Returns 2d numpy array with data from this file, data will be transposed to original dimensions

Return type 2d numpy array

3.1.4 brainflow.ml_model

```
class brainflow.ml_model.BrainFlowMetrics (value)
    Bases: enum.Enum

    Enum to store all supported metrics

    RELAXATION = 0

    CONCENTRATION = 1

class brainflow.ml_model.BrainFlowClassifiers (value)
    Bases: enum.Enum

    Enum to store all supported classifiers

    REGRESSION = 0

    KNN = 1

    SVM = 2

    LDA = 3

class brainflow.ml_model.BrainFlowModelParams (metric, classifier)
    Bases: object

    inputs parameters for prepare_session method

    Parameters
        • metric (int) – metric to calculate
        • classifier (int) – classifier to use
        • file (str) – file to load model
        • other_info (int) – additional information

class brainflow.ml_model.MLModel (model_params: brainflow.ml_model.BrainFlowModelParams)
    Bases: object

    MLModel class used to calc derivative metrics from raw data

    Parameters model_params (BrainFlowModelParams) – Model Params

    classmethod enable_ml_logger () → None
        enable ML Logger with level INFO, uses stderr for log messages by default

    classmethod disable_ml_logger () → None
        disable BrainFlow Logger

    classmethod enable_dev_ml_logger () → None
        enable ML Logger with level TRACE, uses stderr for log messages by default

    classmethod set_log_file (log_file: str) → None
        redirect logger from stderr to file, can be called any time

    Parameters log_file (str) – log file name

    prepare () → None
        prepare classifier

    release () → None
        release classifier
```

predict (*data: nptyping.types._ndarray.NDArray*) → float
calculate metric from data

Parameters *data* (*NDArray*) – input array

Returns metric value

Return type float

3.2 C++ API Reference

3.2.1 BoardShim class

class BoardShim

BoardShim class to communicate with a board.

Public Functions

BoardShim (int *board_id*, **struct** BrainFlowInputParams *params*)

~BoardShim ()

void **prepare_session** ()

prepare BrainFlow's streaming session, should be called first

void **start_stream** (int *buffer_size* = 450000, char **streamer_params* = NULL)

start streaming thread and store data in ringbuffer

Parameters

- *buffer_size*: size of internal ring buffer
- *streamer_params*: use it to pass data packages further or store them directly during streaming, supported values: "file://%file_name%.w", "file://%file_name%.a", "streaming_board://%multicast_group_ip%:%port%".

Range for multicast addresses is from "224.0.0.0" to "239.255.255.255"

bool **is_prepared** ()

check if session is ready or not

void **stop_stream** ()

stop streaming thread, doesnt release other resources

void **release_session** ()

release streaming session

double ****get_current_board_data** (int *num_samples*, int **num_data_points*)

get latest collected data, doesnt remove it from ringbuffer

int **get_board_id** ()

Get board id, for some boards can be different than provided (playback, streaming)

int **get_board_data_count** ()

get number of packages in ringbuffer

double ****get_board_data** (int **num_data_points*)

get all collected data and flush it from internal buffer

std::string **config_board** (char **config*)
send string to a board, use it carefully and only if you understand what you are doing

Public Members

int **board_id**

Public Static Functions

void **disable_board_logger** ()
disable BrainFlow loggers

void **enable_board_logger** ()
enable BrainFlow logger with LEVEL_INFO

void **enable_dev_board_logger** ()
enable BrainFlow logger with LEVEL_TRACE

void **set_log_file** (char **log_file*)
redirect BrainFlow logger from stderr to file

void **set_log_level** (int *log_level*)
use set_log_level only if you want to write your own log messages to BrainFlow logger

void **log_message** (int *log_level*, const char **format*, ...)
write user defined string to BrainFlow logger

int **get_sampling_rate** (int *board_id*)
get sampling rate for this board

Parameters

- *board_id*: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int **get_package_num_channel** (int *board_id*)
get row index which holds package nums

Parameters

- *board_id*: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int **get_timestamp_channel** (int *board_id*)
get row index which holds timestamps

Parameters

- *board_id*: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int **get_battery_channel** (int *board_id*)
 get row index which holds battery level info

Parameters

- *board_id*: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int **get_num_rows** (int *board_id*)
 get number of rows in returned from *get_board_data()* 2d array

Parameters

- *board_id*: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

std::string ***get_eeg_names** (int *board_id*, int **len*)
 get eeg channel names in 10-20 system for devices with fixed electrode locations

Parameters

- *board_id*: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int ***get_eeg_channels** (int *board_id*, int **len*)
 get row indices which hold EEG data, for some board we can not split EEG...

Parameters

- *board_id*: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int ***get_emg_channels** (int *board_id*, int **len*)
 get row indices which hold EMG data, for some board we can not split EEG...

Parameters

- *board_id*: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int ***get_ecg_channels** (int *board_id*, int **len*)
 get row indices which hold ECG data, for some board we can not split EEG...

Parameters

- *board_id*: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int ***get_eog_channels** (int *board_id*, int **len*)
get row indices which hold EOG data, for some board we can not split EEG...

Parameters

- *board_id*: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int ***get_exg_channels** (int *board_id*, int **len*)
get row indices which hold EXG data

Parameters

- *board_id*: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int ***get_ppg_channels** (int *board_id*, int **len*)
get row indices which hold PPG data

Parameters

- *board_id*: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int ***get_eda_channels** (int *board_id*, int **len*)
get row indices which hold EDA data

Parameters

- *board_id*: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int ***get_accel_channels** (int *board_id*, int **len*)
get row indices which hold accel data

Parameters

- *board_id*: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int ***get_analog_channels** (int *board_id*, int **len*)
get row indices which hold analog data

Parameters

- `board_id`: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is `UNSUPPORTED_BOARD_ERROR`

int ***get_gyro_channels** (int *board_id*, int **len*)
get row indices which hold gyro data

Parameters

- `board_id`: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is `UNSUPPORTED_BOARD_ERROR`

int ***get_other_channels** (int *board_id*, int **len*)
get row indices which hold other information

Parameters

- `board_id`: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is `UNSUPPORTED_BOARD_ERROR`

int ***get_temperature_channels** (int *board_id*, int **len*)
get row indices which hold temperature data

Parameters

- `board_id`: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is `UNSUPPORTED_BOARD_ERROR`

int ***get_resistance_channels** (int *board_id*, int **len*)
get row indices which hold resistance data

Parameters

- `board_id`: board id of your device

Exceptions

- *BrainFlowException*: If this board has no such data exit code is `UNSUPPORTED_BOARD_ERROR`

3.2.2 DataFilter class

class DataFilter

DataFilter class to perform signal processing.

Public Static Functions

void **enable_data_logger** ()

enable Data logger with LEVEL_INFO

void **disable_data_logger** ()

disable Data loggers

void **enable_dev_data_logger** ()

enable Data logger with LEVEL_TRACE

void **set_log_file** (char **log_file*)

void **perform_lowpass** (double **data*, int *data_len*, int *sampling_rate*, double *cutoff*, int *order*, int *filter_type*, double *ripple*)

perform low pass filter in-place

void **perform_highpass** (double **data*, int *data_len*, int *sampling_rate*, double *cutoff*, int *order*, int *filter_type*, double *ripple*)

perform high pass filter in-place

void **perform_bandpass** (double **data*, int *data_len*, int *sampling_rate*, double *center_freq*, double *band_width*, int *order*, int *filter_type*, double *ripple*)

perform bandpass filter in-place

void **perform_bandstop** (double **data*, int *data_len*, int *sampling_rate*, double *center_freq*, double *band_width*, int *order*, int *filter_type*, double *ripple*)

perform bandstop filter in-place

void **perform_rolling_filter** (double **data*, int *data_len*, int *period*, int *agg_operation*)

perform moving average or moving median filter in-place

double ***perform_downsampling** (double **data*, int *data_len*, int *period*, int *agg_operation*, int **filtered_size*)

perform data downsampling, it just aggregates several data points

std::pair<double*, int*> **perform_wavelet_transform** (double **data*, int *data_len*, char **wavelet*, int *decomposition_level*)

perform wavelet transform

Return std::pair of wavelet coeffs array in format [A(J) D(J) D(J-1) D(1)] where J is decomposition level A - app coeffs, D - detailed coeffs, and array of lengths for each block in wavelet coeffs array, length of this array is decomposition_level + 1

Parameters

- *data*: input array, any size
- *data_len*: length of input array
- *wavelet*: supported vals: db1..db15,haar,sym2..sym10,coif1..coif5,bior1.1,bior1.3,bior1.5,bior2.2,bior2.4,bior2.6,bior3.7,bior3.9,bior4.4,bior5.5,bior6.8
- *decomposition_level*: level of decomposition in wavelet transform

double ***perform_inverse_wavelet_transform** (std::pair<double*, int*> *wavelet_output*, int *original_data_len*, char **wavelet*, int *decomposition_level*)

performs inverse wavelet transform

void **perform_wavelet_denoising** (double **data*, int *data_len*, char **wavelet*, int *decomposition_level*)

perform wavelet denoising

std::complex<double> ***perform_fft** (double **data*, int *data_len*, int *window*)

perform direct fft

Return complex array with size $\text{data_len} / 2 + 1$, it holds only positive im values

Parameters

- *data*: input array
- *data_len*: must be power of 2
- *window*: window function

double ***perform_iftft** (std::complex<double> **data*, int *data_len*)

perform inverse fft

Return restored data

Parameters

- *data*: complex array from `perform_fft`
- *data_len*: len of original array, must be power of 2

int **get_nearest_power_of_two** (int *value*)

calculate nearest power of 2

Return nearest power of 2

Parameters

- *value*: input value

std::pair<double*, double*> **get_psd** (double **data*, int *data_len*, int *sampling_rate*, int *window*)

calculate PSD

Return pair of amplitude and freq arrays of size $\text{data_len} / 2 + 1$

Parameters

- *data*: input array
- *data_len*: must be power of 2
- *sampling_rate*: sampling rate
- *window*: window function

void **detrend** (double **data*, int *data_len*, int *detrend_operation*)

subtract trend from data

Parameters

- *data*: input array
- *data_len*:
- *detrend_operation*: use `DetrendOperations` enum

```
std::pair<double*, double*> get_psd_welch (double *data, int data_len, int nfft, int overlap, int sampling_rate, int window)
```

```
double get_band_power (std::pair<double*, double*> psd, int data_len, double freq_start, double freq_end)
```

calculate band power

Return band power

Parameters

- psd: psd calculated using get_psd
- data_len: len of ampl and freq arrays: $N / 2 + 1$ where N is FFT size
- freq_start: lowest frequency
- freq_end: highest frequency

```
std::pair<double*, double*> get_avg_band_powers (double **data, int cols, int *channels, int channels_len, int sampling_rate, bool apply_filters)
```

calculate avg and stddev of BandPowers across all channels

Return pair of double arrays of size 5, first of them - avg band powers, second stddev

Parameters

- data: input 2d array
- cols: number of cols in 2d array - number of datapoints
- channels: array of rows - eeg channels which should be used
- channels_len: - len of channels array
- sampling_rate: sampling rate
- apply_filters: set to true to apply filters before band power calculations

```
void write_file (double **data, int num_rows, int num_cols, char *file_name, char *file_mode)
```

write file, in file data will be transposed

```
double **read_file (int *num_rows, int *num_cols, char *file_name)
```

read data from file, data will be transposed to original format

3.2.3 BrainFlowException class

```
class BrainFlowException : public exception
```

BrainFlowException class to notify about errors.

Public Functions

```
BrainFlowException (const char *msg, int exit_code)
```

```
~BrainFlowException ()
```

```
const char *what () const
```

Public Members

int **exit_code**
exit code returned from low level API

3.2.4 BrainFlowModelParams class

Warning: doxygenclass: Cannot find class “BrainFlowModelParams” in doxygen xml output for project “BrainFlowCpp” from directory: build-cpp/xml

3.2.5 MLModel class

class MLModel
Calculates different metrics from raw data.

Public Functions

MLModel (**struct** BrainFlowModelParams *params*)

~MLModel ()

void **prepare** ()
initialize classifier, should be called first

double **predict** (double **data*, int *data_len*)
calculate metric from data

void **release** ()
release classifier

Public Static Functions

void **set_log_file** (char **log_file*)
redirect logger to a file

void **enable_ml_logger** ()
enable ML logger with LEVEL_INFO

void **disable_ml_logger** ()
disable ML loggers

void **enable_dev_ml_logger** ()
enable ML logger with LEVEL_TRACE

3.2.6 BrainFlow constants

```
#pragma once

enum class BrainFlowExitCodes : int
{
    STATUS_OK = 0,
    PORT_ALREADY_OPEN_ERROR = 1,
    UNABLE_TO_OPEN_PORT_ERROR = 2,
    SET_PORT_ERROR = 3,
    BOARD_WRITE_ERROR = 4,
    INCOMING_MSG_ERROR = 5,
    INITIAL_MSG_ERROR = 6,
    BOARD_NOT_READY_ERROR = 7,
    STREAM_ALREADY_RUN_ERROR = 8,
    INVALID_BUFFER_SIZE_ERROR = 9,
    STREAM_THREAD_ERROR = 10,
    STREAM_THREAD_IS_NOT_RUNNING = 11,
    EMPTY_BUFFER_ERROR = 12,
    INVALID_ARGUMENTS_ERROR = 13,
    UNSUPPORTED_BOARD_ERROR = 14,
    BOARD_NOT_CREATED_ERROR = 15,
    ANOTHER_BOARD_IS_CREATED_ERROR = 16,
    GENERAL_ERROR = 17,
    SYNC_TIMEOUT_ERROR = 18,
    JSON_NOT_FOUND_ERROR = 19,
    NO_SUCH_DATA_IN_JSON_ERROR = 20,
    CLASSIFIER_IS_NOT_PREPARED_ERROR = 21,
    ANOTHER_CLASSIFIER_IS_PREPARED_ERROR = 22,
    UNSUPPORTED_CLASSIFIER_AND_METRIC_COMBINATION_ERROR = 23,
};

enum class BoardIds : int
{
    PLAYBACK_FILE_BOARD = -3,
    STREAMING_BOARD = -2,
    SYNTHETIC_BOARD = -1,
    CYTON_BOARD = 0,
    GANGLION_BOARD = 1,
    CYTON_DAISSY_BOARD = 2,
    GALEA_BOARD = 3,
    GANGLION_WIFI_BOARD = 4,
    CYTON_WIFI_BOARD = 5,
    CYTON_DAISSY_WIFI_BOARD = 6,
    BRAINBIT_BOARD = 7,
    UNICORN_BOARD = 8,
    CALLIBRI_EEG_BOARD = 9,
    CALLIBRI_EMG_BOARD = 10,
    CALLIBRI_ECG_BOARD = 11,
    FASCIA_BOARD = 12,
    NOTION_1_BOARD = 13,
    NOTION_2_BOARD = 14,
    IRONBCI_BOARD = 15,
    FREEEEG32_BOARD = 17
};

enum class FilterTypes : int
```

(continues on next page)

(continued from previous page)

```

{
    BUTTERWORTH = 0,
    CHEBYSHEV_TYPE_1 = 1,
    BESSEL = 2
};

enum class AggOperations : int
{
    MEAN = 0,
    MEDIAN = 1,
    EACH = 2
};

enum class WindowFunctions : int
{
    NO_WINDOW = 0,
    HANNING = 1,
    HAMMING = 2,
    BLACKMAN_HARRIS = 3
};

enum class DetrendOperations : int
{
    NONE = 0,
    CONSTANT = 1,
    LINEAR = 2
};

enum class BrainFlowMetrics : int
{
    RELAXATION = 0,
    CONCENTRATION = 1
};

enum class BrainFlowClassifiers : int
{
    REGRESSION = 0,
    KNN = 1,
    SVM = 2,
    LDA = 3
};

/// LogLevels enum to store all possible log levels
enum class LogLevels : int
{
    LEVEL_TRACE = 0,    /// TRACE
    LEVEL_DEBUG = 1,    /// DEBUG
    LEVEL_INFO = 2,     /// INFO
    LEVEL_WARN = 3,     /// WARN
    LEVEL_ERROR = 4,    /// ERROR
    LEVEL_CRITICAL = 5, /// CRITICAL
    LEVEL_OFF = 6       // OFF
};

```

3.3 Java API Reference

Content of Brainflow Package:

enum AggOperations

enum to store all supported aggregation operations

Public Functions

```
int get_code ()
```

```
brainflow::AggOperations (final int code)
```

Public Members

```
brainflow::MEAN    =(0)
```

```
brainflow::MEDIAN  =(1)
```

```
brainflow::EACH    =(2)
```

Public Static Functions

```
String brainflow::string_from_code (final int code)
```

```
AggOperations brainflow::from_code (final int code)
```

```
brainflow::[static initializer]
```

Private Members

```
final int brainflow::agg_operation
```

Private Static Attributes

```
final Map< Integer, AggOperations > brainflow::ao_map    = new HashMap<Integer, AggOperations>()
```

enum BoardIds

enum to store all supported boards

Public Functions

```
int get_code ()
```

```
brainflow::BoardIds (final int code)
```

Public Members

```
brainflow::PLAYBACK_FILE_BOARD    = (-3)
brainflow::STREAMING_BOARD        = (-2)
brainflow::SYNTHETIC_BOARD        = (-1)
brainflow::CYTON_BOARD            = (0)
brainflow::GANGLION_BOARD         = (1)
brainflow::CYTON_DAISSY_BOARD     = (2)
brainflow::GALEA_BOARD            = (3)
brainflow::GANGLION_WIFI_BOARD    = (4)
brainflow::CYTON_WIFI_BOARD       = (5)
brainflow::CYTON_DAISSY_WIFI_BOARD = (6)
brainflow::BRAINBIT_BOARD         = (7)
brainflow::UNICORN_BOARD          = (8)
brainflow::CALLIBRI_EEG_BOARD     = (9)
brainflow::CALLIBRI_EMG_BOARD     = (10)
brainflow::CALLIBRI_ECG_BOARD     = (11)
brainflow::FASCIA_BOARD           = (12)
brainflow::NOTION_1_BOARD         = (13)
brainflow::NOTION_2_BOARD         = (14)
brainflow::IRONBCI_BOARD          = (15)
brainflow::FREEEEEG32_BOARD       = (17)
```

Public Static Functions

```
String brainflow::string_from_code (final int code)
BoardIds brainflow::from_code (final int code)
brainflow::[static initializer]
```

Private Members

```
final int brainflow::board_id
```

Private Static Attributes

```
final Map< Integer, BoardIds > brainflow::bi_map = new HashMap<Integer, BoardIds> ()  
  
class brainflow::brainflow::BoardShim  
    BoardShim class to communicate with a board
```

Public Functions

```
BoardShim (int board_id, BrainFlowInputParams params)  
    Create BoardShim object  
  
void prepare_session ()  
    prepare steaming session, allocate resources  
  
int get_board_id ()  
    Get Board Id, can be different than provided (playback or streaming board)  
  
String config_board (String config)  
    send string to a board, use this method carefully and only if you understand what you are doing  
  
void start_stream (int buffer_size, String streamer_params)  
    start streaming thread, store data in internal ringbuffer and stream them from brainflow at the same time
```

Parameters

- `buffer_size`: size of internal ringbuffer
- `streamer_params`: supported vals: “file://%file_name%:w”, “file://%file_name%:a”, “streaming_board://%multicast_group_ip%:%port%”. Range for multicast addresses is from “224.0.0.0” to “239.255.255.255”

```
void start_stream ()  
    start streaming thread, store data in internal ringbuffer  
  
void start_stream (int buffer_size)  
    start streaming thread, store data in internal ringbuffer  
  
void stop_stream ()  
    stop streaming thread  
  
void release_session ()  
    release all resources  
  
int get_board_data_count ()  
    get number of packages in ringbuffer  
  
boolean is_prepared ()  
    check session status  
  
double [][] get_current_board_data (int num_samples)  
    get latest collected data, can return less than “num_samples”, doesnt flush it from ringbuffer  
  
double [][] get_board_data ()  
    get all data from ringbuffer and flush it
```


Public Members

int **board_id**
BrainFlow's board id

Public Static Functions

void **enable_board_logger** ()
enable BrainFlow logger with level INFO

void **enable_dev_board_logger** ()
enable BrainFlow logger with level TRACE

void **disable_board_logger** ()
disable BrainFlow logger

void **set_log_file** (String *log_file*)
redirect logger from stderr to a file

void **set_log_level** (int *log_level*)
set log level

void **log_message** (int *log_level*, String *message*)
send user defined strings to BrainFlow logger

int **get_sampling_rate** (int *board_id*)
get sampling rate for this board

int **get_timestamp_channel** (int *board_id*)
get row index in returned by *get_board_data()* 2d array which contains timestamps

int **get_num_rows** (int *board_id*)
get number of rows in returned by *get_board_data()* 2d array

int **get_package_num_channel** (int *board_id*)
get row index in returned by *get_board_data()* 2d array which contains package nums

int **get_battery_channel** (int *board_id*)
get row index in returned by *get_board_data()* 2d array which contains battery level

String [] **get_eeg_names** (int *board_id*)
Get names of EEG electrodes in 10-20 system. Only if electrodes have freezed locations

int [] **get_eeg_channels** (int *board_id*)
get row indices in returned by *get_board_data()* 2d array which contain EEG data, for some boards we can not split EEG... and return the same array

int [] **get_emg_channels** (int *board_id*)
get row indices in returned by *get_board_data()* 2d array which contain EMG data, for some boards we can not split EEG... and return the same array

int [] **get_ecg_channels** (int *board_id*)
get row indices in returned by *get_board_data()* 2d array which contain ECG data, for some boards we can not split EEG... and return the same array

int [] **get_temperature_channels** (int *board_id*)
get row indices in returned by *get_board_data()* 2d array which contain temperature data

int [] **get_resistance_channels** (int *board_id*)
get row indices in returned by *get_board_data()* 2d array which contain resistance data

```
int [] get_eog_channels (int board_id)
    get row indices in returned by get_board_data() 2d array which contain EOG data, for some boards we
    can not split EEG... and return the same array

int [] get_exg_channels (int board_id)
    get row indices in returned by get_board_data() 2d array which contain EXG data

int [] get_eda_channels (int board_id)
    get row indices in returned by get_board_data() 2d array which contain EDA data, for some boards we
    can not split EEG... and return the same array

int [] get_ppg_channels (int board_id)
    get row indices in returned by get_board_data() 2d array which contain PPG data, for some boards we can
    not split EEG... and return the same array

int [] get_accel_channels (int board_id)
    get row indices in returned by get_board_data() 2d array which contain accel data

int [] get_analog_channels (int board_id)
    get row indices in returned by get_board_data() 2d array which contain analog data

int [] get_gyro_channels (int board_id)
    get row indices in returned by get_board_data() 2d array which contain gyro data

int [] get_other_channels (int board_id)
    get row indices in returned by get_board_data() 2d array which contain other data
```

```
enum BrainFlowClassifiers
```

Public Functions

```
int get_code ()

brainflow::BrainFlowClassifiers (final int code)
```

Public Members

```
brainflow::REGRESSION    = (0)
brainflow::KNN           = (1)
brainflow::SVM           = (2)
brainflow::LDA           = (3)
```

Public Static Functions

```
String brainflow::string_from_code (final int code)

BrainFlowClassifiers brainflow::from_code (final int code)

brainflow::[static initializer]
```

Private Members

```
final int brainflow::protocol
```

Private Static Attributes

```
final Map< Integer, BrainFlowClassifiers > brainflow::cl_map = new HashMap<Integer, I
class brainflow::brainflow::BrainFlowError : public Exception
    BrainFlowError exception to notify about errors
```

Public Functions

```
BrainFlowError (String message, int ec)
```

Public Members

```
String msg
```

```
int exit_code
    exit code returned from low level API
```

```
class brainflow::brainflow::BrainFlowInputParams
    to get fields which are required for your board check SupportedBoards section
```

Public Functions

```
BrainFlowInputParams ()
```

```
String to_json ()
```

```
String get_ip_address ()
```

```
void set_ip_address (String ip_address)
```

```
String get_mac_address ()
```

```
void set_mac_address (String mac_address)
```

```
String get_serial_port ()
```

```
void set_serial_port (String serial_port)
```

```
int get_ip_port ()
```

```
void set_ip_port (int ip_port)
```

```
int get_ip_protocol ()
```

```
void set_ip_protocol (int ip_protocol)
```

```
String get_other_info ()
```

```
void set_other_info (String other_info)
```

```
void set_timeout (int timeout)
```

```
int get_timeout ()
```

```
String get_serial_number ()
```

```
void set_serial_number (String serial_number)  
String get_file ()  
void set_file (String file)
```

Public Members

```
String ip_address  
String mac_address  
String serial_port  
int ip_port  
int ip_protocol  
String other_info  
int timeout  
String serial_number  
String file  
enum BrainFlowMetrics
```

Public Functions

```
int get_code ()  
brainflow::BrainFlowMetrics (final int code)
```

Public Members

```
brainflow::RELAXATION    = (0)  
brainflow::CONCENTRATION = (1)
```

Public Static Functions

```
String brainflow::string_from_code (final int code)  
BrainFlowMetrics brainflow::from_code (final int code)  
brainflow::[static initializer]
```

Private Members

```
final int brainflow::protocol
```

Private Static Attributes

```

    final Map< Integer, BrainFlowMetrics > brainflow::metr_map    = new HashMap<Integer, Br
class brainflow::brainflow::BrainFlowModelParams
    describe model parameters

```

Public Functions

```

BrainFlowModelParams (int metric, int classifier)

int get_metric ()
void set_metric (int metric)
int get_classifier ()
void set_classifier (int classifier)
String get_file ()
void set_file (String file)
String get_other_info ()
void set_other_info (String other_info)
String to_json ()

```

Public Members

```

int metric
int classifier
String file
String other_info
class brainflow::brainflow::DataFilter
    DataFilter class to perform signal processing

```

Public Static Functions

```

void enable_data_logger ()
    enable Data logger with level INFO
void enable_dev_data_logger ()
    enable Data logger with level TRACE
void disable_data_logger ()
    disable Data logger
void set_log_file (String log_file)
    redirect logger from stderr to a file

void perform_lowpass (double[] data, int sampling_rate, double cutoff, int order, int
    perform lowpass filter in-place

void perform_highpass (double[] data, int sampling_rate, double cutoff, int order, int
    perform highpass filter in-place

```

void perform_bandpass (double[] data, int sampling_rate, double center_freq, double band_width)
perform bandpass filter in-place

void perform_bandstop (double[] data, int sampling_rate, double center_freq, double band_width)
perform bandstop filter in-place

void perform_rolling_filter (double[] data, int period, int operation)
perform moving average or moving median filter in-place

void detrend (double[] data, int operation)
subtract trend from data in-place

double [] perform_downsampling (double[] data, int period, int operation)
perform data downsampling, it doesnt apply lowpass filter for you, it just aggregates several data points

void perform_wavelet_denoising (double[] data, String wavelet, int decomposition_level)
perform wavelet based denoising in-place

Parameters

- wavelet: supported vals: db1..db15,haar,sym2..sym10,coif1..coif5,bior1.1,bior1.3,bior1.5,bior2.2,bior2.4,bior2.6,bior3.7,bior3.9,bior4.4,bior5.5,bior6.8
- decomposition_level: level of decomposition of wavelet transform

Pair< double[], int[]> perform_wavelet_transform (double[] data, String wavelet, int decomposition_level)
perform wavelet transform

Parameters

- wavelet: supported vals: db1..db15,haar,sym2..sym10,coif1..coif5,bior1.1,bior1.3,bior1.5,bior2.2,bior2.4,bior2.6,bior3.7,bior3.9,bior4.4,bior5.5,bior6.8

double [] perform_inverse_wavelet_transform (Pair< double[], int[]> wavelet_output, int decomposition_level)
perform inverse wavelet transform

Complex [] perform_fft (double[] data, int start_pos, int end_pos, int window)
perform direct fft

Return array of complex values with size $N / 2 + 1$

Parameters

- data: data for fft transform
- start_pos: starting position to calc fft
- end_pos: end position to calc fft, total_len must be a power of two
- window: window function

double [] perform_ifft (Complex[] data)
perform inverse fft

Return restored data

Parameters

- data: data from fft transform(array of complex values)

Pair< double[], double[]> get_avg_band_powers (double[][] data, int[] channels, int sampling_rate)
calc average and stddev of band powers across all channels

Return pair of avgs and stddevs for bandpowers

Parameters

- `data`: data to process
- `channels`: rows of data arrays which should be used in calculation
- `sampling_rate`: sampling rate
- `apply_filters`: apply bandpass and bandstop filters before calculation

Pair< double[], double[]> get_psd (double[] data, int start_pos, int end_pos, int samp
get PSD

Return pair of ampl and freq arrays with len $N / 2 + 1$

Parameters

- `data`: data to process
- `start_pos`: starting position to calc PSD
- `end_pos`: end position to calc PSD, total_len must be a power of two
- `sampling_rate`: sampling rate
- `window`: window function

Pair< double[], double[]> get_psd_welch (double[] data, int nfft, int overlap, int samp
get PSD using Welch Method

Return pair of ampl and freq arrays

Parameters

- `data`: data to process
- `nfft`: size of FFT, must be power of two
- `overlap`: overlap between FFT Windows, must be between 0 and nfft
- `sampling_rate`: sampling rate
- `window`: window function

double get_band_power (Pair<double[], double[]> psd, double freq_start, double freq_end)
get band power

Return band power

Parameters

- `psd`: PSD from `get_psd` or `get_log_psd`
- `freq_start`: lowest frequency of band
- `freq_end`: highest frequency of band

int get_nearest_power_of_two (int value)
calculate nearest power of two

void write_file (double[][] data, String file_name, String file_mode)
write data to csv file, in file data will be transposed

double [][] read_file (String file_name)
read data from file, transpose it back to original format

enum DetrendOperations

enum to store all supported detrend operations

Public Functions

```
int get_code ()
```

```
brainflow::DetrendOperations (final int code)
```

Public Members

```
brainflow::NONE      = (0)
```

```
brainflow::CONSTANT   = (1)
```

```
brainflow::LINEAR     = (2)
```

Public Static Functions

```
String brainflow::string_from_code (final int code)
```

```
DetrendOperations brainflow::from_code (final int code)
```

```
brainflow::[static initializer]
```

Private Members

```
final int brainflow::detrend_operation
```

Private Static Attributes

```
final Map< Integer, DetrendOperations > brainflow::dt_map = new HashMap<Integer, Det.
```

enum ExitCode

Public Functions

```
int get_code ()
```

```
brainflow::ExitCode (final int code)
```

Public Members

```
brainflow::STATUS_OK      = (0)
```

```
brainflow::PORT_ALREADY_OPEN_ERROR = (1)
```

```
brainflow::UNABLE_TO_OPEN_PORT_ERROR = (2)
```

```
brainflow::SET_PORT_ERROR = (3)
```

```
brainflow::BOARD_WRITE_ERROR = (4)
```

```
brainflow::INCOMING_MSG_ERROR = (5)
```

```
brainflow::INITIAL_MSG_ERROR = (6)
```



```

brainflow::BOARD_NOT_READY_ERROR    =(7)
brainflow::STREAM_ALREADY_RUN_ERROR  =(8)
brainflow::INVALID_BUFFER_SIZE_ERROR =(9)
brainflow::STREAM_THREAD_ERROR       =(10)
brainflow::STREAM_THREAD_IS_NOT_RUNNING =(11)
brainflow::EMPTY_BUFFER_ERROR        =(12)
brainflow::INVALID_ARGUMENTS_ERROR   =(13)
brainflow::UNSUPPORTED_BOARD_ERROR   =(14)
brainflow::BOARD_NOT_CREATED_ERROR   =(15)
brainflow::ANOTHER_BOARD_IS_CREATED_ERROR =(16)
brainflow::GENERAL_ERROR              =(17)
brainflow::SYNC_TIMEOUT_ERROR         =(18)
brainflow::JSON_NOT_FOUND_ERROR       =(19)
brainflow::NO_SUCH_DATA_IN_JSON_ERROR =(20)
brainflow::CLASSIFIER_IS_NOT_PREPARED_ERROR =(21)
brainflow::ANOTHER_CLASSIFIER_IS_PREPARED_ERROR =(22)
brainflow::UNSUPPORTED_CLASSIFIER_AND_METRIC_COMBINATION_ERROR =(23)

```

Public Static Functions

```

String brainflow::string_from_code (final int code)
ExitCode brainflow::from_code (final int code)
brainflow::[static initializer]

```

Private Members

```
final int brainflow::exit_code
```

Private Static Attributes

```
final Map< Integer, ExitCode > brainflow::ec_map    = new HashMap<Integer, ExitCode> ()
enum FilterTypes
    enum to store all possible filter types

```

Public Functions

```
int get_code ()
brainflow::FilterTypes (final int code)
```

Public Members

```
brainflow::BUTTERWORTH    =(0)
brainflow::CHEBYSHEV_TYPE_1  =(1)
brainflow::BESSEL         =(2)
```

Public Static Functions

```
String brainflow::string_from_code (final int code)
FilterTypes brainflow::from_code (final int code)
brainflow::[static initializer]
```

Private Members

```
final int brainflow::filter_type
```

Private Static Attributes

```
final Map< Integer, FilterTypes > brainflow::ft_map    = new HashMap<Integer, FilterTypes>()
enum IpProtocolType
```

Public Functions

```
int get_code ()
brainflow::IpProtocolType (final int code)
```

Public Members

```
brainflow::NONE    =(0)
brainflow::UDP      =(1)
brainflow::TCP      =(2)
```

Public Static Functions

```
String brainflow::string_from_code (final int code)
IpProtocolType brainflow::from_code (final int code)
brainflow::[static initializer]
```

Private Members

```
final int brainflow::protocol
```

Private Static Attributes

```
final Map< Integer, IpProtocolType > brainflow::ip_map    = new HashMap<Integer, IpProt
enum LogLevels
```

Public Functions

```
int get_code ()
brainflow::LogLevels (final int code)
```

Public Members

```
brainflow::LEVEL_TRACE    =(0)
brainflow::LEVEL_DEBUG    =(1)
brainflow::LEVEL_INFO     =(2)
brainflow::LEVEL_WARN     =(3)
brainflow::LEVEL_ERROR    =(4)
brainflow::LEVEL_CRITICAL  =(5)
brainflow::LEVEL_OFF      =(6)
```

Public Static Functions

```
String brainflow::string_from_code (final int code)
LogLevels brainflow::from_code (final int code)
brainflow::[static initializer]
```

Private Members

```
final int brainflow::log_level
```

Private Static Attributes

```
final Map< Integer, LogLevels > brainflow::ll_map    = new HashMap<Integer, LogLevels>  
class brainflow::brainflow::MLModel
```

Public Functions

```
MLModel (BrainFlowModelParams params)  
    Create MLModel object
```

```
void prepare ()  
    Prepare classifier
```

Exceptions

- *BrainFlowError:*

```
void release ()  
    Release classifier
```

Exceptions

- *BrainFlowError:*

```
double predict (double[] data)  
    Get score of classifier
```

Exceptions

- *BrainFlowError:*

Public Static Functions

```
void enable_ml_logger ()  
    enable ML logger with level INFO
```

```
void enable_dev_ml_logger ()  
    enable ML logger with level TRACE
```

```
void disable_ml_logger ()  
    disable BrainFlow logger
```

```
void set_log_file (String log_file)  
    redirect logger from stderr to a file
```

```
enum WindowFunctions
```

Public Functions

```
int get_code ()
brainflow::WindowFunctions (final int code)
```

Public Members

```
brainflow::NO_WINDOW    =(0)
brainflow::HANNING      =(1)
brainflow::HAMMING      =(2)
brainflow::BLACKMAN_HARRIS =(3)
```

Public Static Functions

```
String brainflow::string_from_code (final int code)
WindowFunctions brainflow::from_code (final int code)
brainflow::[static initializer]
```

Private Members

```
final int brainflow::window
```

Private Static Attributes

```
final Map< Integer, WindowFunctions > brainflow::window_map = new HashMap<Integer, W
```

3.4 C# API Reference

To simplify 2D array manipulation we use [Accord Library](#).

Content of brainflow namespace:

```
enum brainflow::LogLevels
    Values:
        enumerator LEVEL_TRACE = 0
        enumerator LEVEL_DEBUG = 1
        enumerator LEVEL_INFO = 2
        enumerator LEVEL_WARN = 3
        enumerator LEVEL_ERROR = 4
        enumerator LEVEL_CRITICAL = 5
        enumerator LEVEL_OFF = 6
enum brainflow::CustomExitCodes
    Values:
```

```
enumerator STATUS_OK = 0
enumerator PORT_ALREADY_OPEN_ERROR = 1
enumerator UNABLE_TO_OPEN_PORT_ERROR = 2
enumerator SET_PORT_ERROR = 3
enumerator BOARD_WRITE_ERROR = 4
enumerator INCOMING_MSG_ERROR = 5
enumerator INITIAL_MSG_ERROR = 6
enumerator BOARD_NOT_READY_ERROR = 7
enumerator STREAM_ALREADY_RUN_ERROR = 8
enumerator INVALID_BUFFER_SIZE_ERROR = 9
enumerator STREAM_THREAD_ERROR = 10
enumerator STREAM_THREAD_IS_NOT_RUNNING = 11
enumerator EMPTY_BUFFER_ERROR = 12
enumerator INVALID_ARGUMENTS_ERROR = 13
enumerator UNSUPPORTED_BOARD_ERROR = 14
enumerator BOARD_NOT_CREATED_ERROR = 15
enumerator ANOTHER_BOARD_IS_CREATED_ERROR = 16
enumerator GENERAL_ERROR = 17
enumerator SYNC_TIMEOUT_ERROR = 18
enumerator JSON_NOT_FOUND_ERROR = 19
enumerator NO_SUCH_DATA_IN_JSON_ERROR = 20
enumerator CLASSIFIER_IS_NOT_PREPARED_ERROR = 21
enumerator ANOTHER_CLASSIFIER_IS_PREPARED_ERROR = 22
enumerator UNSUPPORTED_CLASSIFIER_AND_METRIC_COMBINATION_ERROR = 23

enum brainflow::BoardIds
    Values:
        enumerator PLAYBACK_FILE_BOARD = -3
        enumerator STREAMING_BOARD = -2
        enumerator SYNTHETIC_BOARD = -1
        enumerator CYTON_BOARD = 0
        enumerator GANGLION_BOARD = 1
        enumerator CYTON_DAISY_BOARD = 2
        enumerator GALEA_BOARD = 3
        enumerator GANGLION_WIFI_BOARD = 4
        enumerator CYTON_WIFI_BOARD = 5
        enumerator CYTON_DAISY_WIFI_BOARD = 6
```

```
    enumerator BRAINBIT_BOARD = 7
    enumerator UNICORN_BOARD = 8
    enumerator CALLIBRI_EEG_BOARD = 9
    enumerator CALLIBRI_EMG_BOARD = 10
    enumerator CALLIBRI_ECG_BOARD = 11
    enumerator FASCIA_BOARD = 12
    enumerator NOTION_1_BOARD = 13
    enumerator NOTION_2_BOARD = 14
    enumerator IRONBCI_BOARD = 15
    enumerator FREEEEG32_BOARD = 17

enum brainflow::IpProtocolType
    Values:
    enumerator NONE = 0
    enumerator UDP = 1
    enumerator TCP = 2

enum brainflow::FilterTypes
    Values:
    enumerator BUTTERWORTH = 0
    enumerator CHEBYSHEV_TYPE_1 = 1
    enumerator BESSEL = 2

enum brainflow::AggOperations
    Values:
    enumerator MEAN = 0
    enumerator MEDIAN = 1
    enumerator EACH = 2

enum brainflow::WindowFunctions
    Values:
    enumerator NO_WINDOW = 0
    enumerator HANNING = 1
    enumerator HAMMING = 2
    enumerator BLACKMAN_HARRIS = 3

enum brainflow::DetrendOperations
    Values:
    enumerator NONE = 0
    enumerator CONSTANT = 1
    enumerator LINEAR = 2

enum brainflow::BrainFlowMetrics
    Values:
```

```
    enumerator RELAXATION = 0
    enumerator CONCENTRATION = 1
enum brainflow::BrainFlowClassifiers
    Values:
    enumerator REGRESSION = 0
    enumerator KNN = 1
    enumerator SVM = 2
    enumerator LDA = 3
class brainflow::brainflow::BoardShim
    BoardShim class to communicate with a board
```

Public Functions

BoardShim (int *board_id*, *BrainFlowInputParams* *input_params*)
Create an instance of BoardShim class

Parameters

- *board_id*
- *input_params*

void **prepare_session** ()
prepare BrainFlow's streaming session, allocate required resources

string **config_board** (string *config*)
send string to a board, use this method carefully and only if you understand what you are doing

Parameters

- *config*

void **start_stream** (int *buffer_size* = 3600 * 250, string *streamer_params* = "")
start streaming thread, store data in internal ringbuffer

Parameters

- *buffer_size*: size of internal ringbuffer
- *streamer_params*: supported values: `file://%file_name%:w`, `file://%file_name%:a`, `streaming_board://multicast_group_ip%:port%`

void **stop_stream** ()
stop streaming thread, doesnt release other resources

void **release_session** ()
release BrainFlow's session

bool **is_prepared** ()
check session status

summary> Get Board Id, for some boards can be different than provided /summary>

Return session status

Return Master board id

int **get_board_id** ()

int **get_board_data_count** ()

get number of packages in ringbuffer

Return number of packages

double [,] **get_current_board_data** (int num_samples)

get latest collected data, doesnt remove it from ringbuffer

Return latest collected data, can be less than “num_samples”

Parameters

- num_samples

double [,] **get_board_data** ()

get all collected data and remove it from ringbuffer

Return all collected data

Public Members

int **board_id**

BrainFlow’s board id

Public Static Functions

int **get_sampling_rate** (int board_id)

get sampling rate for this board id

Return sampling rate

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int **get_package_num_channel** (int board_id)

get row index in returned by *get_board_data()* 2d array which hold package nums

Return row num in 2d array

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int **get_timestamp_channel** (int board_id)

get row index which hold timestamps

Return row num in 2d array

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int **get_battery_channel** (int board_id)
get row undex which holds battery level

Return row num in 2d array

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int **get_num_rows** (int board_id)
get number of rows in returned by *get_board_data()* 2d array

Return number of rows in 2d array

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

string [] **get_eeg_names** (int board_id)
get names of EEG channels in 10-20 system. Only if electrodes have fixed locations

Return array of 10-20 locations

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int [] **get_eeg_channels** (int board_id)
get row indices of EEG channels for this board, for some board we can not split EMG.. data and return the same array for all of them

Return array of row nums

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int [] get_emg_channels (int board_id)

get row indices of EMG channels for this board, for some board we can not split EMG.. data and return the same array for all of them

Return array of row nums

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int [] get_ecg_channels (int board_id)

get row indices of ECG channels for this board, for some board we can not split EMG.. data and return the same array for all of them

Return array of row nums

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int [] get_eog_channels (int board_id)

get row indices of EOG channels for this board, for some board we can not split EMG.. data and return the same array for all of them

Return array of row nums

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int [] get_exg_channels (int board_id)

get row indices of EXG channels for this board

Return array of row nums

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

```
int [] get_eda_channels (int board_id)
```

get row indices of EDA channels for this board, for some board we can not split EMG.. data and return the same array for all of them

Return array of row nums

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

```
int [] get_ppg_channels (int board_id)
```

get row indeces which hold ppg data

Return array of row nums

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

```
int [] get_accel_channels (int board_id)
```

get row indices which hold accel data

Return array of row nums

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

```
int [] get_analog_channels (int board_id)
```

get row indices which hold analog data

Return array of row nums

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

```
int [] get_gyro_channels (int board_id)
```

get row indices which hold gyro data

Return array of row nums

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int [] get_other_channels (int board_id)
get other channels for this board

Return array of row nums

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int [] get_temperature_channels (int board_id)
get temperature channels for this board

Return array of row nums

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

int [] get_resistance_channels (int board_id)
get resistance channels for this board

Return array of row nums

Parameters

- board_id

Exceptions

- *BrainFlowException*: If this board has no such data exit code is UNSUPPORTED_BOARD_ERROR

void set_log_level (int log_level)
set log level, logger is disabled by default

Parameters

- log_level

void enable_board_logger ()
enable BrainFlow's logger with level INFO

void disable_board_logger ()
disable BrainFlow's logger

void enable_dev_board_logger ()
enable BrainFlow's logger with level TRACE

void **set_log_file** (string *log_file*)
redirect BrainFlow's logger from stderr to file

Parameters

- *log_file*

void **log_message** (int *log_level*, string *message*)
send your own log message to BrainFlow's logger

Parameters

- *log_level*
- *message*

class `brainflow::brainflow::BrainFlowException` : **public** Exception
BrainFlowException class to notify about errors

Public Functions

BrainFlowException (int *code*)

Public Members

int **exit_code**
exit code returned from low level API

class `brainflow::brainflow::BrainFlowInputParams`
Check SupportedBoards to get information about fields which are required for specific board

Public Functions

BrainFlowInputParams ()

string **to_json** ()

Public Members

string **serial_port**
serial port name

string **mac_address**
MAC address

string **ip_address**
IP address

int **ip_port**
PORT

int **ip_protocol**
IP protocol, use IpProtocolType

string **other_info**
you can provide additional info to low level API using this field

```
int timeout
    timeout for device discovery or connection

string serial_number
    serial number

string file
    file
```

```
class brainflow::brainflow::BrainFlowModelParams
    Describe model
```

Public Functions

```
BrainFlowModelParams (int metric, int classifier)

string to_json ()
```

Public Members

```
int metric
    metric to calculate

int classifier
    classifier to use

string file
    path to model file

string other_info
    other info
```

```
class brainflow::brainflow::DataFilter
    DataFilter class to perform signal processing
```

Public Static Functions

```
void enable_data_logger ()
    enable Data logger with level INFO

void disable_data_logger ()
    disable Data logger

void enable_dev_data_logger ()
    enable Data logger with level TRACE

void set_log_file (string log_file)
    redirect BrainFlow's logger from stderr to file
```

Parameters

- *log_file*

```
double [] perform_lowpass (double[] data, int sampling_rate, double cutoff, int order,
    perform lowpass filter, unlike other bindings instead in-place calculation it returns new array
```

Return filtered data

Parameters

- data
- sampling_rate
- cutoff
- order
- filter_type
- ripple

double [] perform_highpass (double[] data, int sampling_rate, double cutoff, int order
perform highpass filter, unlike other bindings instead in-place calculation it returns new array

Return filtered data

Parameters

- data
- sampling_rate
- cutoff
- order
- filter_type
- ripple

double [] perform_bandpass (double[] data, int sampling_rate, double center_freq, doub
perform bandpass filter, unlike other bindings instead in-place calculation it returns new array

Return filtered data

Parameters

- data
- sampling_rate
- center_freq
- band_width
- order
- filter_type
- ripple

double [] perform_bandstop (double[] data, int sampling_rate, double center_freq, doub
perform bandstop filter, unlike other bindings instead in-place calculation it returns new array

Return filtered data

Parameters

- data
- sampling_rate
- center_freq
- band_width
- order

- filter_type
- ripple

double [] perform_rolling_filter (double[] data, int period, int operation)
 perform moving average or moving median filter, unlike other bindings instead in-place calculation it returns new array

Return filtered data

Parameters

- data
- period
- operation

double [] detrend (double[] data, int operation)
 detrend, unlike other bindings instead in-place calculation it returns new array

Return data with removed trend

Parameters

- data
- operation

double [] perform_downsampling (double[] data, int period, int operation)
 perform data downsampling, it just aggregates data without applying lowpass filter

Return data after downsampling

Parameters

- data
- period
- operation

Tuple< double[], int[]> perform_wavelet_transform (double[] data, string wavelet, int decomposition_level)
 perform wavelet transform

Return tuple of wavelet coeffs in format [A(J) D(J) D(J-1) ... D(1)] where J is decomposition level, A - app coeffs, D - detailed coeffs, and array with lengths for each block

Parameters

- data: data for wavelet transform
- wavelet: db1..db15,haar,sym2..sym10,coif1..coif5,bior1.1,bior1.3,bior1.5,bior2.2,bior2.4,bior2.6,bior2.8,bior3.1,bior3.3,bior3.5,bior3.7,bior3.9,bior4.4,bior5.5,bior6.8
- decomposition_level: decomposition level

double [] perform_inverse_wavelet_transform (Tuple< double[], int[]> wavelet_data, int decomposition_level)
 perform inverse wavelet transform

Return restored data

Parameters

- wavelet_data: tuple returned by perform_wavelet_transform

- `original_data_len`: size of original data before direct wavelet transform
- `wavelet`: db1..db15,haar,sym2..sym10,coif1..coif5,bior1.1,bior1.3,bior1.5,bior2.2,bior2.4,bior2.6,bior2.8,bior3.1,bior3.7,bior3.9,bior4.4,bior5.5,bior6.8
- `decomposition_level`: level of decomposition

double [] perform_wavelet_denoising (double[] data, string wavelet, int decomposition_level)
perform wavelet based denoising

Return denoised data

Parameters

- `data`: data for denoising
- `wavelet`: db1..db15,haar,sym2..sym10,coif1..coif5,bior1.1,bior1.3,bior1.5,bior2.2,bior2.4,bior2.6,bior2.8,bior3.1,bior3.7,bior3.9,bior4.4,bior5.5,bior6.8
- `decomposition_level`: level of decomposition in wavelet transform

Complex [] perform_fft (double[] data, int start_pos, int end_pos, int window)
perform direct fft

Return complex array of size $N / 2 + 1$ of fft data

Parameters

- `data`: data for fft
- `start_pos`: start pos
- `end_pos`: end pos, `end_pos - start_pos` must be a power of 2
- `window`: window function

double [] perform_ifft (Complex[] data)
perform inverse fft

Return restored data

Parameters

- `data`: data from `perform_fft`

void write_file (double[,] data, string file_name, string file_mode)
write data to csv file, data will be transposed

Parameters

- `data`
- `file_name`
- `file_mode`

double [,] read_file (string file_name)
read data from file, data will be transposed back to original format

Return

Parameters

- `file_name`

int `get_nearest_power_of_two` (int *value*)
calculate nearest power of two

Return nearest power of two

Parameters

- *value*

Tuple< double[], double[]> `get_avg_band_powers` (double[,] *data*, int[] *channels*, int *sampling_rate*)
calculate avg and stddev bandpowers across channels

Return Tuple of avgs and stddev arrays

Parameters

- *data*: 2d array with values
- *channels*: rows of data array which should be used for calculation
- *sampling_rate*: sampling rate
- *apply_filters*: apply bandpass and bandstop filters before calculation

Tuple< double[], double[]> `get_psd` (double[] *data*, int *start_pos*, int *end_pos*, int *sampling_rate*)
calculate PSD

Return Tuple of ampls and freqs arrays of size $N / 2 + 1$

Parameters

- *data*: data for PSD
- *start_pos*: start pos
- *end_pos*: end pos, *end_pos* - *start_pos* must be a power of 2
- *sampling_rate*: sampling rate
- *window*: window function

Tuple< double[], double[]> `get_psd_welch` (double[] *data*, int *nfft*, int *overlap*, int *sampling_rate*)
calculate PSD using Welch method

Return Tuple of ampls and freqs arrays

Parameters

- *data*: data for log PSD
- *nfft*: FFT Size
- *overlap*: FFT Window overlap, must be between 0 and *nfft*
- *sampling_rate*: sampling rate
- *window*: window function

double `get_band_power` (Tuple<double[], double[]> *psd*, double *start_freq*, double *stop_freq*)
calculate band power

Return band power

Parameters

- *psd*: psd data returned by `get_psd` or `get_psd_welch`

- `start_freq`: lowest frequency of band
- `stop_freq`: highest frequency of band

class `brainflow::brainflow::MLModel`

Public Functions

MLModel (*BrainFlowModelParams* `input_params`)
Create an instance of MLModel class

Parameters

- `input_params`

void **prepare** ()
Prepare classifier

void **release** ()
Release classifier

double **predict** (double[] `data`)
Get score of classifier

Public Static Functions

void **enable_ml_logger** ()
enable ML logger with level INFO

void **disable_ml_logger** ()
disable ML logger

void **enable_dev_ml_logger** ()
enable ML logger with level TRACE

void **set_log_file** (string *log_file*)
redirect BrainFlow's logger from stderr to file

Parameters

- `log_file`

3.5 R API Reference

R binding is a wrapper on top of Python binding. It is implemented using [reticulate](#).

Check R samples to see how to use it.

Full code for R binding:

```
#' @import reticulate
NULL

#' @export
brainflow_python <- NULL
#' @export
```

(continues on next page)

(continued from previous page)

```

np <- NULL
#' @export
pandas <- NULL
sys <- NULL
type_map <- NULL

.onLoad <- function (libname, pkgname)
{
  brainflow_python <-< import ('brainflow', delay_load = TRUE)
  np <-< import ('numpy', delay_load = TRUE)
  pandas <-< import ('pandas', delay_load = TRUE)
  sys <-< import ('sys', delay_load = TRUE)
  type_map <-< function (type)
  {
    if (is.character (type))
    {
      return (list (
        'float32' = np$float32,
        'float64' = np$float64,
        'auto' = NULL
      )[[type]])
    }
    type
  }
}

```

3.6 Matlab API Reference

Matlab binding calls CC++ code as any other binding, it's not compatible with Octave. Use Matlab examples and API reference for other languages as a starting point.

A few general rules to keep in mind:

- Use char arrays instead strings to work with BrainFlow API, it means 'my_string' instead "my_string", otherwise you will get calllib error
- Use int32 values instead enums, it means int32 (BoardIDs.SYNTHETIC_BOARD) instead BoardIDs.SYNTHETIC_BOARD, the same is true for all enums in BrainFlow API

Like here:

```

BoardShim.set_log_file ('brainflow.log');
BoardShim.enable_dev_board_logger ();

params = BrainFlowInputParams ();
board_shim = BoardShim (int32 (BoardIDs.SYNTHETIC_BOARD), params);
board_shim.prepare_session ();
board_shim.start_stream (45000, '');
pause (2)
board_shim.stop_stream ()
data = board_shim.get_current_board_data (20);
board_shim.release_session ();

DataFilter.write_file (data, 'data.csv', 'w');
restored_data = DataFilter.read_file ('data.csv');

```

3.7 Julia API Reference

Julia binding calls CC++ code as any other binding. Use Julia examples and API reference for other languages as a starting point.

Since Julia is not Object-Oriented language, there is no DataFilter class. BoardShim class exists but all BoardShim class methods were moved to BrainFlow package and you need to pass BoardShim object to them.

Like here:

```
import brainflow

# specify logging library to use
brainflow.enable_dev_brainflow_logger(Integer(brainflow.BOARD_CONTROLLER))

params = brainflow.BrainFlowInputParams()
board_shim = brainflow.BoardShim(Integer(brainflow.SYNTHETIC_BOARD), params)

brainflow.prepare_session(board_shim)
brainflow.start_stream(board_shim)
sleep(5)
brainflow.stop_stream(board_shim)
data = brainflow.get_current_board_data(32, board_shim)
brainflow.release_session(board_shim)

brainflow.write_file(data, "test.csv", "w")
restored_data = brainflow.read_file("test.csv")

println("Original Data")
println(data)
println("Restored Data")
println(restored_data)
```

DATA FORMAT DESCRIPTION

4.1 Units of Measure

For EEG, EMG, etc BrainFlow returns uV.

For timestamps BrainFlow uses UNIX timestamp, this count starts at the Unix Epoch on January 1st, 1970 at UTC. Precision is microsecond, but for some boards timestamps are generated on PC side as soon as package was received.

You can compare BrainFlow's timestamp with time returned by code like this:

```
import time
print (time.time ())
```

4.2 Generic Format Description

Methods like:

```
get_board_data ()
get_current_board_data (max_num_packages)
```

Return 2d double array [num_channels x num_data_points], rows of this array represent different channels like EEG channels, EMG channels, Accel channels, Timesteps and so on, while columns in this array represent actual packages from a board.

Exact format for this array is board specific. To keep the API uniform, we have methods like:

```
# these methods return an array of rows in this 2d array containing eeg\emg\ecg\accel_
↪data
get_eeg_channels (board_id)
get_emg_channels (board_id)
get_ecg_channels (board_id)
get_accel_channels (board_id)
# and so on, check docs for full list
# also we have methods to get sampling rate from board id, get number of timestamp_
↪channel and others
get_sampling_rate (board_id)
get_timestamp_channel (board_id)
# and so on
```

For some boards like OpenBCI Cyton, OpenBCI Ganglion, etc we cannot separate EMG, EEG, EDA and ECG and in this case we return exactly the same array for all these methods but for some devices EMG and EEG channels will differ.

If board has no such data these methods throw an exception with `UNSUPPORTED_BOARD_ERROR` exit code.

Using the methods above, you can write completely board agnostic code and switch boards using a single parameter! Even if you have only one board using these methods you can easily switch to Synthetic Board or Streaming Board.

4.3 OpenBCI Specific Data

4.3.1 Special Channels for OpenBCI Cyton Based Boards

Cyton-based boards from OpenBCI support different output formats, described [here](#).

For Cyton based boards, we add Cyton End byte to a first channel from:

```
get_other_channels (board_id)
```

If Cyton End byte is equal to 0xC0 we add accel data. To get rows which contain accel data use:

```
get_accel_channels (board_id)
```

If Cyton End byte is equal to 0xC1 we add analog data. To get rows which contain analog data use:

```
get_analog_channels (board_id)
```

For analog data, we return int32 values. But since from low level API, we return double array, these values are converted to double without any changes.

Also we add raw unprocessed bytes to the second and next channels returned by:

```
get_other_channels (board_id)
```

If Cyton End Byte is outside [this range](#), we drop the entire package.

Check this example for details:

```
import argparse
import time
import numpy as np

import brainflow
from brainflow.board_shim import BoardShim, BrainFlowInputParams, BoardIds
from brainflow.data_filter import DataFilter, FilterTypes

def main ():
    parser = argparse.ArgumentParser ()
    # use docs to check which parameters are required for specific board, e.g. for
    ↪ Cyton - set serial port
    parser.add_argument ('--ip-port', type = int, help = 'ip port', required = False,
    ↪ default = 0)
    parser.add_argument ('--ip-protocol', type = int, help = 'ip protocol, check
    ↪ IpProtocolType enum', required = False, default = 0)
    parser.add_argument ('--ip-address', type = str, help = 'ip address', required =
    ↪ False, default = '')
    parser.add_argument ('--serial-port', type = str, help = 'serial port', required
    ↪ = False, default = '')
    parser.add_argument ('--mac-address', type = str, help = 'mac address', required
    ↪ = False, default = '')
```

(continues on next page)

(continued from previous page)

```

    parser.add_argument ('--other-info', type = str, help = 'other info', required =
↪False, default = '')
    parser.add_argument ('--streamer-params', type = str, help = 'other info',
↪required = False, default = '')
    parser.add_argument ('--board-id', type = int, help = 'board id, check docs to
↪get a list of supported boards', required = True)
    parser.add_argument ('--log', action = 'store_true')
    args = parser.parse_args ()

    params = BrainFlowInputParams ()
    params.ip_port = args.ip_port
    params.serial_port = args.serial_port
    params.mac_address = args.mac_address
    params.other_info = args.other_info
    params.ip_address = args.ip_address
    params.ip_protocol = args.ip_protocol

    if (args.log):
        BoardShim.enable_dev_board_logger ()
    else:
        BoardShim.disable_board_logger ()

    board = BoardShim (args.board_id, params)
    board.prepare_session ()

    board.start_stream ()
    time.sleep (5)
    board.config_board ('/2') # enable analog mode only for Cyton Based Boards!
    time.sleep (5)
    data = board.get_board_data ()
    board.stop_stream ()
    board.release_session ()

    """
    data[BoardShim.get_other_channels(args.board_id)[0]] contains cyton end byte
    data[BoardShim.get_other_channels(args.board_id)[1...]] contains unprocessed
↪bytes
    if end byte is 0xC0 there are accel data in data[BoardShim.get_accel_
↪channels(args.board_id)[...]] else there are zeros
    if end byte is 0xC1 there are analog data in data[BoardShim.get_analog_
↪channels(args.board_id)[...]] else there are zeros
    """
    print (data[BoardShim.get_other_channels(args.board_id)[0]][0:5]) # should be
↪standard end byte 0xC0
    print (data[BoardShim.get_other_channels(args.board_id)[0]][-5:]) # should be
↪analog and byte 0xC1

    DataFilter.write_file (data, 'cyton_data.csv', 'w')

if __name__ == "__main__":
    main ()

```


CODE SAMPLES

Make sure that you've installed BrainFlow package before running the code samples below.

See *Installation Instructions* for details.

5.1 Python

To run some signal processing samples, you may need to install:

- matplotlib
- pandas
- mne

BrainFlow doesn't use these packages and doesn't install them, but the packages will be used in demos below.

5.1.1 Python Get Data from a Board

```
import argparse
import time
import numpy as np

import brainflow
from brainflow.board_shim import BoardShim, BrainFlowInputParams
from brainflow.data_filter import DataFilter, FilterTypes, AggOperations

def main ():
    parser = argparse.ArgumentParser ()
    # use docs to check which parameters are required for specific board, e.g. for
    ↪ Cyton - set serial port
    parser.add_argument ('--timeout', type = int, help = 'timeout for device
    ↪ discovery or connection', required = False, default = 0)
    parser.add_argument ('--ip-port', type = int, help = 'ip port', required = False,
    ↪ default = 0)
    parser.add_argument ('--ip-protocol', type = int, help = 'ip protocol, check
    ↪ IpProtocolType enum', required = False, default = 0)
    parser.add_argument ('--ip-address', type = str, help = 'ip address', required =
    ↪ False, default = '')
    parser.add_argument ('--serial-port', type = str, help = 'serial port', required
    ↪ = False, default = '')
    parser.add_argument ('--mac-address', type = str, help = 'mac address', required
    ↪ = False, default = '')
```

(continues on next page)

(continued from previous page)

```

    parser.add_argument ('--other-info', type = str, help = 'other info', required =
↪False, default = '')
    parser.add_argument ('--streamer-params', type = str, help = 'streamer params',
↪required = False, default = '')
    parser.add_argument ('--serial-number', type = str, help = 'serial number',
↪required = False, default = '')
    parser.add_argument ('--board-id', type = int, help = 'board id, check docs to
↪get a list of supported boards', required = True)
    parser.add_argument ('--file', type = str, help = 'file', required = False,
↪default = '')
    parser.add_argument ('--log', action = 'store_true')
    args = parser.parse_args ()

    params = BrainFlowInputParams ()
    params.ip_port = args.ip_port
    params.serial_port = args.serial_port
    params.mac_address = args.mac_address
    params.other_info = args.other_info
    params.serial_number = args.serial_number
    params.ip_address = args.ip_address
    params.ip_protocol = args.ip_protocol
    params.timeout = args.timeout
    params.file = args.file

    if (args.log):
        BoardShim.enable_dev_board_logger ()
    else:
        BoardShim.disable_board_logger ()

    board = BoardShim (args.board_id, params)
    board.prepare_session ()

    # board.start_stream () # use this for default options
    board.start_stream (45000, args.streamer_params)
    time.sleep (10)
    # data = board.get_current_board_data (256) # get latest 256 packages or less,
↪doesn't remove them from internal buffer
    data = board.get_board_data () # get all data and remove it from internal buffer
    board.stop_stream ()
    board.release_session ()

    print (data)

if __name__ == "__main__":
    main ()

```

5.1.2 Python Read Write File

```
import argparse
import time
import numpy as np
import pandas as pd

import brainflow
from brainflow.board_shim import BoardShim, BrainFlowInputParams, LogLevels, BoardIds
from brainflow.data_filter import DataFilter, FilterTypes, AggOperations

def main ():
    BoardShim.enable_dev_board_logger ()

    # use synthetic board for demo
    params = BrainFlowInputParams ()
    board = BoardShim (BoardIds.SYNTHETIC_BOARD.value, params)
    board.prepare_session ()
    board.start_stream ()
    BoardShim.log_message (LogLevels.LEVEL_INFO.value, 'start sleeping in the main_
↳thread')
    time.sleep (10)
    data = board.get_current_board_data (20) # get 20 latest data points dont remove_
↳them from internal buffer
    board.stop_stream ()
    board.release_session ()

    # demo how to convert it to pandas DF and plot data
    eeg_channels = BoardShim.get_eeg_channels (BoardIds.SYNTHETIC_BOARD.value)
    df = pd.DataFrame (np.transpose (data))
    print ('Data From the Board')
    print (df.head (10))

    # demo for data serialization using brainflow API, we recommend to use it instead_
↳pandas.to_csv()
    DataFilter.write_file (data, 'test.csv', 'w') # use 'a' for append mode
    restored_data = DataFilter.read_file ('test.csv')
    restored_df = pd.DataFrame (np.transpose (restored_data))
    print ('Data From the File')
    print (restored_df.head (10))

if __name__ == "__main__":
    main ()
```

5.1.3 Python Downsample Data

```
import time
import numpy as np

import brainflow
from brainflow.board_shim import BoardShim, BrainFlowInputParams, LogLevels, BoardIds
from brainflow.data_filter import DataFilter, FilterTypes, AggOperations
```

(continues on next page)

(continued from previous page)

```

def main ():
    BoardShim.enable_dev_board_logger ()

    # use synthetic board for demo
    params = BrainFlowInputParams ()
    board = BoardShim (BoardIds.SYNTHETIC_BOARD.value, params)
    board.prepare_session ()
    board.start_stream ()
    BoardShim.log_message (LogLevels.LEVEL_INFO.value, 'start sleeping in the main_
↳thread')
    time.sleep (10)
    data = board.get_current_board_data (20) # get 20 latest data points dont remove_
↳them from internal buffer
    board.stop_stream ()
    board.release_session ()

    eeg_channels = BoardShim.get_eeg_channels (BoardIds.SYNTHETIC_BOARD.value)
    # demo for downsampling, it just aggregates data
    for count, channel in enumerate (eeg_channels):
        print ('Original data for channel %d:' % channel)
        print (data[channel])
        if count == 0:
            downsampled_data = DataFilter.perform_downsampling (data[channel], 3,
↳AggOperations.MEDIAN.value)
            elif count == 1:
                downsampled_data = DataFilter.perform_downsampling (data[channel], 2,
↳AggOperations.MEAN.value)
            else:
                downsampled_data = DataFilter.perform_downsampling (data[channel], 2,
↳AggOperations.EACH.value)
            print ('Downsampled data for channel %d:' % channel)
            print (downsampled_data)

if __name__ == "__main__":
    main ()

```

5.1.4 Python Transforms

```

import argparse
import time
import brainflow
import numpy as np

from brainflow.board_shim import BoardShim, BrainFlowInputParams, LogLevels, BoardIds
from brainflow.data_filter import DataFilter, FilterTypes, AggOperations,
↳WindowFunctions

def main ():
    BoardShim.enable_dev_board_logger ()

    # use synthetic board for demo
    params = BrainFlowInputParams ()

```

(continues on next page)

(continued from previous page)

```

board_id = BoardIds.SYNTHETIC_BOARD.value
sampling_rate = BoardShim.get_sampling_rate (board_id)
board = BoardShim (board_id, params)
board.prepare_session ()
board.start_stream ()
BoardShim.log_message (LogLevels.LEVEL_INFO.value, 'start sleeping in the main_
↪thread')
time.sleep (10)
data = board.get_current_board_data (DataFilter.get_nearest_power_of_two_
↪(sampling_rate))
board.stop_stream ()
board.release_session ()

eeg_channels = BoardShim.get_eeg_channels (board_id)
# demo for transforms
for count, channel in enumerate (eeg_channels):
    print ('Original data for channel %d:' % channel)
    print (data[channel])
    # demo for wavelet transforms
    # wavelet_coeffs format is[A(J) D(J) D(J-1) ..... D(1)] where J is_
↪decomposition level, A - app coeffs, D - detailed coeffs
    # lengths array stores lengths for each block
    wavelet_coeffs, lengths = DataFilter.perform_wavelet_transform (data[channel],
↪ 'db5', 3)
    app_coefs = wavelet_coeffs[0: lengths[0]]
    detailed_coeffs_first_block = wavelet_coeffs[lengths[0] : lengths[1]]
    # you can do smth with wavelet coeffs here, for example denoising works via_
↪thresholds
    # for wavelets coefficients
    restored_data = DataFilter.perform_inverse_wavelet_transform ((wavelet_coeffs,
↪ lengths), data[channel].shape[0], 'db5', 3)
    print ('Restored data after wavelet transform for channel %d:' % channel)
    print (restored_data)

    # demo for fft, len of data must be a power of 2
    fft_data = DataFilter.perform_fft (data[channel], WindowFunctions.NO_WINDOW.
↪value)
    # len of fft_data is N / 2 + 1
    restored_fft_data = DataFilter.perform_ifft (fft_data)
    print ('Restored data after fft for channel %d:' % channel)
    print (restored_fft_data)

if __name__ == "__main__":
    main ()

```

5.1.5 Python Signal Filtering

```

import argparse
import time
import brainflow
import numpy as np

import pandas as pd
import matplotlib
matplotlib.use ('Agg')
import matplotlib.pyplot as plt

from brainflow.board_shim import BoardShim, BrainFlowInputParams, LogLevels, BoardIds
from brainflow.data_filter import DataFilter, FilterTypes, AggOperations

def main ():
    BoardShim.enable_dev_board_logger ()

    # use synthetic board for demo
    params = BrainFlowInputParams ()
    board_id = BoardIds.SYNTHETIC_BOARD.value
    board = BoardShim (board_id, params)
    board.prepare_session ()
    board.start_stream ()
    BoardShim.log_message (LogLevels.LEVEL_INFO.value, 'start sleeping in the main_
↪thread')
    time.sleep (10)
    data = board.get_board_data ()
    board.stop_stream ()
    board.release_session ()

    # demo how to convert it to pandas DF and plot data
    eeg_channels = BoardShim.get_eeg_channels (board_id)
    df = pd.DataFrame (np.transpose (data))
    plt.figure ()
    df[eeg_channels].plot (subplots = True)
    plt.savefig ('before_processing.png')

    # for demo apply different filters to different channels, in production choose one
    for count, channel in enumerate (eeg_channels):
        # filters work in-place
        if count == 0:
            DataFilter.perform_bandpass (data[channel], BoardShim.get_sampling_rate_
↪(board_id), 15.0, 6.0, 4, FilterTypes.BESSEL.value, 0)
        elif count == 1:
            DataFilter.perform_bandstop (data[channel], BoardShim.get_sampling_rate_
↪(board_id), 30.0, 1.0, 3, FilterTypes.BUTTERWORTH.value, 0)
        elif count == 2:
            DataFilter.perform_lowpass (data[channel], BoardShim.get_sampling_rate_
↪(board_id), 20.0, 5, FilterTypes.CHEBYSHEV_TYPE_1.value, 1)
        elif count == 3:
            DataFilter.perform_highpass (data[channel], BoardShim.get_sampling_rate_
↪(board_id), 3.0, 4, FilterTypes.BUTTERWORTH.value, 0)
        elif count == 4:
            DataFilter.perform_rolling_filter (data[channel], 3, AggOperations.MEAN.
↪value)

```

(continues on next page)

(continued from previous page)

```

        elif count == 5:
            DataFilter.perform_rolling_filter (data[channel], 3, AggOperations.MEDIAN.
↪value)

    df = pd.DataFrame (np.transpose (data))
    plt.figure ()
    df[eeg_channels].plot (subplots = True)
    plt.savefig ('after_processing.png')

if __name__ == "__main__":
    main ()

```

5.1.6 Python Denoising

```

import argparse
import time
import brainflow
import numpy as np

import pandas as pd
import matplotlib
matplotlib.use ('Agg')
import matplotlib.pyplot as plt

from brainflow.board_shim import BoardShim, BrainFlowInputParams, LogLevels, BoardIds
from brainflow.data_filter import DataFilter, FilterTypes, AggOperations

def main ():
    BoardShim.enable_dev_board_logger ()

    # use synthetic board for demo
    params = BrainFlowInputParams ()
    board_id = BoardIds.SYNTHETIC_BOARD.value
    board = BoardShim (board_id, params)
    board.prepare_session ()
    board.start_stream ()
    BoardShim.log_message (LogLevels.LEVEL_INFO.value, 'start sleeping in the main_
↪thread')
    time.sleep (20)
    data = board.get_board_data ()
    board.stop_stream ()
    board.release_session ()

    # demo how to convert it to pandas DF and plot data
    eeg_channels = BoardShim.get_eeg_channels (board_id)
    df = pd.DataFrame (np.transpose (data))
    plt.figure ()
    df[eeg_channels].plot (subplots = True)
    plt.savefig ('before_processing.png')

    # demo for denoising, apply different methods to different channels for demo
    for count, channel in enumerate (eeg_channels):
        # first of all you can try simple moving median or moving average with_
↪different window size

```

(continues on next page)

(continued from previous page)

```

        if count == 0:
            DataFilter.perform_rolling_filter (data[channel], 3, AggOperations.MEAN.
↪value)
        elif count == 1:
            DataFilter.perform_rolling_filter (data[channel], 3, AggOperations.MEDIAN.
↪value)
        # if methods above dont work for your signal you can try wavelet based_
↪denoising
        # feel free to try different functions and decomposition levels
        elif count == 2:
            DataFilter.perform_wavelet_denoising (data[channel], 'db6', 3)
        elif count == 3:
            DataFilter.perform_wavelet_denoising (data[channel], 'bior3.9', 3)
        elif count == 4:
            DataFilter.perform_wavelet_denoising (data[channel], 'sym7', 3)
        elif count == 5:
            # with synthetic board this one looks like the best option, but it_
↪depends on many circumstances
            DataFilter.perform_wavelet_denoising (data[channel], 'coif3', 3)

df = pd.DataFrame (np.transpose (data))
plt.figure ()
df[eeg_channels].plot (subplots = True)
plt.savefig ('after_processing.png')

if __name__ == "__main__":
    main ()

```

5.1.7 Python MNE Integration

```

import time
import numpy as np
import matplotlib
matplotlib.use ('Agg')
import matplotlib.pyplot as plt
import pandas as pd

import brainflow
from brainflow.board_shim import BoardShim, BrainFlowInputParams, BoardIds

import mne
from mne.channels import read_layout

def main():
    BoardShim.enable_dev_board_logger ()
    # use synthetic board for demo
    params = BrainFlowInputParams ()
    board = BoardShim (BoardIds.SYNTHETIC_BOARD.value, params)
    board.prepare_session ()
    board.start_stream ()
    time.sleep (10)
    data = board.get_board_data ()

```

(continues on next page)

(continued from previous page)

```

board.stop_stream ()
board.release_session ()

eeg_channels = BoardShim.get_eeg_channels (BoardIds.SYNTHETIC_BOARD.value)
eeg_data = data[eeg_channels, :]
eeg_data = eeg_data / 1000000 # BrainFlow returns uV, convert to V for MNE

# Creating MNE objects from brainflow data arrays
ch_types = ['eeg'] * len (eeg_channels)
ch_names = BoardShim.get_eeg_names (BoardIds.SYNTHETIC_BOARD.value)
sfreq = BoardShim.get_sampling_rate (BoardIds.SYNTHETIC_BOARD.value)
info = mne.create_info (ch_names = ch_names, sfreq = sfreq, ch_types = ch_types)
raw = mne.io.RawArray (eeg_data, info)
# its time to plot something!
raw.plot_psd (average = True)
plt.savefig ('psd.png')

if __name__ == '__main__':
    main ()

```

5.1.8 Python Band Power

```

import argparse
import time
import brainflow
import numpy as np

from brainflow.board_shim import BoardShim, BrainFlowInputParams, LogLevels, BoardIds
from brainflow.data_filter import DataFilter, FilterTypes, AggOperations, \
↳ WindowFunctions, DetrendOperations

def main ():
    BoardShim.enable_dev_board_logger ()

    # use synthetic board for demo
    params = BrainFlowInputParams ()
    board_id = BoardIds.SYNTHETIC_BOARD.value
    sampling_rate = BoardShim.get_sampling_rate (board_id)
    board = BoardShim (board_id, params)
    board.prepare_session ()
    board.start_stream ()
    BoardShim.log_message (LogLevels.LEVEL_INFO.value, 'start sleeping in the main_
↳ thread')
    time.sleep (10)
    nfft = DataFilter.get_nearest_power_of_two (sampling_rate)
    data = board.get_board_data ()
    board.stop_stream ()
    board.release_session ()

    eeg_channels = BoardShim.get_eeg_channels (board_id)
    # second eeg channel of synthetic board is a sine wave at 10Hz, should see huge_
↳ alpha
    eeg_channel = eeg_channels[1]

```

(continues on next page)

(continued from previous page)

```

# optional detrend
DataFilter.detrend (data[eeg_channel], DetrendOperations.LINEAR.value)
psd = DataFilter.get_psd_welch (data[eeg_channel], nfft, nfft // 2, sampling_rate,
↪ WindowFunctions.BLACKMAN_HARRIS.value)

band_power_alpha = DataFilter.get_band_power (psd, 7.0, 13.0)
band_power_beta = DataFilter.get_band_power (psd, 14.0, 30.0)
print ("alpha/beta:%f", band_power_alpha / band_power_beta)

# fail test if ratio is not smth we expect
if (band_power_alpha / band_power_beta < 100):
    raise ValueError ('Wrong Ratio')

if __name__ == "__main__":
    main ()

```

5.1.9 Python EEG Metrics

```

import argparse
import time
import brainflow
import numpy as np

from brainflow.board_shim import BoardShim, BrainFlowInputParams, LogLevels, BoardIds,
↪ BrainFlowError
from brainflow.data_filter import DataFilter, FilterTypes, AggOperations,
↪ WindowFunctions, DetrendOperations
from brainflow.ml_model import MLModel, BrainFlowMetrics, BrainFlowClassifiers,
↪ BrainFlowModelParams
from brainflow.exit_codes import *

def main ():
    BoardShim.enable_board_logger ()
    DataFilter.enable_data_logger ()
    MLModel.enable_ml_logger ()

    parser = argparse.ArgumentParser ()
    # use docs to check which parameters are required for specific board, e.g. for
↪ Cyton - set serial port
    parser.add_argument ('--timeout', type = int, help = 'timeout for device
↪ discovery or connection', required = False, default = 0)
    parser.add_argument ('--ip-port', type = int, help = 'ip port', required = False,
↪ default = 0)
    parser.add_argument ('--ip-protocol', type = int, help = 'ip protocol, check
↪ IpProtocolType enum', required = False, default = 0)
    parser.add_argument ('--ip-address', type = str, help = 'ip address', required =
↪ False, default = '')
    parser.add_argument ('--serial-port', type = str, help = 'serial port', required
↪ = False, default = '')
    parser.add_argument ('--mac-address', type = str, help = 'mac address', required
↪ = False, default = '')
    parser.add_argument ('--other-info', type = str, help = 'other info', required =
↪ False, default = '')

```

(continues on next page)

(continued from previous page)

```

    parser.add_argument ('--streamer-params', type = str, help = 'streamer params',
↳required = False, default = '')
    parser.add_argument ('--serial-number', type = str, help = 'serial number',
↳required = False, default = '')
    parser.add_argument ('--board-id', type = int, help = 'board id, check docs to
↳get a list of supported boards', required = True)
    parser.add_argument ('--file', type = str, help = 'file', required = False,
↳default = '')
    args = parser.parse_args ()

    params = BrainFlowInputParams ()
    params.ip_port = args.ip_port
    params.serial_port = args.serial_port
    params.mac_address = args.mac_address
    params.other_info = args.other_info
    params.serial_number = args.serial_number
    params.ip_address = args.ip_address
    params.ip_protocol = args.ip_protocol
    params.timeout = args.timeout
    params.file = args.file

    board = BoardShim (args.board_id, params)
    master_board_id = board.get_board_id ()
    sampling_rate = BoardShim.get_sampling_rate (master_board_id)
    board.prepare_session ()
    board.start_stream (45000, args.streamer_params)
    BoardShim.log_message (LogLevels.LEVEL_INFO.value, 'start sleeping in the main
↳thread')
    time.sleep (5) # recommended window size for eeg metric calculation is at least 4
↳seconds, bigger is better
    data = board.get_board_data ()
    board.stop_stream ()
    board.release_session ()

    eeg_channels = BoardShim.get_eeg_channels (int (master_board_id))
    bands = DataFilter.get_avg_band_powers (data, eeg_channels, sampling_rate, True)
    feature_vector = np.concatenate ((bands[0], bands[1]))
    print (feature_vector)

    # calc concentration
    concentration_params = BrainFlowModelParams (BrainFlowMetrics.CONCENTRATION.value,
↳BrainFlowClassifiers.KNN.value)
    concentration = MLModel (concentration_params)
    concentration.prepare ()
    print ('Concentration: %f' % concentration.predict (feature_vector))
    concentration.release ()

    # calc relaxation
    relaxation_params = BrainFlowModelParams (BrainFlowMetrics.RELAXATION.value,
↳BrainFlowClassifiers.REGRESSION.value)
    relaxation = MLModel (relaxation_params)
    relaxation.prepare ()
    print ('Relaxation: %f' % relaxation.predict (feature_vector))
    relaxation.release ()

if __name__ == "__main__":
    main ()

```

5.2 Java

5.2.1 Java Get Data from a Board

```

package brainflow.examples;

import java.util.Arrays;

import brainflow.BoardShim;
import brainflow.BrainFlowInputParams;
import brainflow.LogLevels;

public class BrainFlowGetData
{
    public static void main (String[] args) throws Exception
    {
        BoardShim.enable_board_logger ();
        BrainFlowInputParams params = new BrainFlowInputParams ();
        int board_id = parse_args (args, params);
        BoardShim board_shim = new BoardShim (board_id, params);

        board_shim.prepare_session ();
        // board_shim.start_stream (); // use this for default options
        board_shim.start_stream (450000, "file://file_stream.csv:w");
        BoardShim.log_message (LogLevels.LEVEL_INFO.get_code (), "Start sleeping in_
↪the main thread");
        Thread.sleep (5000);
        board_shim.stop_stream ();
        System.out.println (board_shim.get_board_data_count ());
        double[][] data = board_shim.get_current_board_data (30); // doesnt flush it_
↪from ring buffer
        // double[][] data = board_shim.get_board_data (); // get all data and flush
        // from ring buffer
        for (int i = 0; i < data.length; i++)
        {
            System.out.println (Arrays.toString (data[i]));
        }
        board_shim.release_session ();
    }

    private static int parse_args (String[] args, BrainFlowInputParams params)
    {
        int board_id = -1;
        for (int i = 0; i < args.length; i++)
        {
            if (args[i].equals ("--ip-address"))
            {
                params.ip_address = args[i + 1];
            }
            if (args[i].equals ("--serial-port"))
            {
                params.serial_port = args[i + 1];
            }
            if (args[i].equals ("--ip-port"))
            {

```

(continues on next page)

(continued from previous page)

```

        params.ip_port = Integer.parseInt (args[i + 1]);
    }
    if (args[i].equals ("--ip-protocol"))
    {
        params.ip_protocol = Integer.parseInt (args[i + 1]);
    }
    if (args[i].equals ("--other-info"))
    {
        params.other_info = args[i + 1];
    }
    if (args[i].equals ("--board-id"))
    {
        board_id = Integer.parseInt (args[i + 1]);
    }
    if (args[i].equals ("--timeout"))
    {
        params.timeout = Integer.parseInt (args[i + 1]);
    }
    if (args[i].equals ("--serial-number"))
    {
        params.serial_number = args[i + 1];
    }
    if (args[i].equals ("--file"))
    {
        params.file = args[i + 1];
    }
    }
    return board_id;
}
}

```

5.2.2 Java Read Write File

```

package brainflow.examples;

import java.util.Arrays;

import brainflow.BoardIds;
import brainflow.BoardShim;
import brainflow.BrainFlowInputParams;
import brainflow.DataFilter;
import brainflow.LogLevels;

public class Serialization
{
    public static void main (String[] args) throws Exception
    {
        // use Synthetic board for demo
        BoardShim.enable_board_logger ();
        BrainFlowInputParams params = new BrainFlowInputParams ();
        int board_id = BoardIds.SYNTHETIC_BOARD.get_code ();
        BoardShim board_shim = new BoardShim (board_id, params);

        board_shim.prepare_session ();
    }
}

```

(continues on next page)

(continued from previous page)

```

        board_shim.start_stream (3600);
        BoardShim.log_message (LogLevels.LEVEL_INFO.get_code (), "Start sleeping in_
↪the main thread");
        Thread.sleep (5000);
        board_shim.stop_stream ();
        System.out.println (board_shim.get_board_data_count ());
        int num_rows = BoardShim.get_num_rows (board_id);
        double[][] data = board_shim.get_current_board_data (30);
        for (int i = 0; i < num_rows; i++)
        {
            System.out.println (Arrays.toString (data[i]));
        }
        board_shim.release_session ();

        // demo for serialization
        DataFilter.write_file (data, "test.csv", "w");
        double[][] restored_data = DataFilter.read_file ("test.csv");
        System.out.println ("After Serialization:");
        for (int i = 0; i < num_rows; i++)
        {
            System.out.println (Arrays.toString (restored_data[i]));
        }
    }
}

```

5.2.3 Java Downsample Data

```

package brainflow.examples;

import java.util.Arrays;

import brainflow.AggOperations;
import brainflow.BoardIds;
import brainflow.BoardShim;
import brainflow.BrainFlowInputParams;
import brainflow.DataFilter;
import brainflow.LogLevels;

public class Downsampling
{
    public static void main (String[] args) throws Exception
    {
        // use synthetic board for demo
        BoardShim.enable_board_logger ();
        BrainFlowInputParams params = new BrainFlowInputParams ();
        int board_id = BoardIds.SYNTHETIC_BOARD.get_code ();

        BoardShim board_shim = new BoardShim (board_id, params);
        board_shim.prepare_session ();
        board_shim.start_stream (3600);
        BoardShim.log_message (LogLevels.LEVEL_INFO.get_code (), "Start sleeping in_
↪the main thread");
        Thread.sleep (5000);
    }
}

```

(continues on next page)

(continued from previous page)

```

board_shim.stop_stream ();
System.out.println (board_shim.get_board_data_count ());
double[][] data = board_shim.get_current_board_data (30);
board_shim.release_session ();

int[] eeg_channels = BoardShim.get_eeg_channels (board_id);
for (int i = 0; i < eeg_channels.length; i++)
{
    System.out.println ("Original data:");
    System.out.println (Arrays.toString (data[i]));
    // keep each second element, you can use MEAN and MEDIAN as well
    double[] downsampled_data = DataFilter.perform_downsampling (data[eeg_
↪channels[i]], 2,
        AggOperations.EACH.get_code ());
    System.out.println ("Downsampled data:");
    System.out.println (Arrays.toString (downsampled_data));
}
}
}

```

5.2.4 Java Transforms

```

package brainflow.examples;

import java.util.Arrays;

import org.apache.commons.lang3.tuple.Pair;
import org.apache.commons.math3.complex.Complex;

import brainflow.BoardIds;
import brainflow.BoardShim;
import brainflow.BrainFlowInputParams;
import brainflow.DataFilter;
import brainflow.LogLevels;
import brainflow.WindowFunctions;

public class Transforms
{
    public static void main (String[] args) throws Exception
    {
        // use synthetic board for demo
        BoardShim.enable_board_logger ();
        BrainFlowInputParams params = new BrainFlowInputParams ();
        int board_id = BoardIds.SYNTHETIC_BOARD.get_code ();
        BoardShim board_shim = new BoardShim (board_id, params);

        board_shim.prepare_session ();
        board_shim.start_stream (3600);
        BoardShim.log_message (LogLevels.LEVEL_INFO.get_code (), "Start sleeping in_
↪the main thread");
        Thread.sleep (10000);
        board_shim.stop_stream ();
        System.out.println (board_shim.get_board_data_count ());
        int num_rows = BoardShim.get_num_rows (board_id);
    }
}

```

(continues on next page)

(continued from previous page)

```

double[][] data = board_shim.get_current_board_data (64);
for (int i = 0; i < num_rows; i++)
{
    System.out.println (Arrays.toString (data[i]));
}
board_shim.release_session ();

int[] eeg_channels = BoardShim.get_eeg_channels (board_id);
for (int i = 0; i < eeg_channels.length; i++)
{
    System.out.println ("Original data:");
    System.out.println (Arrays.toString (data[eeg_channels[i]]));
    // demo for wavelet transform
    // Pair of coeffs array in format[A(J) D(J) D(J-1) ..... D(1)] where J is_
↪a
    // decomposition level, A - app coeffs, D - detailed coeffs, and array_
↪which
    // stores
    // length for each block, len of this array is decomposition_length + 1
    Pair<double[], int[]> wavelet_data = DataFilter.perform_wavelet_transform_
↪(data[eeg_channels[i]], "db4", 3);
    // print approximation coeffs
    for (int j = 0; j < wavelet_data.getRight ()[0]; j++)
    {
        System.out.print (wavelet_data.getLeft ()[j] + " ");
    }
    System.out.println ();
    // you can do smth with these coeffs here, for example denoising works via
    // thresholds for wavelet coeffs
    double[] restored_data = DataFilter.perform_inverse_wavelet_transform_
↪(wavelet_data,
        data[eeg_channels[i]].length, "db4", 3);
    System.out.println ("Restored data after wavelet:");
    System.out.println (Arrays.toString (restored_data));

    // demo for fft works only for power of 2
    // len of fft_data is N / 2 + 1
    Complex[] fft_data = DataFilter.perform_fft (data[eeg_channels[i]], 0, 64,
        WindowFunctions.NO_WINDOW.get_code ());
    double[] restored_fft_data = DataFilter.perform_ifft (fft_data);
    System.out.println ("Restored data after fft:");
    System.out.println (Arrays.toString (restored_fft_data));
}
}
}

```

5.2.5 Java Signal Filtering

```

package brainflow.examples;

import java.util.Arrays;

import brainflow.BoardIds;
import brainflow.BoardShim;
import brainflow.BrainFlowInputParams;
import brainflow.DataFilter;
import brainflow.FilterTypes;
import brainflow.LogLevels;

public class SignalFiltering
{
    public static void main (String[] args) throws Exception
    {
        // use synthetic board for demo
        BoardShim.enable_board_logger ();
        BrainFlowInputParams params = new BrainFlowInputParams ();
        int board_id = BoardIds.SYNTHETIC_BOARD.get_code ();
        BoardShim board_shim = new BoardShim (board_id, params);

        board_shim.prepare_session ();
        board_shim.start_stream (3600);
        BoardShim.log_message (LogLevels.LEVEL_INFO.get_code (), "Start sleeping in_
↪the main thread");
        Thread.sleep (5000);
        board_shim.stop_stream ();
        System.out.println (board_shim.get_board_data_count ());
        int num_rows = BoardShim.get_num_rows (board_id);
        double[][] data = board_shim.get_current_board_data (30);
        for (int i = 0; i < num_rows; i++)
        {
            System.out.println (Arrays.toString (data[i]));
        }
        board_shim.release_session ();

        int[] eeg_channels = BoardShim.get_eeg_channels (board_id);
        for (int i = 0; i < eeg_channels.length; i++)
        {
            // just for demo - apply different filters to different eeg channels
            switch (i)
            {
                case 0:
                    DataFilter.perform_lowpass (data[eeg_channels[i]], BoardShim.get_
↪sampling_rate (board_id), 20.0, 4,
                    FilterTypes.BESSEL.get_code (), 0.0);
                    break;
                case 1:
                    DataFilter.perform_highpass (data[eeg_channels[i]], BoardShim.get_
↪sampling_rate (board_id), 5.0, 4,
                    FilterTypes.BUTTERWORTH.get_code (), 0.0);
                    break;
                case 2:
                    DataFilter.perform_bandpass (data[eeg_channels[i]], BoardShim.get_
↪sampling_rate (board_id), 15.0,

```

(continues on next page)

(continued from previous page)

```

                    5.0, 4, FilterTypes.CHEBYSHEV_TYPE_1.get_code (), 1.0);
                break;
            case 3:
                DataFilter.perform_bandstop (data[eeg_channels[i]], BoardShim.get_
→sampling_rate (board_id), 50.0,
                    1.0, 4, FilterTypes.CHEBYSHEV_TYPE_1.get_code (), 1.0);
                break;
            }
        }
        System.out.println ("After signal processing:");
        for (int i = 0; i < num_rows; i++)
        {
            System.out.println (Arrays.toString (data[i]));
        }
    }
}

```

5.2.6 Java Denoising

```

package brainflow.examples;

import brainflow.Aggregations;
import brainflow.BoardIds;
import brainflow.BoardShim;
import brainflow.BrainFlowInputParams;
import brainflow.DataFilter;
import brainflow.LogLevels;

import java.util.Arrays;

public class Denoising
{
    public static void main (String[] args) throws Exception
    {
        // use synthetic board for demo
        BoardShim.enable_board_logger ();
        BrainFlowInputParams params = new BrainFlowInputParams ();
        int board_id = BoardIds.SYNTHETIC_BOARD.get_code ();
        BoardShim board_shim = new BoardShim (board_id, params);

        board_shim.prepare_session ();
        board_shim.start_stream (3600);
        BoardShim.log_message (LogLevels.LEVEL_INFO.get_code (), "Start sleeping in_
→the main thread");
        Thread.sleep (5000);
        board_shim.stop_stream ();
        System.out.println (board_shim.get_board_data_count ());
        int num_rows = BoardShim.get_num_rows (board_id);
        double[][] data = board_shim.get_current_board_data (64);
        for (int i = 0; i < num_rows; i++)
        {
            System.out.println (Arrays.toString (data[i]));
        }
        board_shim.release_session ();
    }
}

```

(continues on next page)

(continued from previous page)

```

int[] eeg_channels = BoardShim.get_eeg_channels (board_id);
for (int i = 0; i < eeg_channels.length; i++)
{
    // just for demo - apply different methods to different eeg channels
    switch (i)
    {
        // first of all you can try simple moving average or moving median
        case 0:
            DataFilter.perform_rolling_filter (data[eeg_channels[i]], 3,
↪AggOperations.MEAN.get_code ());
            break;
        case 1:
            DataFilter.perform_rolling_filter (data[eeg_channels[i]], 3,
↪AggOperations.MEDIAN.get_code ());
            break;
        // if methods above dont work good for you you should try wavelet_
↪based
        // denoising
        default:
            // try different functions and different decomposition levels here
            DataFilter.perform_wavelet_denoising (data[eeg_channels[i]], "db4
↪", 3);
            break;
    }
}
System.out.println ("After signal processing:");
for (int i = 0; i < num_rows; i++)
{
    System.out.println (Arrays.toString (data[i]));
}
}

```

5.2.7 Java Band Power

```

package brainflow.examples;

import java.util.Arrays;

import org.apache.commons.lang3.tuple.Pair;
import org.apache.commons.math3.complex.Complex;

import brainflow.BoardIds;
import brainflow.BoardShim;
import brainflow.BrainFlowInputParams;
import brainflow.DataFilter;
import brainflow.DetrendOperations;
import brainflow.LogLevels;
import brainflow.WindowFunctions;

public class BandPower
{
    public static void main (String[] args) throws Exception

```

(continues on next page)

(continued from previous page)

```

{
    // use synthetic board for demo
    BoardShim.enable_board_logger ();
    BrainFlowInputParams params = new BrainFlowInputParams ();
    int board_id = BoardIds.SYNTHETIC_BOARD.get_code ();
    BoardShim board_shim = new BoardShim (board_id, params);
    int sampling_rate = BoardShim.get_sampling_rate (board_id);
    int nfft = DataFilter.get_nearest_power_of_two (sampling_rate);

    board_shim.prepare_session ();
    board_shim.start_stream (3600);
    BoardShim.log_message (LogLevels.LEVEL_INFO.get_code (), "Start sleeping in_
↪the main thread");
    Thread.sleep (10000);
    board_shim.stop_stream ();
    double[][] data = board_shim.get_board_data ();
    board_shim.release_session ();

    int[] eeg_channels = BoardShim.get_eeg_channels (board_id);
    // seconds channel of synthetic board has big 'alpha' use it for test
    int eeg_channel = eeg_channels[1];
    // optional: detrend before psd
    DataFilter.detrend (data[eeg_channel], DetrendOperations.LINEAR.get_code ());
    Pair<double[], double[]> psd = DataFilter.get_psd_welch (data[eeg_channel],_
↪nfft, nfft / 2, sampling_rate,
        WindowFunctions.HANNING.get_code ());
    double band_power_alpha = DataFilter.get_band_power (psd, 7.0, 13.0);
    double band_power_beta = DataFilter.get_band_power (psd, 14.0, 30.0);
    System.out.println ("Alpha/Beta Ratio: " + (band_power_alpha / band_power_
↪beta));
}
}

```

5.2.8 Java EEG Metrics

```

package brainflow.examples;

import org.apache.commons.lang3.ArrayUtils;
import org.apache.commons.lang3.tuple.Pair;

import brainflow.BoardIds;
import brainflow.BoardShim;
import brainflow.BrainFlowClassifiers;
import brainflow.BrainFlowInputParams;
import brainflow.BrainFlowMetrics;
import brainflow.BrainFlowModelParams;
import brainflow.DataFilter;
import brainflow.LogLevels;
import brainflow.MLModel;

public class EEGMetrics
{
    public static void main (String[] args) throws Exception
    {

```

(continues on next page)

(continued from previous page)

```

BoardShim.enable_board_logger ();
BrainFlowInputParams params = new BrainFlowInputParams ();
int board_id = parse_args (args, params);
BoardShim board_shim = new BoardShim (board_id, params);
int master_board_id = board_shim.get_board_id ();
int sampling_rate = BoardShim.get_sampling_rate (master_board_id);
int[] eeg_channels = BoardShim.get_eeg_channels (master_board_id);

board_shim.prepare_session ();
board_shim.start_stream (3600);
BoardShim.log_message (LogLevels.LEVEL_INFO.get_code (), "Start sleeping in_
↳the main thread");
    // recommended window size for eeg metric calculation is at least 4 seconds,
    // bigger is better
    Thread.sleep (5000);
board_shim.stop_stream ();
double[][] data = board_shim.get_board_data ();
board_shim.release_session ();

Pair<double[], double[]> bands = DataFilter.get_avg_band_powers (data, eeg_
↳channels, sampling_rate, true);
double[] feature_vector = ArrayUtils.addAll (bands.getLeft (), bands.getRight_
↳());
BrainFlowModelParams model_params = new BrainFlowModelParams_
↳(BrainFlowMetrics.CONCENTRATION.get_code (),
    BrainFlowClassifiers.REGRESSION.get_code ());
MLModel concentration = new MLModel (model_params);
concentration.prepare ();
System.out.print ("Concentration: " + concentration.predict (feature_vector));
concentration.release ();
}

private static int parse_args (String[] args, BrainFlowInputParams params)
{
    int board_id = -1;
    for (int i = 0; i < args.length; i++)
    {
        if (args[i].equals ("--ip-address"))
        {
            params.ip_address = args[i + 1];
        }
        if (args[i].equals ("--serial-port"))
        {
            params.serial_port = args[i + 1];
        }
        if (args[i].equals ("--ip-port"))
        {
            params.ip_port = Integer.parseInt (args[i + 1]);
        }
        if (args[i].equals ("--ip-protocol"))
        {
            params.ip_protocol = Integer.parseInt (args[i + 1]);
        }
        if (args[i].equals ("--other-info"))
        {
            params.other_info = args[i + 1];
        }
    }
}

```

(continues on next page)

(continued from previous page)

```

        if (args[i].equals ("--board-id"))
        {
            board_id = Integer.parseInt (args[i + 1]);
        }
        if (args[i].equals ("--timeout"))
        {
            params.timeout = Integer.parseInt (args[i + 1]);
        }
        if (args[i].equals ("--serial-number"))
        {
            params.serial_number = args[i + 1];
        }
        if (args[i].equals ("--file"))
        {
            params.file = args[i + 1];
        }
    }
    return board_id;
}
}

```

5.3 C#

5.3.1 C# Read Data from a Board

```

using System;
using brainflow;

using Accord.Math;

namespace test
{
    class GetBoardData
    {
        static void Main (string[] args)
        {
            BoardShim.enable_dev_board_logger ();

            BrainFlowInputParams input_params = new BrainFlowInputParams ();
            int board_id = parse_args (args, input_params);

            BoardShim board_shim = new BoardShim (board_id, input_params);
            board_shim.prepare_session ();
            // board_shim.start_stream (); // use this for default options
            board_shim.start_stream (450000, "file://file_stream.csv:w");
            System.Threading.Thread.Sleep (5000);
            board_shim.stop_stream ();
            double[,] unprocessed_data = board_shim.get_current_board_data (20);
            int[] eeg_channels = BoardShim.get_eeg_channels (board_id);
            foreach (var index in eeg_channels)
                Console.WriteLine ("[{0}]", string.Join (" ", unprocessed_data.
↪GetRow (index)));
            board_shim.release_session ();
        }
    }
}

```

(continues on next page)

(continued from previous page)

```

    }

    static int parse_args (string[] args, BrainFlowInputParams input_params)
    {
        int board_id = (int)BoardIds.SYNTHETIC_BOARD; //assume synthetic board by_
↪default
        // use docs to get params for your specific board, e.g. set serial_port_
↪for Cyton
        for (int i = 0; i < args.Length; i++)
        {
            if (args[i].Equals ("--ip-address"))
            {
                input_params.ip_address = args[i + 1];
            }
            if (args[i].Equals ("--mac-address"))
            {
                input_params.mac_address = args[i + 1];
            }
            if (args[i].Equals ("--serial-port"))
            {
                input_params.serial_port = args[i + 1];
            }
            if (args[i].Equals ("--other-info"))
            {
                input_params.other_info = args[i + 1];
            }
            if (args[i].Equals ("--ip-port"))
            {
                input_params.ip_port = Convert.ToInt32 (args[i + 1]);
            }
            if (args[i].Equals ("--ip-protocol"))
            {
                input_params.ip_protocol = Convert.ToInt32 (args[i + 1]);
            }
            if (args[i].Equals ("--board-id"))
            {
                board_id = Convert.ToInt32 (args[i + 1]);
            }
            if (args[i].Equals ("--timeout"))
            {
                input_params.timeout = Convert.ToInt32(args[i + 1]);
            }
            if (args[i].Equals ("--serial-number"))
            {
                input_params.serial_number = args[i + 1];
            }
            if (args[i].Equals ("--file"))
            {
                input_params.file = args[i + 1];
            }
        }
        return board_id;
    }
}

```

5.3.2 C# Read Write File

```

using System;
using brainflow;

using Accord.Math;

namespace test
{
    class Serialization
    {
        static void Main (string[] args)
        {
            // use synthetic board for demo
            BoardShim.enable_dev_board_logger ();
            BrainFlowInputParams input_params = new BrainFlowInputParams ();
            int board_id = (int)BoardIds.SYNTHETIC_BOARD;

            BoardShim board_shim = new BoardShim (board_id, input_params);
            board_shim.prepare_session ();
            board_shim.start_stream (3600);
            System.Threading.Thread.Sleep (5000);
            board_shim.stop_stream ();
            double[,] unprocessed_data = board_shim.get_current_board_data (20);
            int[] eeg_channels = BoardShim.get_eeg_channels (board_id);
            Console.WriteLine ("Before serialization:");
            foreach (var index in eeg_channels)
                Console.WriteLine ("[{0}]", string.Join (" ", unprocessed_data.
↪GetRow (index)));
            board_shim.release_session ();

            // demo for data serialization
            DataFilter.write_file (unprocessed_data, "test.csv", "w");
            double[,] restored_data = DataFilter.read_file ("test.csv");
            Console.WriteLine ("After serialization:");
            foreach (var index in eeg_channels)
                Console.WriteLine ("[{0}]", string.Join (" ", restored_data.GetRow_
↪(index)));
        }
    }
}

```

5.3.3 C# Downsample Data

```

using System;
using brainflow;

using Accord.Math;

namespace test
{
    class Downsampling
    {
        static void Main (string[] args)
        {

```

(continues on next page)

(continued from previous page)

```

// use synthetic board for demo
BoardShim.enable_dev_board_logger ();
BrainFlowInputParams input_params = new BrainFlowInputParams ();
int board_id = (int)BoardIds.SYNTHETIC_BOARD;

BoardShim board_shim = new BoardShim (board_id, input_params);
board_shim.prepare_session ();
board_shim.start_stream (3600);
System.Threading.Thread.Sleep (5000);
board_shim.stop_stream ();
double[,] unprocessed_data = board_shim.get_current_board_data (20);
int[] eeg_channels = BoardShim.get_eeg_channels (board_id);
board_shim.release_session ();

for (int i = 0; i < eeg_channels.Length; i++)
{
    Console.WriteLine ("Before processing:");
    Console.WriteLine ("[{0}]", string.Join (" ", unprocessed_data.
↪GetRow(eeg_channels[i])));
    // you can use MEAN, MEDIAN or EACH for downsampling
    double[] filtered = DataFilter.perform_downsampling (unprocessed_data.
↪GetRow (eeg_channels[i]), 3, (int)AggOperations.MEDIAN);
    Console.WriteLine ("Before processing:");
    Console.WriteLine ("[{0}]", string.Join (" ", filtered));
}
}
}

```

5.3.4 C# Transforms

```

using System;
using System.Numerics;
using brainflow;

using Accord.Math;

namespace test
{
    class Transforms
    {
        static void Main (string[] args)
        {
            // use synthetic board for demo
            BoardShim.enable_dev_board_logger ();
            BrainFlowInputParams input_params = new BrainFlowInputParams ();
            int board_id = (int)BoardIds.SYNTHETIC_BOARD;

            BoardShim board_shim = new BoardShim (board_id, input_params);
            board_shim.prepare_session ();
            board_shim.start_stream (3600);
            System.Threading.Thread.Sleep (5000);
            board_shim.stop_stream ();
            double[,] unprocessed_data = board_shim.get_current_board_data (64);
            int[] eeg_channels = BoardShim.get_eeg_channels (board_id);

```

(continues on next page)

(continued from previous page)

```

        board_shim.release_session ();

        for (int i = 0; i < eeg_channels.Length; i++)
        {
            Console.WriteLine ("Original data:");
            Console.WriteLine ("[{0}]", string.Join (" ", unprocessed_data.
↪GetRow (eeg_channels[i])));
            // demo for wavelet transform
            // tuple of coeffs array in format[A(J) D(J) D(J-1) ..... D(1)] where_
↪J is a
            // decomposition level, A - app coeffs, D - detailed coeffs, and_
↪array which stores
            // length for each block, len of this array is decomposition_length +_
↪1
            Tuple<double[], int[]> wavelet_data = DataFilter.perform_wavelet_
↪transform(unprocessed_data.GetRow (eeg_channels[i]), "db4", 3);
            // print app coeffs
            for (int j = 0; j < wavelet_data.Item2[0]; j++)
            {
                Console.Write (wavelet_data.Item1[j] + " ");
            }
            Console.WriteLine ();
            // you can do smth with wavelet coeffs here, for example denoising_
↪works via thresholds for wavelets coeffs
            double[] restored_data = DataFilter.perform_inverse_wavelet_transform_
↪(wavelet_data, unprocessed_data.GetRow (eeg_channels[i]).Length, "db4", 3);
            Console.WriteLine ("Restored wavelet data:");
            Console.WriteLine ("[{0}]", string.Join (" ", restored_data));

            // demo for fft
            // end_pos - start_pos must be a power of 2
            Complex[] fft_data = DataFilter.perform_fft (unprocessed_data.GetRow_
↪(eeg_channels[i]), 0, 64, (int)WindowFunctions.HAMMING);
            // len of fft_data is N / 2 + 1
            double[] restored_fft_data = DataFilter.perform_ifft (fft_data);
            Console.WriteLine ("Restored fft data:");
            Console.WriteLine ("[{0}]", string.Join (" ", restored_fft_data));
        }
    }
}

```

5.3.5 C# Signal Filtering

```

using System;
using brainflow;

using Accord.Math;

namespace test
{
    class SignalFiltering
    {
        static void Main (string[] args)
        {

```

(continues on next page)

(continued from previous page)

```

// use synthetic board for demo
BoardShim.enable_dev_board_logger ();
BrainFlowInputParams input_params = new BrainFlowInputParams ();
int board_id = (int)BoardIds.SYNTHETIC_BOARD;

BoardShim board_shim = new BoardShim (board_id, input_params);
board_shim.prepare_session ();
board_shim.start_stream (3600);
System.Threading.Thread.Sleep (5000);
board_shim.stop_stream ();
double[,] unprocessed_data = board_shim.get_current_board_data (20);
int[] eeg_channels = BoardShim.get_eeg_channels (board_id);
board_shim.release_session ();

// for demo apply different filters to different channels
double[] filtered;
for (int i = 0; i < eeg_channels.Length; i++)
{
    Console.WriteLine ("Before processing:");
    Console.WriteLine ("[{0}]", string.Join (" ", unprocessed_data.
↪GetRow (eeg_channels[i])));
    switch (i)
    {
        case 0:
            filtered = DataFilter.perform_lowpass (unprocessed_data.
↪GetRow(eeg_channels[i]), BoardShim.get_sampling_rate (board_id), 20.0, 4,
↪(int)FilterTypes.BESSEL, 0.0);
            Console.WriteLine ("Filtered channel " + eeg_channels[i]);
            Console.WriteLine ("[{0}]", string.Join (" ", filtered));
            break;
        case 1:
            filtered = DataFilter.perform_highpass (unprocessed_data.
↪GetRow (eeg_channels[i]), BoardShim.get_sampling_rate (board_id), 2.0, 4,
↪(int)FilterTypes.BUTTERWORTH, 0.0);
            Console.WriteLine ("Filtered channel " + eeg_channels[i]);
            Console.WriteLine ("[{0}]", string.Join (" ", filtered));
            break;
        case 2:
            filtered = DataFilter.perform_bandpass (unprocessed_data.
↪GetRow (eeg_channels[i]), BoardShim.get_sampling_rate (board_id), 15.0, 5.0, 2,
↪(int)FilterTypes.BUTTERWORTH, 0.0);
            Console.WriteLine ("Filtered channel " + eeg_channels[i]);
            Console.WriteLine ("[{0}]", string.Join (" ", filtered));
            break;
        case 3:
            filtered = DataFilter.perform_bandstop (unprocessed_data.
↪GetRow (eeg_channels[i]), BoardShim.get_sampling_rate (board_id), 50.0, 1.0, 6,
↪(int)FilterTypes.CHEBYSHEV_TYPE_1, 1.0);
            Console.WriteLine ("Filtered channel " + eeg_channels[i]);
            Console.WriteLine ("[{0}]", string.Join (" ", filtered));
            break;
    }
}
}
}
}
}

```

5.3.6 C# Denoising

```

using System;
using brainflow;

using Accord.Math;

namespace test
{
    class Denoising
    {
        static void Main (string[] args)
        {
            // use synthetic board for demo
            BoardShim.enable_dev_board_logger ();
            BrainFlowInputParams input_params = new BrainFlowInputParams ();
            int board_id = (int)BoardIds.SYNTHETIC_BOARD;

            BoardShim board_shim = new BoardShim (board_id, input_params);
            board_shim.prepare_session ();
            board_shim.start_stream (3600);
            System.Threading.Thread.Sleep (5000);
            board_shim.stop_stream ();
            double[,] unprocessed_data = board_shim.get_current_board_data (64);
            int[] eeg_channels = BoardShim.get_eeg_channels (board_id);
            foreach (var index in eeg_channels)
                Console.WriteLine ("[{0}]", string.Join (" ", unprocessed_data.
↪GetRow (index)));
            board_shim.release_session ();

            // for demo apply different methods to different channels
            double[] filtered;
            for (int i = 0; i < eeg_channels.Length; i++)
            {
                switch (i)
                {
                    // first of all you can try simple moving average or moving median
                    case 0:
                        filtered = DataFilter.perform_rolling_filter (unprocessed_
↪data.GetRow (eeg_channels[i]), 3, (int)AggOperations.MEAN);
                        Console.WriteLine ("Filtered channel " + eeg_channels[i]);
                        Console.WriteLine ("[{0}]", string.Join (" ", filtered));
                        break;
                    case 1:
                        filtered = DataFilter.perform_rolling_filter (unprocessed_
↪data.GetRow (eeg_channels[i]), 3, (int)AggOperations.MEDIAN);
                        Console.WriteLine ("Filtered channel " + eeg_channels[i]);
                        Console.WriteLine ("[{0}]", string.Join (" ", filtered));
                        break;
                    // if for your signal these methods dont work good you can try_
↪wavelet based denoising
                    default:
                        // feel free to try different functions and different_
↪decomposition levels
                        filtered = DataFilter.perform_wavelet_denoising (unprocessed_
↪data.GetRow (eeg_channels[i]), "db4", 3);
                        Console.WriteLine ("Filtered channel " + eeg_channels[i]);
                }
            }
        }
    }
}

```

(continues on next page)

(continued from previous page)

```

        Console.WriteLine ("[{0}]", string.Join (" ", filtered));
        break;
    }
}
}
}
}
}

```

5.3.7 C# Band Power

```

using System;
using System.Numerics;
using brainflow;

using Accord.Math;

namespace test
{
    class BandPower
    {
        static void Main (string[] args)
        {
            // use synthetic board for demo
            BoardShim.enable_dev_board_logger ();
            BrainFlowInputParams input_params = new BrainFlowInputParams ();
            int board_id = (int)BoardIds.SYNTHETIC_BOARD;
            int sampling_rate = BoardShim.get_sampling_rate (board_id);
            int nfft = DataFilter.get_nearest_power_of_two(sampling_rate);

            BoardShim board_shim = new BoardShim (board_id, input_params);
            board_shim.prepare_session ();
            board_shim.start_stream (3600);
            System.Threading.Thread.Sleep (10000);
            board_shim.stop_stream ();
            double[,] data = board_shim.get_board_data ();
            int[] eeg_channels = BoardShim.get_eeg_channels (board_id);
            // use second channel of synthetic board to see 'alpha'
            int channel = eeg_channels[1];
            board_shim.release_session ();
            double[] detrend = DataFilter.detrend(data.GetRow(channel),
↪(int)DetrendOperations.LINEAR);
            Tuple<double[], double[]> psd = DataFilter.get_psd_welch (detrend, nfft,
↪nfft / 2, sampling_rate, (int)WindowFunctions.HANNING);
            double band_power_alpha = DataFilter.get_band_power (psd, 7.0, 13.0);
            double band_power_beta = DataFilter.get_band_power (psd, 14.0, 30.0);
            Console.WriteLine ("Alpha/Beta Ratio:" + (band_power_alpha/ band_power_
↪beta));
        }
    }
}

```

5.3.8 C# EEG Metrics

```

using System;
using System.Numerics;
using brainflow;

using Accord.Math;

namespace test
{
    class EEGMetrics
    {
        static void Main (string[] args)
        {
            // use synthetic board for demo
            BoardShim.enable_dev_board_logger ();
            BrainFlowInputParams input_params = new BrainFlowInputParams ();
            int board_id = parse_args (args, input_params);
            BoardShim board_shim = new BoardShim (board_id, input_params);
            int sampling_rate = BoardShim.get_sampling_rate (board_shim.get_board_id_
↪());

            int nfft = DataFilter.get_nearest_power_of_two (sampling_rate);
            int[] eeg_channels = BoardShim.get_eeg_channels (board_shim.get_board_id_
↪());

            board_shim.prepare_session ();
            board_shim.start_stream (3600);
            System.Threading.Thread.Sleep (10000);
            board_shim.stop_stream ();
            double[,] data = board_shim.get_board_data ();
            board_shim.release_session ();

            Tuple<double[], double[]> bands = DataFilter.get_avg_band_powers (data,
↪eeg_channels, sampling_rate, true);
            double[] feature_vector = bands.Item1.Concatenate (bands.Item2);
            BrainFlowModelParams model_params = new BrainFlowModelParams_
↪((int)BrainFlowMetrics.CONCENTRATION, (int)BrainFlowClassifiers.REGRESSION);
            MLModel concentration = new MLModel (model_params);
            concentration.prepare ();
            Console.WriteLine ("Concentration: " + concentration.predict (feature_
↪vector));
            concentration.release ();
        }

        static int parse_args (string[] args, BrainFlowInputParams input_params)
        {
            int board_id = (int)BoardIds.SYNTHETIC_BOARD; //assume synthetic board by_
↪default
            // use docs to get params for your specific board, e.g. set serial_port_
↪for Cyton
            for (int i = 0; i < args.Length; i++)
            {
                if (args[i].Equals ("--ip-address"))
                {
                    input_params.ip_address = args[i + 1];
                }
                if (args[i].Equals ("--mac-address"))

```

(continues on next page)

(continued from previous page)

```

        {
            input_params.mac_address = args[i + 1];
        }
        if (args[i].Equals ("--serial-port"))
        {
            input_params.serial_port = args[i + 1];
        }
        if (args[i].Equals ("--other-info"))
        {
            input_params.other_info = args[i + 1];
        }
        if (args[i].Equals ("--ip-port"))
        {
            input_params.ip_port = Convert.ToInt32 (args[i + 1]);
        }
        if (args[i].Equals ("--ip-protocol"))
        {
            input_params.ip_protocol = Convert.ToInt32 (args[i + 1]);
        }
        if (args[i].Equals ("--board-id"))
        {
            board_id = Convert.ToInt32 (args[i + 1]);
        }
        if (args[i].Equals ("--timeout"))
        {
            input_params.timeout = Convert.ToInt32 (args[i + 1]);
        }
        if (args[i].Equals ("--serial-number"))
        {
            input_params.serial_number = args[i + 1];
        }
        if (args[i].Equals ("--file"))
        {
            input_params.file = args[i + 1];
        }
    }
    return board_id;
}
}
}

```

5.4 C++

To compile examples below for Linux or MacOS run:

```

cd tests/cpp/get_data_demo
mkdir build
cd build
cmake -DCMAKE_PREFIX_PATH=TYPE_FULL_PATH_TO_BRAINFLOW_INSTALLED_FOLDER ..
# e.g. cmake -DCMAKE_PREFIX_PATH=/home/andrey/brainflow/installed_linux ..
make

```

For Windows it's almost the same.

Make sure that compiled dynamic libraries exist in search path before running an executable by doing one of

the following:

- for Linux and MacOS add them to LD_LIBRARY_PATH env variable
- for Windows add them to PATH env variable
- or just copypaste them to the folder where your executable is located

5.4.1 CMake File Example

```
cmake_minimum_required (VERSION 3.10)
project (BRAINFLOW_GET_DATA)

set (CMAKE_CXX_STANDARD 11)
set (CMAKE_VERBOSE_MAKEFILE ON)

macro (configure_msvc_runtime)
    if (MSVC)
        # Default to statically-linked runtime.
        if ("${MSVC_RUNTIME}" STREQUAL "")
            set (MSVC_RUNTIME "static")
        endif ()
        # Set compiler options.
        set (variables
            CMAKE_C_FLAGS_DEBUG
            CMAKE_C_FLAGS_MINSIZEREL
            CMAKE_C_FLAGS_RELEASE
            CMAKE_C_FLAGS_RELWITHDEBINFO
            CMAKE_CXX_FLAGS_DEBUG
            CMAKE_CXX_FLAGS_MINSIZEREL
            CMAKE_CXX_FLAGS_RELEASE
            CMAKE_CXX_FLAGS_RELWITHDEBINFO
        )
        if ("${MSVC_RUNTIME}" STREQUAL "static")
            message (STATUS
                "MSVC -> forcing use of statically-linked runtime."
            )
            foreach (variable ${variables})
                if (${variable} MATCHES "/MD")
                    string (REGEX REPLACE "/MD" "/MT" ${variable} "${${variable}}")
                endif ()
            endforeach ()
        else ()
            message (STATUS
                "MSVC -> forcing use of dynamically-linked runtime."
            )
            foreach (variable ${variables})
                if (${variable} MATCHES "/MT")
                    string (REGEX REPLACE "/MT" "/MD" ${variable} "${${variable}}")
                endif ()
            endforeach ()
        endif ()
    endif ()
endmacro ()

# link msvc runtime statically
configure_msvc_runtime()
```

(continues on next page)

(continued from previous page)

```

find_package (
    brainflow CONFIG REQUIRED
)

add_executable (
    brainflow_get_data
    src/brainflow_get_data.cpp
)

target_include_directories (
    brainflow_get_data PUBLIC
    ${brainflow_INCLUDE_DIRS}
)

target_link_libraries (
    brainflow_get_data PUBLIC
    # for some systems(ubuntu for example) order matters
    ${BrainflowPath}
    ${MLModulePath}
    ${DataHandlerPath}
    ${BoardControllerPath}
)

```

5.4.2 C++ Read Data from a Board

```

#include <iostream>
#include <stdlib.h>
#include <string>

#ifdef _WIN32
#include <windows.h>
#else
#include <unistd.h>
#endif

#include "board_shim.h"

using namespace std;

void print_head (double **data_buf, int num_channels, int num_data_points);
bool parse_args (int argc, char *argv[], struct BrainFlowInputParams *params, int_
↳ *board_id);

int main (int argc, char *argv[])
{
    struct BrainFlowInputParams params;
    int board_id = 0;
    if (!parse_args (argc, argv, &params, &board_id))
    {
        return -1;
    }

    BoardShim::enable_dev_board_logger ();

```

(continues on next page)

(continued from previous page)

```

BoardShim *board = new BoardShim (board_id, params);
double **data = NULL;
int res = 0;
int num_rows = 0;

try
{
    board->prepare_session ();
    board->start_stream ();
    // board->start_stream (45000, (char *) "file://file_stream_test.csv:a"); //
↪store data in a
    // file directly during streaming
    BoardShim::log_message ((int)LogLevels::LEVEL_INFO, "Start sleeping in the
↪main thread");
#ifdef _WIN32
    Sleep (5000);
#else
    sleep (5);
#endif

    board->stop_stream ();
    int data_count = 0;
    data = board->get_board_data (&data_count);
    BoardShim::log_message ((int)LogLevels::LEVEL_INFO, "read %d packages", data_
↪count);
    board->release_session ();
    // for STREAMING_BOARD and PLAYBACK_FILE_BOARD you have to query information
↪using board id
    // for master board because for STREAMING_BOARD data format is determined by
↪master board!
    if ((board_id == (int)BoardIds::STREAMING_BOARD) ||
        (board_id == (int)BoardIds::PLAYBACK_FILE_BOARD))
    {
        board_id = std::stoi (params.other_info);
        BoardShim::log_message ((int)LogLevels::LEVEL_INFO, "Use Board Id %d",
↪board_id);
    }
    num_rows = BoardShim::get_num_rows (board_id);
    std::cout << std::endl << "Data from the board" << std::endl << std::endl;
    print_head (data, num_rows, data_count);
}
catch (const BrainFlowException &err)
{
    BoardShim::log_message ((int)LogLevels::LEVEL_ERROR, err.what ());
    res = err.exit_code;
    if (board->is_prepared ())
    {
        board->release_session ();
    }
}

if (data != NULL)
{
    for (int i = 0; i < num_rows; i++)
    {
        delete[] data[i];
    }
}

```

(continues on next page)

(continued from previous page)

```

    }
}
delete[] data;
delete board;

return res;
}

void print_head (double **data_buf, int num_channels, int num_data_points)
{
    std::cout << "Total Channels for this board: " << num_channels << std::endl;
    int num_points = (num_data_points < 5) ? num_data_points : 5;
    for (int i = 0; i < num_channels; i++)
    {
        std::cout << "Channel " << i << ": ";
        for (int j = 0; j < num_points; j++)
        {
            std::cout << data_buf[i][j] << ",";
        }
        std::cout << std::endl;
    }
}

bool parse_args (int argc, char *argv[], struct BrainFlowInputParams *params, int_
↳ *board_id)
{
    bool board_id_found = false;
    for (int i = 1; i < argc; i++)
    {
        if (std::string (argv[i]) == std::string ("--board-id"))
        {
            if (i + 1 < argc)
            {
                i++;
                board_id_found = true;
                *board_id = std::stoi (std::string (argv[i]));
            }
            else
            {
                std::cerr << "missed argument" << std::endl;
                return false;
            }
        }
        if (std::string (argv[i]) == std::string ("--ip-address"))
        {
            if (i + 1 < argc)
            {
                i++;
                params->ip_address = std::string (argv[i]);
            }
            else
            {
                std::cerr << "missed argument" << std::endl;
                return false;
            }
        }
        if (std::string (argv[i]) == std::string ("--ip-port"))

```

(continues on next page)

(continued from previous page)

```

{
    if (i + 1 < argc)
    {
        i++;
        params->ip_port = std::stoi (std::string (argv[i]));
    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
if (std::string (argv[i]) == std::string ("--serial-port"))
{
    if (i + 1 < argc)
    {
        i++;
        params->serial_port = std::string (argv[i]);
    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
if (std::string (argv[i]) == std::string ("--ip-protocol"))
{
    if (i + 1 < argc)
    {
        i++;
        params->ip_protocol = std::stoi (std::string (argv[i]));
    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
if (std::string (argv[i]) == std::string ("--timeout"))
{
    if (i + 1 < argc)
    {
        i++;
        params->timeout = std::stoi (std::string (argv[i]));
    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
if (std::string (argv[i]) == std::string ("--other-info"))
{
    if (i + 1 < argc)
    {
        i++;
        params->other_info = std::string (argv[i]);
    }
}

```

(continues on next page)

(continued from previous page)

```

    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
if (std::string (argv[i]) == std::string ("--mac-address"))
{
    if (i + 1 < argc)
    {
        i++;
        params->mac_address = std::string (argv[i]);
    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
if (std::string (argv[i]) == std::string ("--serial-number"))
{
    if (i + 1 < argc)
    {
        i++;
        params->serial_number = std::string (argv[i]);
    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
if (std::string (argv[i]) == std::string ("--file"))
{
    if (i + 1 < argc)
    {
        i++;
        params->file = std::string (argv[i]);
    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
}
if (!board_id_found)
{
    std::cerr << "board id is not provided" << std::endl;
    return false;
}
return true;
}

```

5.4.3 C++ Read Write File

```

#include <iostream>
#include <stdlib.h>

#ifdef _WIN32
#include <windows.h>
#else
#include <unistd.h>
#endif

#include "board_shim.h"
#include "data_filter.h"

using namespace std;

void print_head (double **data_buf, int num_channels, int num_data_points);

int main (int argc, char *argv[])
{
    struct BrainFlowInputParams params;
    // use synthetic board for demo
    int board_id = (int)BoardIds::SYNTHETIC_BOARD;

    BoardShim::enable_dev_board_logger ();

    BoardShim *board = new BoardShim (board_id, params);
    double **data = NULL;
    int res = 0;
    int num_rows = 0;

    try
    {
        board->prepare_session ();
        board->start_stream ();
        BoardShim::log_message ((int)LogLevels::LEVEL_INFO, "Start sleeping in the_
↪main thread");
#ifdef _WIN32
        Sleep (5000);
#else
        sleep (5);
#endif

        board->stop_stream ();
        int data_count = 0;
        data = board->get_board_data (&data_count);
        BoardShim::log_message ((int)LogLevels::LEVEL_INFO, "read %d packages", data_
↪count);
        board->release_session ();
        num_rows = BoardShim::get_num_rows (board_id);
        std::cout << std::endl << "Data from the board" << std::endl << std::endl;
        print_head (data, num_rows, data_count);

        // demo for serialization
        DataFilter::write_file (
            data, num_rows, data_count, "test.csv", "w"); // use "a" for append mode
        int restored_num_rows = 0;
    }
}

```

(continues on next page)

(continued from previous page)

```

    int restored_num_cols = 0;
    double **restored_data =
        DataFilter::read_file (&restored_num_rows, &restored_num_cols, "test.csv
↪");
    std::cout << std::endl
        << "Data from the file, num packages is " << restored_num_cols <<_
↪std::endl
        << std::endl;
    print_head (restored_data, restored_num_rows, restored_num_cols);
    for (int i = 0; i < restored_num_rows; i++)
    {
        delete[] restored_data[i];
    }
    delete[] restored_data;
}
catch (const BrainFlowException &err)
{
    BoardShim::log_message ((int)LogLevels::LEVEL_ERROR, err.what ());
    res = err.exit_code;
}

if (data != NULL)
{
    for (int i = 0; i < num_rows; i++)
    {
        delete[] data[i];
    }
}
delete[] data;
delete board;

return res;
}

void print_head (double **data_buf, int num_channels, int num_data_points)
{
    std::cout << "Total Channels for this board: " << num_channels << std::endl;
    int num_points = (num_data_points < 5) ? num_data_points : 5;
    for (int i = 0; i < num_channels; i++)
    {
        std::cout << "Channel " << i << ": ";
        for (int j = 0; j < num_points; j++)
        {
            std::cout << data_buf[i][j] << ", ";
        }
        std::cout << std::endl;
    }
}

```

5.4.4 C++ Downsample Data

```

#include <iostream>
#include <stdlib.h>

#ifdef _WIN32
#include <windows.h>
#else
#include <unistd.h>
#endif

#include "board_shim.h"
#include "data_filter.h"

using namespace std;

void print_one_row (double *data, int num_data_points);

int main (int argc, char *argv[])
{
    struct BrainFlowInputParams params;
    // use synthetic board for demo
    int board_id = (int)BoardIds::SYNTHETIC_BOARD;

    BoardShim::enable_dev_board_logger ();

    BoardShim *board = new BoardShim (board_id, params);
    double **data = NULL;
    int *eeg_channels = NULL;
    int num_rows = 0;
    int res = 0;

    try
    {
        board->prepare_session ();
        board->start_stream ();
        BoardShim::log_message ((int)LogLevels::LEVEL_INFO, "Start sleeping in the_
↪main thread");
#ifdef _WIN32
        Sleep (5000);
#else
        sleep (5);
#endif

        board->stop_stream ();
        int data_count = 0;
        data = board->get_board_data (&data_count);
        BoardShim::log_message ((int)LogLevels::LEVEL_INFO, "read %d packages", data_
↪count);
        board->release_session ();
        num_rows = BoardShim::get_num_rows (board_id);

        // downsample only eeg channels and print them
        int eeg_num_channels = 0;
        eeg_channels = BoardShim::get_eeg_channels (board_id, &eeg_num_channels);
        double *downsampled_data = NULL;
        int filtered_size = 0;

```

(continues on next page)

(continued from previous page)

```

    for (int i = 0; i < eeg_num_channels; i++)
    {
        std::cout << "Data from :" << eeg_channels[i] << " before downsampling " <
        << std::endl;
        print_one_row (data[eeg_channels[i]], data_count);

        // just for demo apply different downsampling algorithms to different_
        // channels
        // downsampling here doesnt apply lowpass filter for you, it just_
        // aggregates data points
        switch (i)
        {
            case 0:
                downsampled_data = DataFilter::perform_downsampling (data[eeg_
                channels[i]],
                    data_count, 2, (int)AggOperations::MEAN, &filtered_size);
                break;
            case 1:
                downsampled_data = DataFilter::perform_downsampling (data[eeg_
                channels[i]],
                    data_count, 3, (int)AggOperations::MEDIAN, &filtered_size);
                break;
            default:
                downsampled_data = DataFilter::perform_downsampling (data[eeg_
                channels[i]],
                    data_count, 2, (int)AggOperations::EACH, &filtered_size);
                break;
        }

        std::cout << "Data from :" << eeg_channels[i] << " after downsampling " <
        << std::endl;
        print_one_row (downsampled_data, filtered_size);
        delete[] downsampled_data;
    }
}

catch (const BrainFlowException &err)
{
    BoardShim::log_message ((int)LogLevels::LEVEL_ERROR, err.what ());
    res = err.exit_code;
}

if (data != NULL)
{
    for (int i = 0; i < num_rows; i++)
    {
        delete[] data[i];
    }
}

delete[] data;
delete[] eeg_channels;
delete board;

return res;
}

void print_one_row (double *data, int num_data_points)
{

```

(continues on next page)

(continued from previous page)

```

// print only first 10 data points
int num_points = (num_data_points < 10) ? num_data_points : 10;
for (int i = 0; i < num_points; i++)
{
    std::cout << data[i] << " ";
}
std::cout << std::endl;
}

```

5.4.5 C++ Transforms

```

#include <iostream>
#include <stdlib.h>

#ifdef _WIN32
#include <windows.h>
#else
#include <unistd.h>
#endif

#include "board_shim.h"
#include "data_filter.h"

using namespace std;

void print_one_row (double *data, int num_data_points);

int main (int argc, char *argv[])
{
    struct BrainFlowInputParams params;
    // use synthetic board for demo
    int board_id = (int)BoardIds::SYNTHETIC_BOARD;

    BoardShim::enable_dev_board_logger ();

    BoardShim *board = new BoardShim (board_id, params);
    double **data = NULL;
    int *eeg_channels = NULL;
    int num_rows = 0;
    int res = 0;

    try
    {
        board->prepare_session ();
        board->start_stream ();
        BoardShim::log_message ((int)LogLevels::LEVEL_INFO, "Start sleeping in the_
↪main thread");
#ifdef _WIN32
        Sleep (10000);
#else
        sleep (10);
#endif

        board->stop_stream ();
        int data_count = 0;
    }
}

```

(continues on next page)

(continued from previous page)

```

data = board->get_current_board_data (128, &data_count);
if (data_count != 128)
{
    BoardShim::log_message ((int)LogLevel::LEVEL_ERROR,
        "read %d packages, for this test we want exactly 128 packages", data_
↪count);
    return (int)BrainFlowExitCodes::GENERAL_ERROR;
}
board->release_session ();
num_rows = BoardShim::get_num_rows (board_id);

int eeg_num_channels = 0;
eeg_channels = BoardShim::get_eeg_channels (board_id, &eeg_num_channels);
for (int i = 0; i < eeg_num_channels; i++)
{
    // demo for wavelet transform
    // std::pair of coeffs array in format[A(J) D(J) D(J-1) ..... D(1)] where_
↪J is a
    // decomposition level, A - app coeffs, D - detailed coeffs, and array_
↪which stores
    // length for each block, len of this array is decomposition_length + 1
    std::pair<double *, int *> wavelet_output =
        DataFilter::perform_wavelet_transform (data[eeg_channels[i]], data_
↪count, "db4", 4);

    // you can do smth with wavelet coeffs here, for example denoising works_
↪via thresholds
    // for wavelet coefficients
    std::cout << "approximation coefficients:" << std::endl;
    for (int i = 0; i < wavelet_output.second[0]; i++)
    {
        std::cout << wavelet_output.first[i] << " ";
    }
    std::cout << std::endl;
    std::cout << "first block of detailed coefficients:" << std::endl;
    for (int i = wavelet_output.second[0];
        i < wavelet_output.second[0] + wavelet_output.second[1]; i++)
    {
        std::cout << wavelet_output.first[i] << " ";
    }
    std::cout << std::endl;

    double *restored_data = DataFilter::perform_inverse_wavelet_transform (
        wavelet_output, data_count, "db4", 4);

    std::cout << "Original data:" << std::endl;
    print_one_row (data[eeg_channels[i]], data_count);
    std::cout << "Restored after inverse wavelet transform data:" << std::
↪endl;
    print_one_row (restored_data, data_count);

    delete[] wavelet_output.first;
    delete[] restored_data;
    delete[] wavelet_output.second;

    // demo for fft
    // data count must be power of 2 for fft!

```

(continues on next page)

(continued from previous page)

```

        std::complex<double> *fft_data = DataFilter::perform_fft (
            data[eeg_channels[i]], data_count, (int)WindowFunctions::NO_WINDOW);
        // len of fft_data array is N / 2 + 1
        std::cout << "FFT coeffs:" << std::endl;
        for (int i = 0; i < data_count / 2 + 1; i++)
        {
            std::cout << fft_data[i] << " ";
        }
        std::cout << std::endl;
        double *restored_from_fft_data = DataFilter::perform_ifft (fft_data, data_
→count);

        std::cout << "Restored after inverse fft transform data:" << std::endl;
        print_one_row (restored_from_fft_data, data_count);

        delete[] fft_data;
        delete[] restored_from_fft_data;
    }
}
catch (const BrainFlowException &err)
{
    BoardShim::log_message ((int)LogLevels::LEVEL_ERROR, err.what ());
    res = err.exit_code;
}

if (data != NULL)
{
    for (int i = 0; i < num_rows; i++)
    {
        delete[] data[i];
    }
}
delete[] data;
delete[] eeg_channels;
delete board;

return res;
}

void print_one_row (double *data, int num_data_points)
{
    for (int i = 0; i < num_data_points; i++)
    {
        std::cout << data[i] << " ";
    }
    std::cout << std::endl;
}

```

5.4.6 C++ Signal Filtering

```

#include <iostream>
#include <stdlib.h>

#ifdef _WIN32
#include <windows.h>
#else
#include <unistd.h>
#endif

#include "board_shim.h"
#include "data_filter.h"

using namespace std;

void print_head (double **data_buf, int num_channels, int num_data_points);

int main (int argc, char *argv[])
{
    struct BrainFlowInputParams params;
    // use synthetic board for demo
    int board_id = (int)BoardIds::SYNTHETIC_BOARD;

    BoardShim::enable_dev_board_logger ();
    DataFilter::enable_dev_data_logger ();

    BoardShim *board = new BoardShim (board_id, params);
    double **data = NULL;
    int *eeg_channels = NULL;
    int num_rows = 0;
    int res = 0;

    try
    {
        board->prepare_session ();
        board->start_stream ();
        BoardShim::log_message ((int)LogLevels::LEVEL_INFO, "Start sleeping in the_
↪main thread");
#ifdef _WIN32
        Sleep (5000);
#else
        sleep (5);
#endif

        board->stop_stream ();
        int data_count = 0;
        data = board->get_board_data (&data_count);
        BoardShim::log_message ((int)LogLevels::LEVEL_INFO, "read %d packages", data_
↪count);
        board->release_session ();
        num_rows = BoardShim::get_num_rows (board_id);
        std::cout << std::endl << "Data from the board" << std::endl << std::endl;
        print_head (data, num_rows, data_count);

        int eeg_num_channels = 0;
        eeg_channels = BoardShim::get_eeg_channels (board_id, &eeg_num_channels);
    }
}

```

(continues on next page)

(continued from previous page)

```

    int filtered_size = 0;
    double *downsampled_data = NULL;
    for (int i = 0; i < eeg_num_channels; i++)
    {
        switch (i)
        {
            // just for test and demo - apply different filters to different eeg_
←channels
            // signal filtering methods work in-place
            case 0:
                DataFilter::perform_lowpass (data[eeg_channels[i]], data_count,
                    BoardShim::get_sampling_rate (board_id), 30.0, 3,
                    (int)FilterTypes::BUTTERWORTH, 0);
                break;
            case 1:
                DataFilter::perform_highpass (data[eeg_channels[i]], data_count,
                    BoardShim::get_sampling_rate (board_id), 5.0, 5,
                    (int)FilterTypes::CHEBYSHEV_TYPE_1, 1);
                break;
            case 2:
                DataFilter::perform_bandpass (data[eeg_channels[i]], data_count,
                    BoardShim::get_sampling_rate (board_id), 15.0, 5.0, 3,
                    (int)FilterTypes::BESSEL, 0);
                break;
            default:
                DataFilter::perform_bandstop (data[eeg_channels[i]], data_count,
                    BoardShim::get_sampling_rate (board_id), 30.0, 1.0, 3,
                    (int)FilterTypes::BUTTERWORTH, 0);
                break;
        }
    }
    std::cout << std::endl << "Data after processing" << std::endl << std::endl;
    print_head (data, num_rows, data_count);
}
catch (const BrainFlowException &err)
{
    BoardShim::log_message ((int)LogLevels::LEVEL_ERROR, err.what ());
    res = err.exit_code;
}

if (data != NULL)
{
    for (int i = 0; i < num_rows; i++)
    {
        delete[] data[i];
    }
}
delete[] data;
delete[] eeg_channels;
delete board;

return res;
}

void print_head (double **data_buf, int num_channels, int num_data_points)
{
    std::cout << "Total Channels for this board: " << num_channels << std::endl;

```

(continues on next page)

(continued from previous page)

```

int num_points = (num_data_points < 5) ? num_data_points : 5;
for (int i = 0; i < num_channels; i++)
{
    std::cout << "Channel " << i << ": ";
    for (int j = 0; j < num_points; j++)
    {
        std::cout << data_buf[i][j] << ",";
    }
    std::cout << std::endl;
}
}

```

5.4.7 C++ Denoising

```

#include <iostream>
#include <stdlib.h>

#ifdef _WIN32
#include <windows.h>
#else
#include <unistd.h>
#endif

#include "board_shim.h"
#include "data_filter.h"

using namespace std;

void print_head (double **data_buf, int num_channels, int num_data_points);

int main (int argc, char *argv[])
{
    struct BrainFlowInputParams params;
    // use synthetic board for demo
    int board_id = (int)BoardIds::SYNTHETIC_BOARD;

    BoardShim::enable_dev_board_logger ();

    BoardShim *board = new BoardShim (board_id, params);
    double **data = NULL;
    int *eeg_channels = NULL;
    int num_rows = 0;
    int res = 0;

    try
    {
        board->prepare_session ();
        board->start_stream ();
        BoardShim::log_message ((int)LogLevels::LEVEL_INFO, "Start sleeping in the_
↪main thread");
#ifdef _WIN32
        Sleep (5000);
#else
        sleep (5);
#endif
    }
}

```

(continues on next page)

(continued from previous page)

```

board->stop_stream ();
int data_count = 0;
data = board->get_board_data (&data_count);
BoardShim::log_message ((int)LogLevels::LEVEL_INFO, "read %d packages", data_
↪count);
board->release_session ();
num_rows = BoardShim::get_num_rows (board_id);
std::cout << std::endl << "Data from the board" << std::endl << std::endl;
print_head (data, num_rows, data_count);

int eeg_num_channels = 0;
eeg_channels = BoardShim::get_eeg_channels (board_id, &eeg_num_channels);
for (int i = 0; i < eeg_num_channels; i++)
{
    switch (i)
    {
        // for demo apply different methods to different channels, in_
↪production you should
        // choose one
        // first of all you can try simple moving average or moving median to_
↪remove noise
        case 0:
            DataFilter::perform_rolling_filter (
                data[eeg_channels[i]], data_count, 3, (int)AggOperations::
↪MEDIAN);
            break;
        case 1:
            DataFilter::perform_rolling_filter (
                data[eeg_channels[i]], data_count, 3, (int)AggOperations::
↪MEAN);
            break;
        case 2:
            DataFilter::perform_rolling_filter (
                data[eeg_channels[i]], data_count, 5, (int)AggOperations::
↪MEDIAN);
            break;
        case 3:
            DataFilter::perform_rolling_filter (
                data[eeg_channels[i]], data_count, 5, (int)AggOperations::
↪MEAN);
            break;
        // if moving average and moving median dont work well for your_
↪signal you can
        // try wavelet based denoising, feel free to try different_
↪wavelet functions and
        // decomposition levels
        case 4:
            DataFilter::perform_wavelet_denoising (
                data[eeg_channels[i]], data_count, "db4", 3);
            break;
        case 5:
            DataFilter::perform_wavelet_denoising (
                data[eeg_channels[i]], data_count, "coif3", 3);
            break;
    }
}

```

(continues on next page)

(continued from previous page)

```

        std::cout << std::endl << "Data after denoising" << std::endl << std::endl;
        print_head (data, num_rows, data_count);
    }
    catch (const BrainFlowException &err)
    {
        BoardShim::log_message ((int)LogLevels::LEVEL_ERROR, err.what ());
        res = err.exit_code;
    }

    if (data != NULL)
    {
        for (int i = 0; i < num_rows; i++)
        {
            delete[] data[i];
        }
    }
    delete[] data;
    delete[] eeg_channels;
    delete board;

    return res;
}

void print_head (double **data_buf, int num_channels, int num_data_points)
{
    std::cout << "Total Channels for this board: " << num_channels << std::endl;
    int num_points = (num_data_points < 5) ? num_data_points : 5;
    for (int i = 0; i < num_channels; i++)
    {
        std::cout << "Channel " << i << ": ";
        for (int j = 0; j < num_points; j++)
        {
            std::cout << data_buf[i][j] << ", ";
        }
        std::cout << std::endl;
    }
}

```

5.4.8 C++ Band Power

```

#include <iostream>
#include <stdlib.h>

#ifdef _WIN32
#include <windows.h>
#else
#include <unistd.h>
#endif

#include "board_shim.h"
#include "data_filter.h"

using namespace std;

int main (int argc, char *argv[])

```

(continues on next page)

(continued from previous page)

```

{
    struct BrainFlowInputParams params;
    // use synthetic board for demo
    int board_id = (int)BoardIds::SYNTHETIC_BOARD;

    BoardShim::enable_dev_board_logger ();

    BoardShim *board = new BoardShim (board_id, params);
    double **data = NULL;
    int *eeg_channels = NULL;
    int num_rows = 0;
    int res = 0;
    int sampling_rate = BoardShim::get_sampling_rate (board_id);

    try
    {
        board->prepare_session ();
        board->start_stream ();
        BoardShim::log_message ((int)LogLevels::LEVEL_INFO, "Start sleeping in the_
↪main thread");
#ifdef _WIN32
        Sleep (10000);
#else
        sleep (10);
#endif

        board->stop_stream ();
        int data_count = 0;
        int fft_len = DataFilter::get_nearest_power_of_two (sampling_rate);
        data = board->get_board_data (&data_count);
        board->release_session ();
        num_rows = BoardShim::get_num_rows (board_id);

        int eeg_num_channels = 0;
        eeg_channels = BoardShim::get_eeg_channels (board_id, &eeg_num_channels);
        // for synthetic board second channel is a sine wave at 10 Hz, should see big_
↪alpha
        int channel = eeg_channels[1];
        // optional - detrend
        DataFilter::detrend (data[channel], data_count, (int)DetrendOperations::
↪LINEAR);
        std::pair<double *, double *> psd = DataFilter::get_psd_welch (data[channel],_
↪data_count,
            fft_len, fft_len / 2, sampling_rate, (int)WindowFunctions::HANNING);
        // calc band power
        double band_power_alpha = DataFilter::get_band_power (psd, fft_len / 2 + 1, 7.
↪0, 13.0);
        double band_power_beta = DataFilter::get_band_power (psd, fft_len / 2 + 1, 14.
↪0, 30.0);
        std::cout << "alpha/beta:" << band_power_alpha / band_power_beta << std::endl;
        // fail test if unexpected ratio
        if (band_power_alpha / band_power_beta < 100)
        {
            res = -1;
        }
        delete[] psd.first;
        delete[] psd.second;
    }
}

```

(continues on next page)

(continued from previous page)

```

    }
    catch (const BrainFlowException &err)
    {
        BoardShim::log_message ((int)LogLevels::LEVEL_ERROR, err.what ());
        res = err.exit_code;
    }

    if (data != NULL)
    {
        for (int i = 0; i < num_rows; i++)
        {
            delete[] data[i];
        }
        delete[] data;
        delete[] eeg_channels;
        delete board;

        return res;
    }
}

```

5.4.9 C++ EEG Metrics

```

#include <chrono>
#include <iostream>
#include <stdlib.h>

#ifdef _WIN32
#include <windows.h>
#else
#include <unistd.h>
#endif

#include "board_shim.h"
#include "data_filter.h"
#include "ml_model.h"

using namespace std;
using namespace std::chrono;

bool parse_args (int argc, char *argv[], struct BrainFlowInputParams *params, int_
↪ *board_id);

int main (int argc, char *argv[])
{
    struct BrainFlowInputParams params;
    int board_id = 0;
    if (!parse_args (argc, argv, &params, &board_id))
    {
        return -1;
    }

    BoardShim::enable_dev_board_logger ();

```

(continues on next page)

(continued from previous page)

```

BoardShim *board = new BoardShim (board_id, params);
double **data = NULL;
int *eeg_channels = NULL;
int num_rows = 0;
int res = 0;
int master_board_id = board->get_board_id ();

int sampling_rate = BoardShim::get_sampling_rate (master_board_id);

try
{
    // Collect data from device
    board->prepare_session ();
    board->start_stream ();
    BoardShim::log_message ((int)LogLevels::LEVEL_INFO, "Start sleeping in the_
↪main thread");
    // recommended window size for eeg metric calculation is at least 4 seconds,
↪bigger is
    // better
#ifdef _WIN32
    Sleep (5000);
#else
    sleep (5);
#endif

    int data_count = 0;
    data = board->get_board_data (&data_count);
    board->stop_stream ();
    std::cout << "Data Count: " << data_count << std::endl;
    board->release_session ();
    num_rows = BoardShim::get_num_rows (master_board_id);

    // Calc bandpowers and build feature vector
    int eeg_num_channels = 0;
    eeg_channels = BoardShim::get_eeg_channels (master_board_id, &eeg_num_
↪channels);
    std::pair<double *, double *> bands = DataFilter::get_avg_band_powers (
        data, data_count, eeg_channels, eeg_num_channels, sampling_rate, true);
    double feature_vector[10];
    for (int i = 0; i < 5; i++)
    {
        feature_vector[i] = bands.first[i];
        feature_vector[i + 5] = bands.second[i];
    }
    for (int i = 0; i < 10; i++)
    {
        std::cout << feature_vector[i] << " ";
    }
    std::cout << std::endl;

    // Testing all classifiers and metric types
    struct BrainFlowModelParams conc_model_params ((int)BrainFlowMetrics::
↪CONCENTRATION, (int)BrainFlowClassifiers::REGRESSION);
    MLModel concentration_model (conc_model_params);
    concentration_model.prepare ();
    std::cout << "Concentration Regression : " << concentration_model.predict_
↪(feature_vector, 10) << std::endl;
    concentration_model.release ();

```

(continues on next page)

(continued from previous page)

```

        struct BrainFlowModelParams relax_model_params ((int)BrainFlowMetrics::
↪RELAXATION, (int)BrainFlowClassifiers::KNN);
        MLModel relaxation_model (relax_model_params);
        relaxation_model.prepare ();
        std::cout << "Relaxation KNN :" << relaxation_model.predict (feature_vector,
↪10) << std::endl;
        relaxation_model.release ();

        delete[] bands.first;
        delete[] bands.second;
    }
    catch (const BrainFlowException &err)
    {
        BoardShim::log_message ((int)LogLevels::LEVEL_ERROR, err.what ());
        res = err.exit_code;
    }

    if (data != NULL)
    {
        for (int i = 0; i < num_rows; i++)
        {
            delete[] data[i];
        }
    }
    delete[] data;
    delete[] eeg_channels;
    delete board;

    return res;
}

bool parse_args (int argc, char *argv[], struct BrainFlowInputParams *params, int_
↪*board_id)
{
    bool board_id_found = false;
    for (int i = 1; i < argc; i++)
    {
        if (std::string (argv[i]) == std::string ("--board-id"))
        {
            if (i + 1 < argc)
            {
                i++;
                board_id_found = true;
                *board_id = std::stoi (std::string (argv[i]));
            }
            else
            {
                std::cerr << "missed argument" << std::endl;
                return false;
            }
        }
        if (std::string (argv[i]) == std::string ("--ip-address"))
        {
            if (i + 1 < argc)
            {

```

(continues on next page)

(continued from previous page)

```

        i++;
        params->ip_address = std::string (argv[i]);
    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
if (std::string (argv[i]) == std::string ("--ip-port"))
{
    if (i + 1 < argc)
    {
        i++;
        params->ip_port = std::stoi (std::string (argv[i]));
    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
if (std::string (argv[i]) == std::string ("--serial-port"))
{
    if (i + 1 < argc)
    {
        i++;
        params->serial_port = std::string (argv[i]);
    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
if (std::string (argv[i]) == std::string ("--ip-protocol"))
{
    if (i + 1 < argc)
    {
        i++;
        params->ip_protocol = std::stoi (std::string (argv[i]));
    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
if (std::string (argv[i]) == std::string ("--timeout"))
{
    if (i + 1 < argc)
    {
        i++;
        params->timeout = std::stoi (std::string (argv[i]));
    }
    else
    {

```

(continues on next page)

(continued from previous page)

```

        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
if (std::string (argv[i]) == std::string ("--other-info"))
{
    if (i + 1 < argc)
    {
        i++;
        params->other_info = std::string (argv[i]);
    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
if (std::string (argv[i]) == std::string ("--mac-address"))
{
    if (i + 1 < argc)
    {
        i++;
        params->mac_address = std::string (argv[i]);
    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
if (std::string (argv[i]) == std::string ("--serial-number"))
{
    if (i + 1 < argc)
    {
        i++;
        params->serial_number = std::string (argv[i]);
    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
if (std::string (argv[i]) == std::string ("--file"))
{
    if (i + 1 < argc)
    {
        i++;
        params->file = std::string (argv[i]);
    }
    else
    {
        std::cerr << "missed argument" << std::endl;
        return false;
    }
}
}

```

(continues on next page)

(continued from previous page)

```
if (!board_id_found)
{
    std::cerr << "board id is not provided" << std::endl;
    return false;
}
return true;
}
```

5.5 R

5.5.1 R Get Data from a Board

```
library (brainflow)

params <- brainflow_python$BrainFlowInputParams ()
board_shim <- brainflow_python$BoardShim (brainflow_python$BoardIds$SYNTHETIC_BOARD
↪$value, params)
board_shim$prepare_session ()
board_shim$start_stream ()
Sys.sleep (time = 5)
board_shim$stop_stream ()
data <- board_shim$get_current_board_data (as.integer (250))
board_shim$release_session ()
```

5.5.2 R Read Write File

```
library (brainflow)

params <- brainflow_python$BrainFlowInputParams ()
board_shim <- brainflow_python$BoardShim (brainflow_python$BoardIds$SYNTHETIC_BOARD
↪$value, params)
board_shim$prepare_session ()
board_shim$start_stream ()
Sys.sleep (time = 5)
board_shim$stop_stream ()
data <- board_shim$get_current_board_data (as.integer (250))
board_shim$release_session ()

brainflow_python$DataFilter$write_file (data, "test.csv", "w")
data_restored <- brainflow_python$DataFilter$read_file ("test.csv")
print (restored_data)
```

5.5.3 R Transforms

```
library (brainflow)

params <- brainflow_python$BrainFlowInputParams ()
board_shim <- brainflow_python$BoardShim (brainflow_python$BoardIds$SYNTHETIC_BOARD
  ↪$value, params)
board_shim$prepare_session ()
board_shim$start_stream ()
Sys.sleep (time = 5)
board_shim$stop_stream ()
data <- board_shim$get_current_board_data (as.integer (250))
board_shim$release_session ()

# need to convert to numpy array manually
numpy_data <- np$array (data[2,])
print (numpy_data)
wavelet_data <- brainflow_python$DataFilter$perform_wavelet_transform (numpy_data,
  ↪"db4", as.integer (3))
restored_data <- brainflow_python$DataFilter$perform_inverse_wavelet_transform_
  ↪(wavelet_data, length (numpy_data), "db4", as.integer (3))
print (restored_data)
```

5.5.4 R Signal Filtering

```
library (brainflow)

params <- brainflow_python$BrainFlowInputParams ()
board_shim <- brainflow_python$BoardShim (brainflow_python$BoardIds$SYNTHETIC_BOARD
  ↪$value, params)
board_shim$prepare_session ()
board_shim$start_stream ()
Sys.sleep (time = 5)
board_shim$stop_stream ()
data <- board_shim$get_current_board_data (as.integer (250))
board_shim$release_session ()

# need to convert to numpy array manually
numpy_data <- np$array (data[2,])
print (numpy_data)
sampling_rate <- board_shim$get_sampling_rate (brainflow_python$BoardIds$SYNTHETIC_
  ↪BOARD$value)
brainflow_python$DataFilter$perform_bandpass (numpy_data, sampling_rate, 10.0, 5.0,
  ↪as.integer (3), brainflow_python$FilterTypes$BESSEL$value, 0)
print (numpy_data)
```

5.5.5 R Denoising

```
library (brainflow)

params <- brainflow_python$BrainFlowInputParams ()
board_shim <- brainflow_python$BoardShim (brainflow_python$BoardIds$SYNTHETIC_BOARD
  ↪$value, params)
board_shim$prepare_session ()
board_shim$start_stream ()
Sys.sleep (time = 5)
board_shim$stop_stream ()
data <- board_shim$get_current_board_data (as.integer (250))
board_shim$release_session ()

# need to convert to numpy array manually
numpy_data <- np$array (data[2,])
print (numpy_data)
brainflow_python$DataFilter$perform_wavelet_denoising (numpy_data, "db4", as.integer_
  ↪(3))
print (numpy_data)
```

5.5.6 R Band Power

```
library (brainflow)

board_id <- brainflow_python$BoardIds$SYNTHETIC_BOARD$value
sampling_rate <- brainflow_python$BoardShim$get_sampling_rate(board_id)
nfft <- brainflow_python$DataFilter$get_nearest_power_of_two(sampling_rate)
params <- brainflow_python$BrainFlowInputParams ()
board_shim <- brainflow_python$BoardShim (board_id, params)
board_shim$prepare_session ()
board_shim$start_stream ()
Sys.sleep (time = 10)
board_shim$stop_stream ()
data <- board_shim$get_board_data ()
board_shim$release_session ()

# need to convert to numpy array manually
numpy_data <- np$array (data[3,])
psd <- brainflow_python$DataFilter$get_psd_welch (numpy_data, as.integer (nfft), as.
  ↪integer (nfft / 2),
  sampling_rate, brainflow_python$WindowFunctions$BLACKMAN_HARRIS$value)
band_power_alpha <- brainflow_python$DataFilter$get_band_power (psd, 7.0, 13.0)
band_power_beta <- brainflow_python$DataFilter$get_band_power (psd, 14.0, 30.0)
ratio <- band_power_alpha / band_power_beta
```

5.5.7 R EEG Metrics

```
library (brainflow)

board_id <- brainflow_python$BoardIds$SYNTHETIC_BOARD$value
sampling_rate <- brainflow_python$BoardShim$get_sampling_rate(board_id)
nfft <- brainflow_python$DataFilter$get_nearest_power_of_two(sampling_rate)
params <- brainflow_python$BrainFlowInputParams ()
board_shim <- brainflow_python$BoardShim (board_id, params)
board_shim$prepare_session ()
board_shim$start_stream ()
Sys.sleep (time = 10)
board_shim$stop_stream ()
data <- board_shim$get_board_data ()
board_shim$release_session ()

eeg_channels <- brainflow_python$BoardShim$get_eeg_channels (board_id)
bands <- brainflow_python$DataFilter$get_avg_band_powers (data, eeg_channels,
↳sampling_rate, TRUE)
feature_vector <- np$array(c(bands[[1]], bands[[2]]))

concentration_params <- brainflow_python$BrainFlowModelParams (brainflow_python
↳$BrainFlowMetrics$CONCENTRATION$value, brainflow_python$BrainFlowClassifiers
↳$REGRESSION$value)
concentration <- brainflow_python$MLModel (concentration_params)
concentration$prepare()
score <- concentration$predict(feature_vector)
concentration$release()
```

5.6 Matlab

5.6.1 Matlab Get Data from a Board

```
BoardShim.set_log_file ('brainflow.log');
BoardShim.enable_dev_board_logger ();

params = BrainFlowInputParams ();
board_shim = BoardShim (int32 (BoardIDs.SYNTHETIC_BOARD), params);
board_shim.prepare_session ();
a = board_shim.config_board ('~6');
board_shim.start_stream (45000, '');
pause (5);
board_shim.stop_stream ();
data = board_shim.get_current_board_data (10);
disp (data);
board_shim.release_session ();
```

5.6.2 Matlab Read Write File

```
BoardShim.set_log_file ('brainflow.log');
BoardShim.enable_dev_board_logger ();

params = BrainFlowInputParams ();
board_shim = BoardShim (int32 (BoardIDs.SYNTHETIC_BOARD), params);
board_shim.prepare_session ();
board_shim.start_stream (45000, '');
pause (2)
board_shim.stop_stream ()
data = board_shim.get_current_board_data (20);
board_shim.release_session ();

DataFilter.write_file (data, 'data.csv', 'w');
restored_data = DataFilter.read_file ('data.csv');
```

5.6.3 Matlab Transforms

```
BoardShim.set_log_file ('brainflow.log');
BoardShim.enable_dev_board_logger ();

params = BrainFlowInputParams ();
board_shim = BoardShim (int32 (BoardIDs.SYNTHETIC_BOARD), params);
sampling_rate = BoardShim.get_sampling_rate (int32 (BoardIDs.SYNTHETIC_BOARD));
board_shim.prepare_session ();
board_shim.start_stream (45000, '');
pause (5);
board_shim.stop_stream ();
data = board_shim.get_current_board_data (DataFilter.get_nearest_power_of_two_
↳(sampling_rate));
board_shim.release_session ();

eeg_channels = BoardShim.get_eeg_channels (int32 (BoardIDs.SYNTHETIC_BOARD));
% wavelet for first eeg channel %
first_eeg_channel = eeg_channels (1);
original_data = data (first_eeg_channel, :);
[wavelet_data, wavelet_lenghts] = DataFilter.perform_wavelet_transform (original_data,
↳ 'db4', 2);
restored_data = DataFilter.perform_inverse_wavelet_transform (wavelet_data, wavelet_
↳lenghts, size (original_data, 2), 'db4', 2);
% fft for first eeg channel %
fft_data = DataFilter.perform_fft (original_data, int32 (WindowFunctions.NO_WINDOW));
restored_fft_data = DataFilter.perform_ifft (fft_data);
```

5.6.4 Matlab Signal Filtering

```
BoardShim.set_log_file ('brainflow.log');
BoardShim.enable_dev_board_logger ();

params = BrainFlowInputParams ();
board_shim = BoardShim (int32 (BoardIDs.SYNTHETIC_BOARD), params);
board_shim.prepare_session ();
board_shim.start_stream (45000, '');
pause (5);
board_shim.stop_stream ();
data = board_shim.get_current_board_data (64);
board_shim.release_session ();

eeg_channels = BoardShim.get_eeg_channels (int32 (BoardIDs.SYNTHETIC_BOARD));
% apply iir filter to the first eeg channel %
first_eeg_channel = eeg_channels (1);
original_data = data (first_eeg_channel, :);
sampling_rate = BoardShim.get_sampling_rate (int32 (BoardIDs.SYNTHETIC_BOARD));
filtered_data = DataFilter.perform_lowpass (original_data, sampling_rate, 10.0, 3,
↳int32 (FilterTypes.BUTTERWORTH), 0.0);
```

5.6.5 Matlab Denoising

```
BoardShim.set_log_file ('brainflow.log');
BoardShim.enable_dev_board_logger ();

params = BrainFlowInputParams ();
board_shim = BoardShim (int32 (BoardIDs.SYNTHETIC_BOARD), params);
board_shim.prepare_session ();
board_shim.start_stream (45000, '');
pause (5);
board_shim.stop_stream ();
data = board_shim.get_current_board_data (64);
board_shim.release_session ();

eeg_channels = BoardShim.get_eeg_channels (int32 (BoardIDs.SYNTHETIC_BOARD));
% apply wavelet denoising to the first eeg channel %
first_eeg_channel = eeg_channels (1);
noisy_data = data (first_eeg_channel, :);
denoised_data = DataFilter.perform_wavelet_denoising (noisy_data, 'db4', 2);
```

5.6.6 Matlab Band Power

```
BoardShim.set_log_file ('brainflow.log');
BoardShim.enable_dev_board_logger ();

params = BrainFlowInputParams ();
board_shim = BoardShim (int32 (BoardIDs.SYNTHETIC_BOARD), params);
sampling_rate = BoardShim.get_sampling_rate (int32 (BoardIDs.SYNTHETIC_BOARD));
board_shim.prepare_session ();
board_shim.start_stream (45000, '');
pause (10);
```

(continues on next page)

(continued from previous page)

```

board_shim.stop_stream ();
nfft = DataFilter.get_nearest_power_of_two (sampling_rate);
data = board_shim.get_board_data ();
board_shim.release_session ();

eeg_channels = BoardShim.get_eeg_channels (int32 (BoardIDs.SYNTHETIC_BOARD));
eeg_channel = eeg_channels (2);
original_data = data (eeg_channel, :);
detrended = DataFilter.detrend (original_data, int32 (DetrendOperations.LINEAR));
[ampls, freqs] = DataFilter.get_psd_welch (detrended, nfft, nfft / 2, sampling_rate,
↳int32 (WindowFunctions.HANNING));
band_power_alpha = DataFilter.get_band_power (ampls, freqs, 7.0, 13.0);
band_power_beta = DataFilter.get_band_power (ampls, freqs, 14.0, 30.0);
ratio = band_power_alpha / band_power_beta;

```

5.6.7 Matlab EEG Metrics

```

BoardShim.set_log_file ('brainflow.log');
BoardShim.enable_dev_board_logger ();

params = BrainFlowInputParams ();
board_shim = BoardShim (int32 (BoardIDs.SYNTHETIC_BOARD), params);
sampling_rate = BoardShim.get_sampling_rate (int32 (BoardIDs.SYNTHETIC_BOARD));
board_shim.prepare_session ();
board_shim.start_stream (45000, '');
pause (5);
board_shim.stop_stream ();
nfft = DataFilter.get_nearest_power_of_two (sampling_rate);
data = board_shim.get_board_data ();
board_shim.release_session ();

eeg_channels = BoardShim.get_eeg_channels (int32 (BoardIDs.SYNTHETIC_BOARD));
[avgs, stddevs] = DataFilter.get_avg_band_powers (data, eeg_channels, sampling_rate,
↳true);
feature_vector = double([avgs, stddevs]);

concentration_params = BrainFlowModelParams (int32 (BrainFlowMetrics.CONCENTRATION),
↳int32 (BrainFlowClassifiers.REGRESSION));
concentration = MLModel (concentration_params);
concentration.prepare ();
score = concentration.predict (feature_vector);
concentration.release ();

```

5.7 Julia

5.7.1 Julia Get Data from a Board

```

import brainflow

# specify logging library to use

```

(continues on next page)

(continued from previous page)

```

brainflow.enable_dev_brainflow_logger(Integer(brainflow.BOARD_CONTROLLER))

params = brainflow.BrainFlowInputParams()
board_shim = brainflow.BoardShim(Integer(brainflow.SYNTHETIC_BOARD), params)

brainflow.prepare_session(board_shim)
brainflow.start_stream(board_shim)
sleep(5)
brainflow.stop_stream(board_shim)
data = brainflow.get_current_board_data(256, board_shim)
brainflow.release_session(board_shim)

```

5.7.2 Julia Read Write File

```

import brainflow

# specify logging library to use
brainflow.enable_dev_brainflow_logger(Integer(brainflow.BOARD_CONTROLLER))

params = brainflow.BrainFlowInputParams()
board_shim = brainflow.BoardShim(Integer(brainflow.SYNTHETIC_BOARD), params)

brainflow.prepare_session(board_shim)
brainflow.start_stream(board_shim)
sleep(5)
brainflow.stop_stream(board_shim)
data = brainflow.get_current_board_data(32, board_shim)
brainflow.release_session(board_shim)

brainflow.write_file(data, "test.csv", "w")
restored_data = brainflow.read_file("test.csv")

println("Original Data")
println(data)
println("Restored Data")
println(restored_data)

```

5.7.3 Julia Transforms

```

import brainflow

# enable logs
brainflow.enable_dev_brainflow_logger(Integer(brainflow.BOARD_CONTROLLER))
brainflow.enable_dev_brainflow_logger(Integer(brainflow.DATA_HANDLER))

params = brainflow.BrainFlowInputParams()
board_shim = brainflow.BoardShim(Integer(brainflow.SYNTHETIC_BOARD), params)
sampling_rate = brainflow.get_sampling_rate(Integer(brainflow.SYNTHETIC_BOARD))

brainflow.prepare_session(board_shim)
brainflow.start_stream(board_shim)

```

(continues on next page)

(continued from previous page)

```

sleep(5)
brainflow.stop_stream(board_shim)
data = brainflow.get_current_board_data(brainflow.get_nearest_power_of_two(sampling_
↪rate), board_shim)
brainflow.release_session(board_shim)

eeg_channels = brainflow.get_eeg_channels(Integer(brainflow.SYNTHETIC_BOARD))
data_first_channel = data[eeg_channels[1], :]

# returns tuple of wavelet coeffs and lengths
wavelet_data = brainflow.perform_wavelet_transform(data_first_channel, "db4", 2)
restored_wavelet_data = brainflow.perform_inverse_wavelet_transform(wavelet_data, ↪
↪length(data_first_channel), "db4", 2)

fft_data = brainflow.perform_fft(data_first_channel, Integer(brainflow.NO_WINDOW))
restored_fft_data = brainflow.perform_ifft(fft_data)

println("Original Data")
println(data_first_channel)
println("Restored from Wavelet Data")
println(restored_wavelet_data)
println("Restored from FFT Data")
println(restored_fft_data)

```

5.7.4 Julia Signal Filtering

```

import brainflow

# specify logging library to use
brainflow.enable_dev_brainflow_logger(Integer(brainflow.BOARD_CONTROLLER))

params = brainflow.BrainFlowInputParams()
board_shim = brainflow.BoardShim(Integer(brainflow.SYNTHETIC_BOARD), params)

brainflow.prepare_session(board_shim)
brainflow.start_stream(board_shim)
sleep(5)
brainflow.stop_stream(board_shim)
data = brainflow.get_current_board_data(32, board_shim)
brainflow.release_session(board_shim)

eeg_channels = brainflow.get_eeg_channels(Integer(brainflow.SYNTHETIC_BOARD))
sampling_rate = brainflow.get_sampling_rate(Integer(brainflow.SYNTHETIC_BOARD))

data_first_channel = data[eeg_channels[1], :]
println("Original Data First Channel")
println(data_first_channel)
brainflow.perform_lowpass(data_first_channel, sampling_rate, 10.0, 3, ↪
↪Integer(brainflow.BUTTERWORTH), 0.0)
println("After LowPass Filter")
println(data_first_channel)

data_second_channel = data[eeg_channels[2], :]
println("Original Data Second Channel")

```

(continues on next page)

(continued from previous page)

```
println(data_second_channel)
brainflow.perform_highpass(data_second_channel, sampling_rate, 5.0, 3,
↳Integer(brainflow.CHEBYSHEV_TYPE_1), 1.0)
println("After HighPass Filter")
println(data_second_channel)

data_third_channel = data[eeg_channels[3], :]
println("Original Data Third Channel")
println(data_third_channel)
brainflow.perform_bandpass(data_third_channel, sampling_rate, 25.0, 20.0, 3,
↳Integer(brainflow.BESSEL), 0.0)
println("After BandPass Filter")
println(data_third_channel)

data_fourth_channel = data[eeg_channels[4], :]
println("Original Data Fourth Channel")
println(data_fourth_channel)
brainflow.perform_bandstop(data_fourth_channel, sampling_rate, 50.0, 2.0, 3,
↳Integer(brainflow.BESSEL), 0.0)
println("After BandStop Filter")
println(data_fourth_channel)
```

5.7.5 Julia Denoising

```
import brainflow

# specify logging library to use
brainflow.enable_dev_brainflow_logger(Integer(brainflow.BOARD_CONTROLLER))

params = brainflow.BrainFlowInputParams()
board_shim = brainflow.BoardShim(Integer(brainflow.SYNTHETIC_BOARD), params)

brainflow.prepare_session(board_shim)
brainflow.start_stream(board_shim)
sleep(5)
brainflow.stop_stream(board_shim)
data = brainflow.get_current_board_data(32, board_shim)
brainflow.release_session(board_shim)

eeg_channels = brainflow.get_eeg_channels(Integer(brainflow.SYNTHETIC_BOARD))
sampling_rate = brainflow.get_sampling_rate(Integer(brainflow.SYNTHETIC_BOARD))

data_first_channel = data[eeg_channels[1], :]
println("Original Data First Channel")
println(data_first_channel)
brainflow.perform_rolling_filter(data_first_channel, 3, Integer(brainflow.MEAN))
println("After Rolling Filter")
println(data_first_channel)

data_second_channel = data[eeg_channels[2], :]
println("Original Data Second Channel")
println(data_second_channel)
brainflow.perform_wavelet_denoising(data_second_channel, "db4", 2)
println("After Wavelet Denoising")
```

(continues on next page)

(continued from previous page)

```
println(data_second_channel)
```

5.7.6 Julia Band Power

```
import brainflow

# specify logging library to use
brainflow.enable_dev_brainflow_logger(Integer(brainflow.BOARD_CONTROLLER))

params = brainflow.BrainFlowInputParams()
board_shim = brainflow.BoardShim(Integer(brainflow.SYNTHETIC_BOARD), params)
sampling_rate = brainflow.get_sampling_rate(Integer(brainflow.SYNTHETIC_BOARD))
nfft = brainflow.get_nearest_power_of_two(sampling_rate)

brainflow.prepare_session(board_shim)
brainflow.start_stream(board_shim)
sleep(5)
brainflow.stop_stream(board_shim)
data = brainflow.get_board_data(board_shim)
brainflow.release_session(board_shim)

eeg_channels = brainflow.get_eeg_channels(Integer(brainflow.SYNTHETIC_BOARD))
# second channel of synthetic board is sine wave at 10 Hz, should see huge 'alpha'
data_second_channel = data[eeg_channels[2], :]

# optional: detrend
brainflow.detrend(data_second_channel, Integer(brainflow.LINEAR))
# psd is a tuple of ampls and freqs
psd = brainflow.get_psd_welch(data_second_channel, nfft, Integer(nfft / 2), sampling_rate, Integer(brainflow.BLACKMAN_HARRIS))
band_power_alpha = brainflow.get_band_power(psd, 7.0, 13.0)
band_power_beta = brainflow.get_band_power(psd, 14.0, 30.0)
println(band_power_alpha / band_power_beta)
```

5.7.7 Julia EEG Metrics

```
import brainflow

# enable all possible logs from all three libs
brainflow.enable_dev_brainflow_logger(Integer(brainflow.BOARD_CONTROLLER))
brainflow.enable_dev_brainflow_logger(Integer(brainflow.DATA_HANDLER))
brainflow.enable_dev_brainflow_logger(Integer(brainflow.ML_MODULE))

params = brainflow.BrainFlowInputParams()
board_shim = brainflow.BoardShim(Integer(brainflow.SYNTHETIC_BOARD), params)
sampling_rate = brainflow.get_sampling_rate(Integer(brainflow.SYNTHETIC_BOARD))
nfft = brainflow.get_nearest_power_of_two(sampling_rate)

brainflow.prepare_session(board_shim)
brainflow.start_stream(board_shim)
```

(continues on next page)

(continued from previous page)

```

sleep(5)
brainflow.stop_stream(board_shim)
data = brainflow.get_board_data(board_shim)
brainflow.release_session(board_shim)

eeg_channels = brainflow.get_eeg_channels(Integer(brainflow.SYNTHETIC_BOARD))

bands = brainflow.get_avg_band_powers(data, eeg_channels, sampling_rate, true)
feature_vector = vcat(bands[1], bands[2])

# calc concentration
model_params = brainflow.BrainFlowModelParams(Integer(brainflow.CONCENTRATION),
↳Integer(brainflow.KNN))
concentration = brainflow.MLModel(model_params)
brainflow.prepare(concentration)
print(brainflow.predict(feature_vector, concentration))
brainflow.release(concentration)

# calc relaxation
model_params = brainflow.BrainFlowModelParams(Integer(brainflow.RELAXATION),
↳Integer(brainflow.REGRESSION))
relaxation = brainflow.MLModel(model_params)
brainflow.prepare(relaxation)
print(brainflow.predict(feature_vector, relaxation))
brainflow.release(relaxation)

```

5.8 Notebooks

5.8.1 BrainFlow to MNE Python Notebook

```

In [1]: import time
import numpy as np
import pandas as pd
import matplotlib.pyplot as plt

import brainflow
from brainflow.board_shim import BoardShim, BrainFlowInputParams, BoardIds

import mne
from mne.channels import read_layout

/home/docs/checkouts/readthedocs.org/user_builds/brainflow-openbci/envs/stable/lib/
↳python3.7/site-packages/traitlets/traitlets.py:3036: FutureWarning: --rc={'figure.
↳dpi': 96} for dict-traits is deprecated in traitlets 5.0. You can pass --rc
↳<key=value> ... multiple times to add items to a dict.
FutureWarning,

```

```

In [2]: # use synthetic board for demo
params = BrainFlowInputParams ()
board = BoardShim (BoardIds.SYNTHETIC_BOARD.value, params)
board.prepare_session ()
board.start_stream ()
time.sleep (10)

```

(continues on next page)

(continued from previous page)

```
data = board.get_board_data ()
board.stop_stream ()
board.release_session ()
```

```
In [3]: eeg_channels = BoardShim.get_eeg_channels (BoardIds.SYNTHETIC_BOARD.value)
eeg_data = data[eeg_channels, :]
eeg_data = eeg_data / 1000000 # BrainFlow returns uV, convert to V for MNE
```

```
In [4]: # Creating MNE objects from brainflow data arrays
ch_types = ['eeg'] * len (eeg_channels)
ch_names = BoardShim.get_eeg_names (BoardIds.SYNTHETIC_BOARD.value)
sfreq = BoardShim.get_sampling_rate (BoardIds.SYNTHETIC_BOARD.value)
info = mne.create_info (ch_names = ch_names, sfreq = sfreq, ch_types = ch_types)
raw = mne.io.RawArray (eeg_data, info)
# its time to plot something!
raw.plot_psd (average = False)
```


Creating RawArray with float64 data, n_channels=16, n_times=3212

Range : 0 ... 3211 = 0.000 ... 12.844 secs


Ready.

Effective window size : 8.192 (s)

<ipython-input-1-9872dda7122d>:8: RuntimeWarning: Channel locations not available.
↳ Disabling spatial colors.
raw.plot_psd (average = False)



Out [4]:



5.8.2 Denoising Notebook

```
In [1]: import argparse
import time
import brainflow
import numpy as np

import pandas as pd
import matplotlib
import matplotlib.pyplot as plt

from brainflow.board_shim import BoardShim, BrainFlowInputParams, LogLevels, BoardIds
from brainflow.data_filter import DataFilter, FilterTypes, AggOperations

/home/docs/checkouts/readthedocs.org/user_builds/brainflow-openbci/envs/stable/lib/
↳ python3.7/site-packages/traitlets/traitlets.py:3036: FutureWarning: --rc={'figure.
↳ dpi': 96} for dict-traits is deprecated in traitlets 5.0. You can pass --rc
↳ <key=value> ... multiple times to add items to a dict.
FutureWarning,
```


```
In [2]: # use synthetic board for demo
params = BrainFlowInputParams ()
board_id = BoardIds.SYNTHETIC_BOARD.value
board = BoardShim (board_id, params)
```

(continues on next page)

(continued from previous page)

```
board.prepare_session ()
board.start_stream ()
time.sleep (20)
data = board.get_board_data ()
board.stop_stream ()
board.release_session ()
```

```
In [3]: # plot original data
eeg_channels = BoardShim.get_eeg_channels (board_id)
df = pd.DataFrame (np.transpose (data))
df[eeg_channels].plot (subplots = True)
plt.show()
```



```
In [4]: # demo for different denoising methods,
# apply different methods to different channels to determine the best one
for count, channel in enumerate (eeg_channels):
    # first of all you can try simple moving median or moving average with different
    ↪ window size
    if count == 0:
        DataFilter.perform_rolling_filter (data[channel], 3, AggOperations.MEAN.value)
    elif count == 1:
        DataFilter.perform_rolling_filter (data[channel], 3, AggOperations.MEDIAN.
    ↪ value)
    # methods above should increase signal to noise ratio but we can do even better
    # using wavelet based denoising, feel free to try different wavelet functions and
    ↪ decomposition levels
    elif count == 2:
        DataFilter.perform_wavelet_denoising (data[channel], 'db6', 5)
    elif count == 3:
        DataFilter.perform_wavelet_denoising (data[channel], 'bior3.9', 5)
```

(continues on next page)

(continued from previous page)

```


elif count == 4:
    DataFilter.perform_wavelet_denoising (data[channel], 'sym7', 5)
elif count == 5:
    DataFilter.perform_wavelet_denoising (data[channel], 'coif3', 5)
elif count == 6:
    DataFilter.perform_wavelet_denoising (data[channel], 'bior6.8', 5)
elif count == 7:
    DataFilter.perform_wavelet_denoising (data[channel], 'db4', 5)

```

```

In [5]: # plot denoised data
df = pd.DataFrame (np.transpose (data))
df[eeg_channels].plot (subplots = True)
plt.show()
# as you can see wavelet based denoising works much better and increases signal to_
↪noise ratio significantly!

```



5.8.3 BrainFlow Band Power Notebook

```

In [1]: import argparse
import time
import brainflow
import numpy as np

import pandas as pd
import matplotlib
import matplotlib.pyplot as plt

from brainflow.board_shim import BoardShim, BrainFlowInputParams, LogLevels, BoardIds

```

(continues on next page)


(continued from previous page)

```
from brainflow.data_filter import DataFilter, FilterTypes, AggOperations,   
↳ WindowFunctions, DetrendOperations  
  
/home/docs/checkouts/readthedocs.org/user_builds/brainflow-openbci/envs/stable/lib/  
↳ python3.7/site-packages/traitlets/traitlets.py:3036: FutureWarning: --rc={'figure.  
↳ dpi': 96} for dict-traits is deprecated in traitlets 5.0. You can pass --rc  
↳ <key=value> ... multiple times to add items to a dict.  
FutureWarning,
```

```
In [2]: # use synthetic board for demo  
params = BrainFlowInputParams ()  
board_id = BoardIds.SYNTHETIC_BOARD.value  
sampling_rate = BoardShim.get_sampling_rate (board_id)  
nfft = DataFilter.get_nearest_power_of_two (sampling_rate)  
board = BoardShim (board_id, params)  
board.prepare_session ()  
board.start_stream ()  
time.sleep (10)  
data = board.get_board_data ()  
board.stop_stream ()  
board.release_session ()  
eeg_channels = BoardShim.get_eeg_channels (board_id)  
# use first eeg channel for demo  
# second channel of synthetic board is a sine wave at 10 Hz, should see big 'alpha'  
eeg_channel = eeg_channels[1]
```

```
In [3]: # optional: detrend  
DataFilter.detrend (data[eeg_channel], DetrendOperations.LINEAR.value)
```

```
In [4]: psd = DataFilter.get_psd_welch (data[eeg_channel], nfft, nfft // 2, sampling_rate,   
↳ WindowFunctions.HANNING.value)  
plt.plot (psd[1][:60], psd[0][:60])  
plt.show ()
```



```
In [5]: # calc band power
alpha = DataFilter.get_band_power (psd, 7.0, 13.0)
beta = DataFilter.get_band_power (psd, 14.0, 30.0)
print ("Alpha/Beta Ratio is: %f" % (alpha / beta))
```

Alpha/Beta Ratio is: 2316.283381

INTEGRATION WITH GAME ENGINES

6.1 Unity


Integration with Unity can be done only using C# binding. And currently it works only for Windows.

You can build C# binding from source or download compiled package directly from [Nuget](#).

Here we will use Nuget to download and install BrainFlow with dependencies.

Download [nuget.exe](#) and run


```
nuget.exe install brainflow -OutputDirectory <OUTPUTDIR>
```




```
Administrator: Command Prompt
D:\Downloads>nuget.exe install brainflow -OutputDirectory D:\BrainFlowInstalled
Feeds used:
  https://api.nuget.org/v3/index.json
  E:\folder\nuget_packages
  C:\Program Files (x86)\Microsoft SDKs\NuGetPackages\
Installing package 'brainflow' to 'D:\BrainFlowInstalled'.
  GET https://api.nuget.org/v3/registration5-gz-semver2/brainflow/index.json
  OK https://api.nuget.org/v3/registration5-gz-semver2/brainflow/index.json 527ms

Attempting to gather dependency information for package 'brainflow.3.1.0' with respect to project 'D:\BrainFlowInstalled', targeting 'Any,Version=v0.0'
Gathering dependency information took 1.22 sec
Attempting to resolve dependencies for package 'brainflow.3.1.0' with DependencyBehavior 'Lowest'
Resolving dependency information took 0 ms
Resolving actions to install package 'brainflow.3.1.0'
Resolved actions to install package 'brainflow.3.1.0'
Retrieving package 'Accord 3.8.0' from 'nuget.org'.
Retrieving package 'Accord.Math 3.8.0' from 'nuget.org'.
```

Open OUTPUTDIR, you will see BrainFlow library and dependencies. For BrainFlow there are *Managed(C#)* and *Unmanaged(C++)* files. For dependencies like Accord there are only *Managed(C#)* libraries.



Open your Unity project and copy **Managed(C#)** libraries from BrainFlow and **all dependencies** to the Assets folder, after that copy **Unmanaged(C++)** libraries from BrainFlow to the root folder of your project.



Now, you are able to use BrainFlow API in your C# scripts!

For demo we will create a simple script to read data and calculate concentration level.

Add a Sphere object to the Scene and attach script below.

```
using System;
using System.Collections;
using System.Collections.Generic;
using UnityEngine;

using Accord;
using Accord.Math;
using brainflow;

public class Sphere : MonoBehaviour
{
    private BoardShim board_shim = null;
    private MLModel concentration = null;
    private int sampling_rate = 0;
    private int[] eeg_channels = null;

    // Start is called before the first frame update
    void Start()
    {
        try
        {
            BoardShim.set_log_file("brainflow_log.txt");
            BoardShim.enable_dev_board_logger();

            BrainFlowInputParams input_params = new BrainFlowInputParams();
            int board_id = (int)BoardIds.SYNTHETIC_BOARD;
            board_shim = new BoardShim(board_id, input_params);
            board_shim.prepare_session();
            board_shim.start_stream(450000, "file://brainflow_data.csv:w");
            BrainFlowModelParams concentration_params = new
            BrainFlowModelParams((int)BrainFlowMetrics.CONCENTRATION, (int)BrainFlowClassifiers.
            REGRESSION);
```

(continues on next page)

```

        concentration = new MLModel(concentration_params);
        concentration.prepare();


        sampling_rate = BoardShim.get_sampling_rate(board_id);
        eeg_channels = BoardShim.get_eeg_channels(board_id);
        Debug.Log("Brainflow streaming was started");
    }
    catch (BrainFlowException e)
    {
        Debug.Log(e);
    }
}

// Update is called once per frame
void Update()
{
    if ((board_shim == null) || (concentration == null))
    {
        return;
    }
    int number_of_data_points = sampling_rate * 4; // 4 second window is
    ↪recommended for concentration and relaxation calculations
    double[,] data = board_shim.get_current_board_data(number_of_data_points);
    if (data.GetRow(0).Length < number_of_data_points)
    {
        // wait for more data
        return;
    }
    // prepare feature vector
    Tuple<double[], double[]> bands = DataFilter.get_avg_band_powers (data, eeg_
    ↪channels, sampling_rate, true);
    double[] feature_vector = bands.Item1.Concatenate (bands.Item2);
    // calc and print concetration level
    // for synthetic board this value should be close to 1, because of sin waves
    ↪ampls and freqs
    Debug.Log("Concentration: " + concentration.predict (feature_vector));
}

private void OnDestroy()
{
    if (board_shim != null)
    {
        try
        {
            board_shim.release_session();
            concentration.release();
        }
        catch (BrainFlowException e)
        {
            Debug.Log(e);
        }
        Debug.Log("Brainflow streaming was stopped");
    }
}
}

```

If everything is fine, you will see Concentration Score in Console.



After building your game you need to copy *Unmanaged(C++)* libraries to a folder where executable is located.

6.2 Unreal Engine

First of all you need to compile BrainFlow from source. For Windows you need to specify an option to link MSVC Runtime *dynamically*. And you need to use the same version of Visual Studio as in your Unreal Project.

Command line example for Windows and MSVC 2017:

```
# install cmake, clone repo and run commands below
cd brainflow
mkdir build_dyn
cd build_dyn
cmake -G "Visual Studio 15 2017 Win64" -DMSVC_RUNTIME=dynamic -DCMAKE_SYSTEM_
↪VERSION=10.0 -DCMAKE_INSTALL_PREFIX=FULL_PATH_TO_FOLDER_FOR_INSTALLATION ..
# e.g. cmake -G "Visual Studio 15 2017 Win64" -DMSVC_RUNTIME=dynamic -DCMAKE_SYSTEM_
↪VERSION=10.0 -DCMAKE_INSTALL_PREFIX=E:\folder\brainflow\installed_temp ..
cmake --build . --target install --config Release -j 2 --parallel 2
```

Add new entry to your *PATH* environment variable to point to a folder *FULL_PATH_TO_FOLDER_FOR_INSTALLATION* in example above it's *E:\folder\brainflow\installed_temp\lib*. If you have Unreal Engine Editor or Visual Studio running at this point you need to restart these processes.

Open your Visual Studio Solution for your Unreal Engine project, here we created a project called *BrainFlowUnreal*.

Edit file named *ProjectName.Build.cs*, in our example this file is called *BrainFlowUnreal.Build.cs*

```
using UnrealBuildTool;
using System.IO;

public class BrainFlowUnreal : ModuleRules
```

(continues on next page)

(continued from previous page)

```

{
    public BrainFlowUnreal(ReadOnlyTargetRules Target) : base(Target)
    {
        PCHUsage = PCHUsageMode.UseExplicitOrSharedPCHs;


        PublicDependencyModuleNames.AddRange(new string[] { "Core", "CoreUObject",
        ↪ "Engine", "InputCore", "HeadMountedDisplay" });

        // here I show it only for Windows, if you wanna target more OSes add similar
        ↪ code for them
        if (Target.Platform == UnrealTargetPlatform.Win64)
        {
            // Add the import library
            PublicLibraryPaths.Add(Path.Combine(ModuleDirectory, "x64"));
            PublicAdditionalLibraries.Add("BrainFlow.lib");
            PublicAdditionalLibraries.Add("DataHandler.lib");
            PublicAdditionalLibraries.Add("BoardController.lib");

            // add headers for static library
            PublicIncludePaths.Add(Path.Combine(ModuleDirectory, "include"));
        }
    }
}

```

After that you need to copy headers and libraries from BrainFlow installation folder to your Unreal Engine project. Here we copied a content of `E:\folder\brainflow\installed_temp\inc` to a folder `E:\gamedev\BrainFlowUnreal\Source\BrainFlowUnreal\include`. Also, you need to copy compiled libraries from `E:\folder\brainflow\installed_temp\lib` to `E:\gamedev\BrainFlowUnreal\Source\BrainFlowUnreal\x64`



Note: in this example we didn't create a new plugin as described [here](#). Also we linked only static libraries and didn't link or load dynamic libraries manually. And we don't recommend to configure it as a plugin.

Finally, you are able to use BrainFlow in your Unreal Engine project.

When you will build a project for production put C++ libraries for BrainFlow in the folder with executable.

BRAINFLOW DEV

7.1 Code style

We use clang-format tool to keep the same code style for all cpp files. You can download clang-format binary from [LLVM Download Page](#) We recommend installing a plugin for your text editor or IDE which will apply clang-format tool during saving. You will need to set code style option to “FILE”

Plugins for text editors and IDEs:

- [Sublime](#)
- [VSCode](#)
- clang-format tool is preinstalled for Visual Studio

Unfortunately clang-format cannot handle naming, so some additional rules are: - methods and variables should be in lower case with underscore - class names should be in camel case - use brackets even for single line if and for statements

We try to keep the same code style for all bindings as well, even if it doesn't match PEP or other standards. For example we add spaces before and after assignment operator to specify default value for method's params and add spaces before brackets.

7.2 CI and tests

If you want to commit to core module of brainflow project please check that all tests are passed, you can enable [Travis CI](#) and [AppVeyour](#) for your fork of BrainFlow to run tests automatically, or check CI status directly in your PR.

Also you can run integration tests manually for any board even if you dont have real hardware using BrainFlow Emulator.

7.3 Pull Requests

Just try to briefly explain a goal of this PR.

7.4 Instructions to add new boards to BrainFlow

- add new board Id to `BoardIds` enum in C code and to the same enum in all bindings
- add new object creation to board controller C interface
- inherit your board from `Board` class and implement all pure virtual methods, store data in `DataBuffer` object, use `synthetic board` as a reference, try to reuse code from `utils` folder
- add information about your board to `brainflow_boards.h`
- add new files to `BOARD_CONTROLLER_SRC` variable in `CmakeLists.txt`, you may also need to add new directory to `target_include_directories` for `BOARD_CONTROLLER_NAME` variable

You've just written Python, Java, C#, R, C++ ... SDKs for your board! Also now you can use your new board with applications and frameworks which use BrainFlow API.

Optional: We use CI to run tests automatically, to add your board to CI pipelines you can develop a simple emulator for your device. Use [emulators for existing boards](#) as a reference and add tests for your device to *Travis* and *Appveyor*.

7.5 Instructions to build docs locally

Don't push changes to Docs without local verification.

- install `pandoc`
- optional: install Doxygen, skip it if you don't understand what it is or don't need to publish your local build

Install requirements:

```
cd docs
python -m pip install -r requirements.txt
```

Build docs:

```
make html
```

7.6 Debug BrainFlow's errors

Since bindings just call methods from dynamic libraries, more likely errors occur in C++ code, it means that you need to use C++ debugger like `gdb`. If there is an error in binding, it should be simple to figure out and resolve the issue using language specific tools.

Steps to get more information about errors in C++ code:

- build BrainFlow's core module and C++ binding in debug mode. In files like `tools/build_linux.sh` default config is `Release`, so you need to change it to `Debug`
- reproduce your issue using C++ binding
- run it with debugger and memory checker

Example for Linux(for MacOS it's the same):

```

vim tools/build_linux.sh
# Change build type to Debug
bash tools/build_linux.sh
# Create a test to reproduce your issue in C++, here we will use get_data_demo
cd tests/cpp/get_data_demo
mkdir build
cd build
cmake -DCMAKE_PREFIX_PATH=TYPE_FULL_PATH_TO_BRAINFLOW_INSTALLED_FOLDER -DCMAKE_BUILD_
↳TYPE=Debug ..
# e.g. cmake -DCMAKE_PREFIX_PATH=/home/andrey/brainflow/installed_linux -DCMAKE_BUILD_
↳TYPE=Debug ..
make
# Run Valgrind to check memory errors
# Here we use command line for Ganglion
sudo valgrind --error-exitcode=1 --leak-check=full ./brainflow_get_data --board-id 1 -
↳-serial-port /dev/ttyACM0 --mac-address e6:73:73:18:09:b1
# Valgrind will print Error Summary and exact line numbers
# Run gdb and get backtrace
sudo gdb --args ./brainflow_get_data --board-id 1 --serial-port /dev/ttyACM0 --mac-
↳address e6:73:73:18:09:b1
# In gdb terminal type 'r' to run the program and as soon as error occurs, type 'bt'
↳to see backtrace with exact lines of code and call stack

```

7.7 BrainFlow Emulator

BrainFlow Emulator allows you to run all integration tests for all supported boards without real hardware. Our CI uses it for test automation. Also, you can run it on your own PC!

7.7.1 Streaming Board

Streaming Board emulator works using Python binding for BrainFlow, so **you need to install Python binding first**.

Install emulator package:

```

cd emulator
python -m pip install -U .

```

Run tests

```

python emulator\brainflow_emulator\streaming_board_emulator.py python tests\python\
↳brainflow_get_data.py --log --board-id -2 --ip-address 225.1.1.1 --ip-port 6677 --
↳other-info -l

```

This emulator uses synthetic board as a master board and the IP address and port are hardcoded.

7.7.2 OpenBCI Cyton

Cyton emulator simulate COM port using:

- `com0com` for Windows
- `pty` for Linux and MacOS

You should pass test command line directly to `cyton_linux.py` or to `cyton_windows.py`. The script will add the port automatically to provided command line and will start an application.

Install emulator package:

```
cd emulator
python -m pip install -U .
```

Run tests for LinuxMacOS and Windows (port argument will be added by Emulator!)

```
python brainflow_emulator/cyton_linux.py python ../tests/python/brainflow_get_data.py -
↳--log --board-id 0 --serial-port
python brainflow_emulator/cyton_windows.py python ../tests/python/brainflow_get_data.
↳py --log --board-id 0 --serial-port
```

7.7.3 OpenBCI NovaXR

NovaXR emulator creates socket server and streams data to BrainFlow like it's a real board, but with much lower sampling rate.

Install emulator package:

```
cd emulator
python -m pip install -U .
```

Run tests:

```
python brainflow_emulator/novaxr_udp.py python ../tests/python/brainflow_get_data.py -
↳-log --ip-address 127.0.0.1 --board-id 3
```

7.7.4 OpenBCI Wifi Shield based boards

Wifi shield emulator starts http server to read commands and creates client socket to stream data.

Install emulator package:

```
cd emulator
python -m pip install -U .
```

Run tests for Ganglion, Cyton and Daisy with Wifi Shield:

```
python brainflow_emulator/wifi_shield_emulator.py python ../tests/python/brainflow_
↳get_data.py --log --ip-address 127.0.0.1 --board-id 4 --ip-protocol 2 --ip-port_
↳17982
python brainflow_emulator/wifi_shield_emulator.py python ../tests/python/brainflow_
↳get_data.py --log --ip-address 127.0.0.1 --board-id 5 --ip-protocol 2 --ip-port_
↳17982
python brainflow_emulator/wifi_shield_emulator.py python ../tests/python/brainflow_
↳get_data.py --log --ip-address 127.0.0.1 --board-id 6 --ip-protocol 2 --ip-port_
↳17982
```

(continues on next page)

(continued from previous page)

--

8.1 Contact Info, Feature Request, Report an Issue

- Join our [slack workspace](#) using [self-invite page](#)
- To report bugs or request features create an issue in our [GitHub Page](#)
- For hardware related questions and issues contact board manufacturer

8.2 Issue format

First of all you need to run your code with:

```
enable_dev_board_logger ()
```

After that, make sure:

- you've specified BrainFlow version and OS version
- you've attached all logs to your issue description
- you've provided steps or a simple example to reproduce your issue

8.3 Contributors

- [Andrey1994](#) - 693 contributions
- [daniellasy](#) - 26 contributions
- [shirleyzhang867](#) - 5 contributions
- [mesca](#) - 5 contributions
- [Fan1117](#) - 5 contributions
- [John42506176Linux](#) - 4 contributions
- [retiutut](#) - 3 contributions
- [alexcastillo](#) - 1 contribution
- [stellarpower](#) - 1 contribution

PARTNERS AND SPONSORS

9.1 OpenBCI

OpenBCI specializes in creating low-cost, high-quality biosensing hardware for brain computer interfacing. Their arduino compatible biosensing boards provide high resolution imaging and recording of EMG, ECG, and EEG signals. Their devices have been used by researchers, makers, and hobbyists in over 60+ countries as brain computer interfaces to power machines and map brain activity. OpenBCI headsets, boards, sensors and electrodes allow anyone interested in biosensing and neurofeedback to purchase high quality equipment at affordable prices.



MIT LICENSE

Copyright (c) 2018 Andrey Parfenov

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the “Software”), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED “AS IS”, WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

PARTNERS AND SPONSORS



SEARCH

- search

B

BoardShim (C++ class), 39
 BoardShim::~BoardShim (C++ function), 39
 BoardShim::board_id (C++ member), 40
 BoardShim::BoardShim (C++ function), 39
 BoardShim::config_board (C++ function), 39
 BoardShim::disable_board_logger (C++ function), 40
 BoardShim::enable_board_logger (C++ function), 40
 BoardShim::enable_dev_board_logger (C++ function), 40
 BoardShim::get_accel_channels (C++ function), 42
 BoardShim::get_analog_channels (C++ function), 42
 BoardShim::get_battery_channel (C++ function), 40
 BoardShim::get_board_data (C++ function), 39
 BoardShim::get_board_data_count (C++ function), 39
 BoardShim::get_board_id (C++ function), 39
 BoardShim::get_current_board_data (C++ function), 39
 BoardShim::get_ecg_channels (C++ function), 41
 BoardShim::get_eda_channels (C++ function), 42
 BoardShim::get_eeg_channels (C++ function), 41
 BoardShim::get_eeg_names (C++ function), 41
 BoardShim::get_emg_channels (C++ function), 41
 BoardShim::get_eog_channels (C++ function), 42
 BoardShim::get_exg_channels (C++ function), 42
 BoardShim::get_gyro_channels (C++ function), 43
 BoardShim::get_num_rows (C++ function), 41
 BoardShim::get_other_channels (C++ function), 43
 BoardShim::get_package_num_channel (C++ function), 40
 BoardShim::get_ppg_channels (C++ function), 42
 BoardShim::get_resistance_channels (C++ function), 43
 BoardShim::get_sampling_rate (C++ function), 40
 BoardShim::get_temperature_channels (C++ function), 43
 BoardShim::get_timestamp_channel (C++ function), 40
 BoardShim::is_prepared (C++ function), 39
 BoardShim::log_message (C++ function), 40
 BoardShim::prepare_session (C++ function), 39
 BoardShim::release_session (C++ function), 39
 BoardShim::set_log_file (C++ function), 40
 BoardShim::set_log_level (C++ function), 40
 BoardShim::start_stream (C++ function), 39
 BoardShim::stop_stream (C++ function), 39
 brainflow::AggOperations (C++ enum), 50, 67
 brainflow::AggOperations::EACH (C++ enumerator), 67
 brainflow::AggOperations::MEAN (C++ enumerator), 67
 brainflow::AggOperations::MEDIAN (C++ enumerator), 67
 brainflow::BoardIds (C++ enum), 50, 66
 brainflow::BoardIds::BRAINBIT_BOARD (C++ enumerator), 66
 brainflow::BoardIds::CALLIBRI_ECG_BOARD (C++ enumerator), 67
 brainflow::BoardIds::CALLIBRI_EEG_BOARD (C++ enumerator), 67
 brainflow::BoardIds::CALLIBRI_EMG_BOARD (C++ enumerator), 67
 brainflow::BoardIds::CYTON_BOARD (C++ enumerator), 66
 brainflow::BoardIds::CYTON_DAISY_BOARD (C++ enumerator), 66

brainflow::BoardIds::CYTON_DAISY_WIFI_BOARD (C++ enumerator), 66
 brainflow::BoardIds::CYTON_WIFI_BOARD (C++ enumerator), 66
 brainflow::BoardIds::FASCIA_BOARD (C++ enumerator), 67
 brainflow::BoardIds::FREEEEG32_BOARD (C++ enumerator), 67
 brainflow::BoardIds::GALEA_BOARD (C++ enumerator), 66
 brainflow::BoardIds::GANGLION_BOARD (C++ enumerator), 66
 brainflow::BoardIds::GANGLION_WIFI_BOARD (C++ enumerator), 66
 brainflow::BoardIds::IRONBCI_BOARD (C++ enumerator), 67
 brainflow::BoardIds::NOTION_1_BOARD (C++ enumerator), 67
 brainflow::BoardIds::NOTION_2_BOARD (C++ enumerator), 67
 brainflow::BoardIds::PLAYBACK_FILE_BOARD (C++ enumerator), 66
 brainflow::BoardIds::STREAMING_BOARD (C++ enumerator), 66
 brainflow::BoardIds::SYNTHETIC_BOARD (C++ enumerator), 66
 brainflow::BoardIds::UNICORN_BOARD (C++ enumerator), 67
 brainflow::brainflow::BoardShim (C++ class), 52, 68
 brainflow::brainflow::BoardShim::board_id (C++ member), 53, 69
 brainflow::brainflow::BoardShim::BoardShim (C++ function), 52, 68
 brainflow::brainflow::BoardShim::config_board (C++ function), 52, 68
 brainflow::brainflow::BoardShim::disable_board_logging (C++ function), 53, 73
 brainflow::brainflow::BoardShim::enable_board_logging (C++ function), 53, 73
 brainflow::brainflow::BoardShim::enable_device_board_logging (C++ function), 53, 73
 brainflow::brainflow::BoardShim::get_battery_level (C++ function), 53, 70
 brainflow::brainflow::BoardShim::get_board_data_rate (C++ function), 52, 69
 brainflow::brainflow::BoardShim::get_board_id (C++ function), 52, 69
 brainflow::brainflow::BoardShim::get_num_channels (C++ function), 53, 70
 brainflow::brainflow::BoardShim::get_package_name (C++ function), 53, 69
 brainflow::brainflow::BoardShim::get_sampling_rate (C++ function), 53, 69
 brainflow::brainflow::BoardShim::get_timestamp_char (C++ function), 53, 69
 brainflow::brainflow::BoardShim::is_prepared (C++ function), 52, 68
 brainflow::brainflow::BoardShim::log_message (C++ function), 53, 74
 brainflow::brainflow::BoardShim::prepare_session (C++ function), 52, 68
 brainflow::brainflow::BoardShim::release_session (C++ function), 52, 68
 brainflow::brainflow::BoardShim::set_log_file (C++ function), 53, 73
 brainflow::brainflow::BoardShim::set_log_level (C++ function), 53, 73
 brainflow::brainflow::BoardShim::start_stream (C++ function), 52, 68
 brainflow::brainflow::BoardShim::stop_stream (C++ function), 52, 68
 brainflow::brainflow::BrainFlowError (C++ class), 55
 brainflow::brainflow::BrainFlowError::BrainFlowError (C++ function), 55
 brainflow::brainflow::BrainFlowError::exit_code (C++ member), 55
 brainflow::brainflow::BrainFlowError::msg (C++ member), 55
 brainflow::brainflow::BrainFlowException (C++ class), 74
 brainflow::brainflow::BrainFlowException::BrainFlowException (C++ function), 74
 brainflow::brainflow::BrainFlowException::exit_code (C++ member), 74
 brainflow::brainflow::BrainFlowInputParams (C++ class), 55, 74
 brainflow::brainflow::BrainFlowInputParams::BrainFlowInputParams (C++ function), 55, 74
 brainflow::brainflow::BrainFlowInputParams::file (C++ member), 56, 75
 brainflow::brainflow::BrainFlowInputParams::get_file_name (C++ function), 56
 brainflow::brainflow::BrainFlowInputParams::get_ip_address (C++ function), 55
 brainflow::brainflow::BrainFlowInputParams::get_ip_port (C++ function), 55
 brainflow::brainflow::BrainFlowInputParams::get_ip_port (C++ function), 55
 brainflow::brainflow::BrainFlowInputParams::get_ip_port (C++ function), 55
 brainflow::brainflow::BrainFlowInputParams::get_mac_address (C++ function), 55
 brainflow::brainflow::BrainFlowInputParams::get_other_params (C++ function), 55
 brainflow::brainflow::BrainFlowInputParams::get_serial_number (C++ function), 55
 brainflow::brainflow::BrainFlowInputParams::get_serial_number (C++ function), 55

brainflow::BrainFlowMetrics (C++ enum), (C++ enumerator), 66
 56, 67
 brainflow::BrainFlowMetrics::CONCENTRATION (C++ enumerator), 68
 brainflow::BrainFlowMetrics::RELAXATION (C++ enumerator), 67
 brainflow::CustomExitCodes (C++ enum), 65
 brainflow::CustomExitCodes::ANOTHER_BOARD_NOT_CREATED_ERROR (C++ enumerator), 66
 brainflow::CustomExitCodes::ANOTHER_CLASSIFIER_IS_NOT_PREPARED_ERROR (C++ enumerator), 66
 brainflow::CustomExitCodes::BOARD_NOT_CREATED_ERROR (C++ enumerator), 66
 brainflow::CustomExitCodes::BOARD_NOT_READY_ERROR (C++ enumerator), 66
 brainflow::CustomExitCodes::BOARD_WRITE_ERROR (C++ enumerator), 66
 brainflow::CustomExitCodes::CLASSIFIER_IS_NOT_PREPARED_ERROR (C++ function), 50, 54, 56, 60, 62, 63, 65
 brainflow::CustomExitCodes::EMPTY_BUFFER_ERROR (C++ enumerator), 66
 brainflow::CustomExitCodes::GENERAL_ERROR (C++ enumerator), 66
 brainflow::CustomExitCodes::INCOMING_MSG_ERROR (C++ enumerator), 66
 brainflow::CustomExitCodes::INITIAL_MSG_ERROR (C++ enumerator), 66
 brainflow::CustomExitCodes::INVALID_ARGUMENTS_ERROR (C++ enumerator), 66
 brainflow::CustomExitCodes::INVALID_BUFFER_SIZE_ERROR (C++ enumerator), 66
 brainflow::CustomExitCodes::JSON_NOT_FOUND_ERROR (C++ enumerator), 65
 brainflow::CustomExitCodes::NO_SUCH_DATA_IN_JSON_ERROR (C++ enumerator), 65
 brainflow::CustomExitCodes::PORT_ALREADY_OPEN_ERROR (C++ enumerator), 65
 brainflow::CustomExitCodes::SET_PORT_ERROR (C++ enumerator), 65
 brainflow::CustomExitCodes::STATUS_OK (C++ enumerator), 65
 brainflow::CustomExitCodes::STREAM_ALREADY_RUN_ERROR (C++ enumerator), 65
 brainflow::CustomExitCodes::STREAM_THREAD_ERROR (C++ enumerator), 67
 brainflow::CustomExitCodes::STREAM_THREAD_IS_NOT_RUNNING_ERROR (C++ enumerator), 66
 brainflow::CustomExitCodes::SYNC_TIMEOUT_ERROR (C++ enumerator), 67
 brainflow::CustomExitCodes::UNABLE_TO_OPEN_PORT_ERROR (C++ enumerator), 66
 brainflow::CustomExitCodes::UNSUPPORTED_BOARD_ERROR (C++ enumerator), 67
 brainflow::CustomExitCodes::UNSUPPORTED_BOARD_FOR_XAND_METRIC_COMBINATION_ERROR (C++ enumerator), 66
 brainflow::DetrendOperations (C++ enum),
 brainflow::DetrendOperations::CONSTANT (C++ enumerator), 67
 brainflow::DetrendOperations::LINEAR (C++ enumerator), 67
 brainflow::DetrendOperations::NONE (C++ enumerator), 67
 brainflow::FilterTypes (C++ enum), 61, 67
 brainflow::FilterTypes::BESSEL (C++ enumerator), 67
 brainflow::FilterTypes::BUTTERWORTH (C++ enumerator), 67
 brainflow::FilterTypes::CHEBYSHEV_TYPE_1 (C++ enumerator), 67
 brainflow::IpProtocolType (C++ enum), 62, 67
 brainflow::IpProtocolType::NONE (C++ enumerator), 67
 brainflow::IpProtocolType::TCP (C++ enumerator), 67
 brainflow::IpProtocolType::UDP (C++ enumerator), 67
 brainflow::LogLevels (C++ enum), 63, 65
 brainflow::LogLevels::LEVEL_CRITICAL (C++ enumerator), 65
 brainflow::LogLevels::LEVEL_DEBUG (C++ enumerator), 65
 brainflow::LogLevels::LEVEL_ERROR (C++ enumerator), 65
 brainflow::LogLevels::LEVEL_INFO (C++ enumerator), 65
 brainflow::LogLevels::LEVEL_OFF (C++ enumerator), 65
 brainflow::LogLevels::LEVEL_TRACE (C++ enumerator), 65
 brainflow::LogLevels::LEVEL_WARN (C++ enumerator), 65
 brainflow::WindowFunctions (C++ enum), 64, 67
 brainflow::WindowFunctions::BLACKMAN_HARRIS (C++ enumerator), 67
 brainflow::WindowFunctions::HAMMING (C++ enumerator), 67
 brainflow::WindowFunctions::HANNING (C++ enumerator), 67
 brainflow::WindowFunctions::NO_WINDOW (C++ enumerator), 67
 BrainFlowException (C++ class), 46

(C++ function), 46
 BrainFlowException::BrainFlowException
 (C++ function), 46
 BrainFlowException::exit_code (C++ member), 47
 BrainFlowException::what (C++ function), 46

D

DataFilter (C++ class), 44
 DataFilter::detrend (C++ function), 45
 DataFilter::disable_data_logger (C++ function), 44
 DataFilter::enable_data_logger (C++ function), 44
 DataFilter::enable_dev_data_logger (C++ function), 44
 DataFilter::get_avg_band_powers (C++ function), 46
 DataFilter::get_band_power (C++ function), 46
 DataFilter::get_nearest_power_of_two (C++ function), 45
 DataFilter::get_psd (C++ function), 45
 DataFilter::get_psd_welch (C++ function), 46
 DataFilter::perform_bandpass (C++ function), 44
 DataFilter::perform_bandstop (C++ function), 44
 DataFilter::perform_downsampling (C++ function), 44
 DataFilter::perform_fft (C++ function), 45
 DataFilter::perform_highpass (C++ function), 44
 DataFilter::perform_ifft (C++ function), 45
 DataFilter::perform_inverse_wavelet_transform (C++ function), 45
 DataFilter::perform_lowpass (C++ function), 44
 DataFilter::perform_rolling_filter (C++ function), 44
 DataFilter::perform_wavelet_denoising (C++ function), 45
 DataFilter::perform_wavelet_transform (C++ function), 44
 DataFilter::read_file (C++ function), 46
 DataFilter::set_log_file (C++ function), 44
 DataFilter::write_file (C++ function), 46

M

MLModel (C++ class), 47
 MLModel::~~MLModel (C++ function), 47
 MLModel::disable_ml_logger (C++ function), 47

MLModel::enable_dev_ml_logger (C++ function), 47
 MLModel::enable_ml_logger (C++ function), 47
 MLModel::MLModel (C++ function), 47
 MLModel::predict (C++ function), 47
 MLModel::prepare (C++ function), 47
 MLModel::release (C++ function), 47
 MLModel::set_log_file (C++ function), 47